



Functional Requirement Specification Document

Version	5.33	
Labguru version	6.37	
Author	Carmel Oren	
Approved by	Gal Haase	
Approved at	12-Aug-2025	

.1	Dashboard (DSH)	9
1.1	'My Daily View' tab	9
1.2	'Recent Activities' tab	9
1.3	'Calendar' tab	9
1.4	'Whiteboard' tab	10
1.5	Section removed	10
1.6	'New Experiment' Popup	11
1.7	'Actions require your attention' boxes	12
1.8	'Dashboards' tab	15
2.	E-Notebook	15
2.1	Labguru project module (PRJ)	15
2.1.1	Projects index page	15
2.1.2	Project show page	16
2.1.3	Archived project	20
2.2	Labguru Project Gantt Chart module(GNT)	21
2.3	Labguru Collaborations module(COL)	25
2.4	Labguru Experiment module(EXP)	27
2.4.1	Experiments Index Page	27
2.4.2	Experiment show page	29
2.4.3	Flagging experiments	32
2.4.4	Actions in experiment sections	33
2.4.5	Experiment with locked format	35
2.4.6	Experiment History Report	36
2.5	Labguru Protocols module (PRT)	37
2.5.1	Protocols index page	37
2.5.2	Protocol show page	38
2.5.3	Archiving protocols	41
2.5.4	Actions in protocol sections	42
2.5.5	External protocols	43
2.6	Data elements(DTL)	43
2.6.1	General	43
2.6.2	Compound element	44
2.6.3	Reaction element	44
2.6.4	Text element	47
2.6.5	Spreadsheet element	47
2.6.6	Steps element	48
2.6.7	Plate element	49
2.6.8	Labguru Form element	54
2.6.9	Samples and reagents element	60
2.6.10	Attachments element	68
2.6.11	Code element	69
2.6.12	Equipment element	70
2.6.13	Date Range element	73

2.6.14	Sketch element	74
2.6.15	Restore element to previous version	74
2.7	Labguru Samples basket module(SMB)	75
2.8	Folders module (FLD)	79
2.8.1	Folder show page	79
2.8.2	Experiments tab	82
2.8.3	Folders tab	83
2.8.4	Figures tab.....	83
3.	Knowledgebase.....	84
3.1	Labguru Datasets module(DTS)	84
3.1.1	'Datasets' index page	84
3.1.2	Dataset show page.....	85
3.2	Labguru SOPs module(SOP)	89
3.2.1	SOP index page.....	89
3.2.2	SOP show page.....	89
3.2.3	Archiving SOP	91
3.2.4	Gantt chart in SOP.....	92
3.3	Labguru Recipes module(RCP)	93
3.3.1	Recipes index page.....	93
3.3.2	Recipe show page	94
3.4	Labguru Papers module(PAP)	97
3.4.1	Papers index page	97
3.4.2	Paper show page.....	97
3.5	Labguru Documents module(DOC)	99
3.5.1	Documents index page.....	99
3.5.2	Document show page	100
3.5.3	Archiving documents	102
3.6	Labguru Reports module(REP)	103
3.6.1	Reports index page	103
3.6.2	Report show page	104
3.6.3	Archiving reports.....	107
3.7	Labguru Image Bank module(IMG)	108
3.7.1	Images index page.....	108
3.7.2	Image show page	108
4.	Inventory	109
4.1	Biocollections(BIO)	109
4.1.1	Index page.....	109
4.1.2	Biocollection's show page.....	114
4.1.3	Biocollection's show pages tabs.....	116
4.1.4	Inventory item 'Add new/Edit item' page.....	121
4.1.5	Labguru antibody module (ANT)	122
4.1.6	Section removed	123
4.1.7	Labguru Bacteria module	123
4.1.8	Labguru Botany module(BTN)	124
4.1.9	Labguru Cell lines module (CEL).....	125
4.1.10	Labguru Compounds Module(COM).....	125
4.1.11	Labguru Flies module (FLY)	127

4.1.12	Labguru Fungi module(FNG)	127
4.1.13	Labguru Genes module(GEN)	128
4.1.14	Labguru Lipids module(LIP)	128
4.1.15	Labguru Plasmids module(PLA)	128
4.1.16	Labguru Primer module(PRM)	129
4.1.17	Labguru Protein module (PRO)	129
4.1.18	Labguru Sequences module(SEQ)	130
4.1.19	Labguru Tissue module(TIS)	130
4.1.20	Labguru Worms module (WRM)	131
4.1.21	Labguru yeasts module (YST)	131
4.1.22	Labguru Viruses module(VIR)	131
4.1.23	Labguru Zebrafishes module(ZEB)	132
4.1.24	Labguru Rodents module	132
4.1.25	Archived item	136
4.2	Labguru Consumables module(CONS)	137
4.2.1	Consumables index page	137
4.2.2	Consumable show page	138
4.2.3	'Add new Material' page	138
4.3	Labguru Shopping List module (SHP)	138
4.3.1	'Add item to shopping list' popup -from inventory item show page	138
4.3.2	'Shopping list' index page	139
4.3.3	Section Removed	141
4.3.4	Order History page	141
4.3.5	Service orders page	142
4.4	Labguru Companies module(CMP)	143
4.5	Customize collection page(CMZ)	145
4.5.1	Pre-defined list	145
4.5.2	Custom fields	146
4.5.3	Derived Collections	147
5.	Storage & Equipment module	148
5.1	Labguru Storages module (STG)	148
5.1.1	Storages page	148
5.1.2	Section removed	150
5.2	Labguru Boxes module(BOX)	150
5.2.1	Boxes index page	150
5.2.2	Box show page	152
5.3	Labguru Stocks module (STK)	156
5.3.1	Stocks index page	156
5.3.2	Consumed stocks index page	158
5.3.3	Stock show page	159
5.3.4	Stocks audit	164
5.3.5	Low stock alert	166
5.3.6	Aliquoting stocks	167
5.3.7	Add/Edit stock page	170
5.3.8	Stocks customization page	173
5.3.9	Stock Pooling	173
5.3.10	Dispatch	176
5.4	Labguru Plates module(PLT)	177
5.4.1	Plates index page	177

5.4.2	Plate show page	178
5.5	Labguru Equipment Module(EQU)	179
5.5.1	Equipment index page	179
5.5.2	Equipment show page.....	182
5.5.3	Archived equipment.....	183
5.5.4	Labguru Maintenance Module.....	183
5.5.5	Labguru Scheduler Module(SCH)	187
5.5.6	Creating equipment events when concurrent events are not allowed	192
5.5.7	Add/Edit Equipment page	193
6.	System functions.....	195
6.1	Labguru Linked Resources module(LRS)	195
6.2	Labguru Tasks module(TSK)	196
6.2.1	Tasks page	196
6.2.2	From inventory item show page	198
6.2.3	From task show page	198
6.3	Labguru Tags module(TAG)	199
6.3.1	From item show page	199
6.3.2	From item index page	199
6.3.3	From tags Index page	200
6.4	Labguru Comments module(CMN)	200
6.5	Labguru Notes module(NTE)	201
6.5.1	Notes in projects and folders	201
6.5.2	Note show page	201
6.6	Global search(GLB)	202
6.6.10	203	
6.7	'Recently viewed' list(RCL)	203
6.8	Labguru Events module(EVT)	203
6.8.1	Calendar Page	203
6.8.2	Adding/editing events	205
6.8.3	Event show page	206
6.9	Section removed	207
6.10	Print to PDF	207
6.11	Label Printing(LBP)	208
6.11.1	General.....	208
6.11.2	Label templates settings	208
6.11.3	Template advanced settings	209
6.12	System notifications	210
7.	Signing and Witnessing (SIG).....	210
7.1	Signed page	210
7.2	Signed experiment	212
7.3	Signed protocol and SOP	212
7.4	Witness a signed page	212
7.5	User authentication in signature actions	214

8.	Account and user management	214
8.1	Labguru Billing module(BIL)	214
8.1.1	Index page.....	214
8.1.2	External credit card link	216
8.2	Account settings module(ACS)	216
8.2.1	General settings	216
8.2.2	Roles management	217
8.2.3	New collections and collections settings	219
8.2.4	Data export	219
8.2.5	Account personalization.....	220
8.2.6	Signature process	220
8.2.7	Account LOVs	221
8.2.8	Forms Settings (Reason for Change)	222
8.2.9	PDF Settings	223
8.2.10	Experiments Settings.....	223
8.2.11	Requests Settings	224
8.3	Privileges in Labguru(PRV)	224
8.4	Labguru Login/Logout module	224
8.4.1	Labguru Login(SGN).....	224
8.4.2	Labguru First Login(SGN).....	226
8.4.3	Labguru Logout(SGN).....	226
8.5	Password management(PAS)	226
8.5.1	Changing password	226
8.5.2	Section removed	227
8.6	Labguru members module(MEM)	227
8.6.1	My profile.....	227
8.6.2	Members list	229
8.6.3	Member show page	231
8.7	Labguru Weekly Digest Email module (WKG)	232
8.8	Labguru Teams module(TMS)	233
8.8.1	General.....	233
8.8.2	Using teams.....	234
8.9	Section removed (SRV)	235
8.10	Section removed(LGT)	235
8.11	Section removed	235
8.12	Support menu	235
9.	Department account(DEP)	236
9.1	General	236
9.1.8	237	
9.2	Sharing equipment	237
9.2.1	On Sharing Lab account	237
9.2.2	On Department account	239
9.2.3	On Subscribing Lab account	242
9.3	Shared protocols	243

10. Labguru import/export (IMP)	244
10.1 Section Removed	244
10.2 Import/Update	244
10.2.1 General	244
10.2.2 Import in collections	246
10.2.3 Import Stocks	248
10.2.4 Import Equipment	250
10.3 Export	251
10.4 Bulk-Update	251
.11 Section removed	251
12. Labguru Workflow Editor module(WFE)	252
12.1 General	252
12.2 Flows index page	252
12.3 Flow page - Editor tab	253
12.4 Flow page - Runs tab	255
12.5 Triggers	256
12.6 Steps	256
12.6.1 Conditions	256
12.6.2 Scripter	257
12.7 Code templates	258
12.8 Exceptions	258
12.9 Flow page - Variables tab	259
13. Labguru Audit trail module(ADT)	259
13.1 General Audit Trail page	259
13.1.9	261
13.2 Diff column	261
14. Components	262
14.1 DNA viewer	262
14.1.1 General	262
14.1.2 Actions menu	264
14.1.3 Editing sequence	265
14.1.4 Quick Search in viewer	265
14.2 EasyEdit - Desktop application	266
14.2.1 EasyEdit setup and access	266
14.2.2 EasyEdit connection	266
14.2.3 Editing and uploading attachments	267
14.2.4 Logs and Information	268
14.2.5 EasyEdit Geneious Prime plug-in	269
14.3 Small Molecule Drawer(SMD)	270
14.3.1 General	270
14.3.2 Drawing Bonds	273

14.3.3	Drawing Atoms.....	273
14.3.4	Marking S-Groups	274
14.3.5	Template toolbar	274
15.	Requests	274
15.1	Labguru Request Templates module	274
15.1.1	Request Templates index page	274
15.1.2	Request Templates show page	276
15.1.3	Actions in Request templates section	277
15.2	Labguru Request module	278
15.2.1	Request Index Page.....	278
15.2.2	Request Show Page.....	279
15.2.3	Navigating Request Flow	281
15.2.4	Actions in Request templates section	282
16.	Batch Management.....	283
16.1	Labguru Processes module	283
16.1.1	Processes index page	283
16.1.2	Processes show page	283
16.1.3	Add/Edit Process	290
16.2	Labguru Batch Records module	296
16.2.1	Batch Records index page	296
16.2.2	Batch Records show page	296
17.	Customer Portal	303
17.1	General	303
17.1.1	Customer Portal First Login.....	303
17.1.2	Customer Portal Login/logout.....	303
17.2	User roles and access controls	304
17.2.1	Super Administrator index page	304
17.2.2	Super Administrator show page.....	305
17.2.3	Customers index page.....	305
17.2.4	Customers show page	306
17.3	Sample management	309
17.3.1	Samples index page.....	309
17.3.2	Samples show page.....	310
17.3.3	Shared Samples	310
Change Log		311
Appendixes		311

1. Dashboard (DSH)

1.1 'My Daily View' tab

Number	1.1.1
Requirement	User will be able to see actions that were made in the last few hours
Note	

Number	1.1.2
Requirement	User will be able to navigate between dates and see other actions
Note	

Number	1.1.3
Requirement	User will be able to add events directly through the event board
Note	By clicking on the event board and filling event's info

1.2 'Recent Activities' tab

Number	1.2.1
Requirement	User will be able to see every action that was made in the last 7 days
Note	• Every action will be followed by member's name, item's name and date • User will be redirected to the item show page, by clicking on the item's name • User will be redirected to member show page, by clicking on the members name

Number	1.2.2
Requirement	User will be able to filter activities by member
Note	By clicking on the 'Filter' link

1.3 'Calendar' tab

Please look at 'Labguru events module'

1.4 'Whiteboard' tab

Number	1.4.1
Requirement	User will be able to write a post on the whiteboard
Note	• By clicking on the 'New' button • By filling the subject and description of the post • By clicking 'Post' • Every post will be followed by member's name and date of the post

Number	1.4.2
Requirement	User will be able to delete a post from the whiteboard
Note	• By clicking on the 'x' icon

Number	1.4.3
Requirement	User will be able to archive a post from the whiteboard
Note	• By clicking on the '√' icon

Number	1.4.4
Requirement	User will be able to edit a post
Note	• By clicking on the post title • By clicking 'Edit'

Number	1.4.5
Requirement	User will be able to enter the post show page
Note	• By clicking on its title

Number	1.4.6
Requirement	User will be able to comment on a post
Note	

Number	1.4.7
Requirement	User will be able to see the number of archived posts
Note	• In the archived posts section

1.5 Section removed

Number	1.5.1
Requirement	Requirement removed
Note	

1.6 'New Experiment' Popup

Number	1.6.1
Requirement	User will be able to add new experiment
Note	• By clicking on the 'New Experiment' button • User will be able to name the experiment • User will be able to select project, folder and protocol from list/search by typing the first letter • User will be able to add new experiment from a protocol • When protocol is not selected, the experiment will be empty

Number	1.6.2
Requirement	Requirement removed
Note	

Number	1.6.3
Requirement	Requirement removed
Note	

Number	1.6.4
Requirement	User will be required to select a name for experiment
Note	By default, the experiment name will be:"{todays_date} - {member name}\s experiment" with daily serial incrementing number

Number	1.6.5
Requirement	Removed requirement
Note	

Number	1.6.6
Requirement	Requirement removed
Note	

Number	1.6.7
Requirement	User will be able to add new protocol
Note	By choosing the option from the drop-down menu

Number	1.6.8
Requirement	User will be able to add new document
Note	By choosing the option from the drop-down menu

Number	1.6.9
Requirement	User will be able to add new project
Note	By choosing the option from the drop-down menu

1.7 'Actions require your attention' boxes

1.7.1 General

Number	1.7.1.1
Requirement	User will be able to see all actions required their attention when clicking on "Actions require your attention"
Note	The total number of actions will be displayed

Number	1.7.1.2
Requirement	User will not be able to see experiments for witness and stock alert boxes if privileges are turned off
Note	

Number	1.7.1.3
Requirement	User will be able to scroll horizontally to view all boxes
Note	

1.7.2 'Tasks' box

Number	1.7.2.1
Requirement	User will see counter zero and "No open tasks" title will be shown if he has no open tasks
Note	

Number	1.7.2.2
Requirement	User will see counter with a number other than zero and "Open tasks" title will be shown if he has any open task
Note	

Number	1.7.2.3
Requirement	User will be able to move to the tasks index page
Note	By clicking on the 'Go to tasks' button

Number	1.7.2.4
Requirement	User will see counter with a number other than zero and "Overdue tasks" title if he has any overdue task
Note	

1.7.3 ‘Unattached images & attachments’ box

Number	1.7.3.1
Requirement	User will see the number of unattached images and unattached attachments
Note	

Number	1.7.3.2
Requirement	User will be able to go to the images index page
Note	By clicking on the ‘Go to images’ button

Number	1.7.3.3
Requirement	User will be able to go to the attachments index page
Note	By clicking on the ‘Go to attachments’ button

1.7.4 ‘Low stocks alert’ box

Number	1.7.4.1
Requirement	User will be able to go to the shopping list page
Note	By clicking on the ‘Go to the shopping list’ button

Number	1.7.4.2
Requirement	User will be able to see a stock alert when low levels conditions are met
Note	

1.7.5 ‘Signed experiments need to be witnessed’ box

Number	1.7.5.1
Requirement	User will see signed experiments that need witnessing
Note	• User with “Can witness” privileges will see notification for all experiments that are part of projects with default signature flow • User who is listed as witness in a custom signature flow will see notification for all experiments that are part of projects with that signature flow • User will see indication to the number of signed experiments that need witnessing that he can witness • User will see notifications only for experiments that are part of projects user has access to

Number	1.7.5.2
Requirement	User will be able to go to the experiments index page
Note	By clicking on the 'Go to experiments'

1.7.6 'New arrived stocks' box

Number	1.7.6.1
Requirement	User will see the number of the arrived stocks
Note	

Number	1.7.6.2
Requirement	User will be able to go to the dispatched events page
Note	By clicking on the 'Go to dispatched/arrived stocks' button

1.7.7 'Experiments with invalid forms' box

Number	1.7.7.1
Requirement	User will see a notification for experiments with invalid forms owned by the user
Note	- User will see indication to number of experiments with invalid forms - When an invalid form is deleted from an experiment, the experiment is excluded from the list

Number	1.7.7.2
Requirement	User will be able to go to the experiments index page
Note	By clicking on the 'Go to invalid experiments' button

1.7.8 'Expired stocks alert' box

Number	1.7.8.1
Requirement	User will be able to see the number of the expired stocks
Note	

Number	1.7.8.2
Requirement	User will be able to view the expired stocks
Note	By clicking on the 'View expired stocks' button

1.7.9 'Equipment required maintenance' box

Number	1.7.9.1
Requirement	User will be able to see the number of the equipment that require maintenance
Note	

Number	1.7.9.2
Requirement	User will be able to view the equipment that require maintenance
Note	By clicking on the 'Go to equipment' page

1.8 'Dashboards' tab

Number	1.8.1
Requirement	User will be able to view all dashboards that were shared to Labguru
Note	

2. E-Notebook

2.1 Labguru project module (PRJ)

2.1.1 Projects index page

Number	2.1.1.1
Requirement	User will see list of projects divided to: • Projects the user is the owner • Projects of other members that user has access to
Note	

Number	2.1.1.2
Requirement	User will be able to filter projects by: • Project owner • Tags • Project creation date • Project last update date
Note	By using the 'More filters' link

Number	2.1.1.3
Requirement	User will be able to add a new project
Note	By clicking on 'New Project'

Number	2.1.1.4
Requirement	User will be able to add a new experiment
Note	• By opening the drop-down menu next to 'New Project' button and selecting new experiment • User will be able to name the experiment • User will be able to select project, folder and protocol from list/search by typing the first letter

Confidential

	When protocol is not selected, the experiment will be empty
--	---

Number	2.1.1.5
Requirement	Requirement removed
Note	

Number	2.1.1.6
Requirement	User will be able to add a new protocol
Note	By opening the drop-down menu next to ‘New Project’ button

Number	2.1.1.7
Requirement	User will be able to add a new document
Note	By opening the drop-down menu next to ‘New Project’ button

Number	2.1.1.8
Requirement	User will not see archived projects in list
Note	

2.1.2 Project show page

Number	2.1.2.1
Requirement	User will be able to edit title of project
Note	

Number	2.1.2.2
Requirement	User will be able to assign the project to another member
Note	

Number	2.1.2.3
Requirement	User will be able to add a new experiment
Note	- By clicking on ‘New Experiment’ - User will be able to name the experiment - User will be able to select project, folder and protocol from list/search by typing the first letter - When protocol is not selected, the experiment will be empty

Number	2.1.2.4
Requirement	User will be able to duplicate a project

Confidential

Note	- By clicking on 'Duplicate' in the 'More' drop-down menu - Linked resources, tags, papers, notes, documents, folders, experiments and other data elements will be duplicated as well
------	--

Number	2.1.2.5
Requirement	User will be able to print the project show page to PDF file
Note	By clicking on 'Print'

Number	2.1.2.6
Requirement	User will be able to export project as zip
Note	Zip file will include all folders, experiments and attachments

Number	2.1.2.7
Requirement	Requirement removed
Note	

Number	2.1.2.8
Requirement	User will be able to view version history of project
Note	By clicking on 'Version history' in the 'More' menu

Number	2.1.2.9
Requirement	User will be able to download attachments added to experiments in project
Note	From 'Attachments' tab

Number	2.1.2.10
Requirement	User will be able to view images added to experiments in project
Note	From 'Attachments' tab

Number	2.1.2.11
Requirement	User will be able to view papers linked to project
Note	By clicking on 'Key Papers' tab

Number	2.1.2.12
Requirement	User will be able to add new papers that are linked to the project
Note	From 'Key Papers' tab

Confidential

Number	2.1.2.13
Requirement	User will be able to view samples added to experiments in project
Note	By clicking on 'Samples' tab

Number	2.1.2.14
Requirement	User will be able to export samples list as a csv file
Note	From 'Samples' tab

Number	2.1.2.15
Requirement	User will be able to view notes added to project
Note	By clicking on 'Notes' tab

Number	2.1.2.16
Requirement	User will be able to add new notes to project
Note	From 'Notes' tab

Number	2.1.2.17
Requirement	User will be able to view documents linked to project
Note	By clicking on 'Documents' tab

Number	2.1.2.18
Requirement	User will be able to add new documents linked to project
Note	From 'Documents' tab, by clicking on 'Add Document'

Number	2.1.2.19
Requirement	User will be able to view the project's calendar
Note	By clicking on 'Calendar' tab

Number	2.1.2.20
Requirement	User will be able to view the project's Gantt chart
Note	By clicking on 'Gantt' tab

Number	2.1.2.21
Requirement	User will be able to add folder to project
Note	- By clicking on 'Add Folder' from: - Top right corner of 'Progress' tab - In 'Completed', 'In progress' or 'Planned' sections

Confidential

Number	2.1.2.22
Requirement	User will be able to re-order folders in the folders tree
Note	By dragging folders up/down

Number	2.1.2.23
Requirement	User will be able to set start date of a project
Note	By clicking on the calendar icon and filling a date in the 'Select start date' popup

Number	2.1.2.24
Requirement	User will be able to set email notification on the project start date
Note	By clicking on the calendar icon and filling 'Email notification' and 'Send notification at' to choose who should receive the notification and when notification should be sent

Number	2.1.2.25
Requirement	User will be able to move experiments within the project
Note	• By drag and drop • By default, experiments order will be by creation date (oldest first)

Number	2.1.2.26
Requirement	User will be able to select signature flow to apply to project
Note	Existing, signed experiment that are part of the project will not be affected from changing the signature flow

Number	2.1.2.27
Requirement	User will be able to open experiment in a new tab
Note	By clicking on icon next to the experiment

Number	2.1.2.28
Requirement	User will be able to set which members have access to the project
Note	• By clicking on 'Permissions' in the 'More' menu • User will be able to select 'All' or 'None' from the offered members list • User will be able to select members and/or teams • User without access will not see the project in the projects' list and will not be able to view any page that is part of the project

Confidential

	- User with access will be able to: - Create new folders - Create new experiments - View all experiments - User added to a team with access to project will inherit the team's permission - Only admins and project owner can edit project access
--	--

Number	2.1.2.29
Requirement	User will be able to set experiment default edit permissions to selected members
Note	- By clicking on 'Permissions' in the 'More' menu - By default, only admins and project owner have edit permissions - User will be able to select 'All' or 'None' from the offered members list - User will be able to select members and/or teams - User with default experiment edit permission will be able to edit any experiment in project - User added to a team with access to project will inherit the team's permission - Only admins and project owner can edit default experiment edit permissions

Number	2.1.2.30
Requirement	User will be able to expand/collapse folders
Note	- User preferences will be saved - New folder will collapse by default

2.1.3 Archived project

Number	2.1.3.1
Requirement	User will be able view archived projects by clicking the 'View archived projects' links
Note	

Number	2.1.3.2
Requirement	User will be able to archive project
Note	By clicking on the 'Archive' link via 'More' menu.

Number	2.1.3.3
Requirement	User will not be able to edit archived project
Note	

Number	2.1.3.4
Requirement	User will not be able to add folders/experiments to archived project
Note	

Number	2.1.3.5
Requirement	User with 'Can administer' privilege will be able to re-activate archived project that is owned by another member
Note	

Number	2.1.3.6
Requirement	User will be able to restore archived project
Note	- By clicking on 'Restore' link via 'More' menu - Only for the owner of the project and user with admin privileges

2.2 Labguru Project Gantt Chart module[\(GNT\)](#)

2.2.1 Chart view

Number	2.2.1.1
Requirement	User will see current date or the last day of the project at the center of the view
Note	Whichever comes first

Number	2.2.1.2
Requirement	User will be able to scroll the chart horizontally
Note	The display should end 2 days before the start day and 2 days after the end date

Number	2.2.1.3
Requirement	In the Gantt view tab User will see the project, it's folders and experiments as Gantt items differentiated to three colors
Note	Note removed

Number	2.2.1.4
Requirement	User will be able to go to project/folder/experiment show page
Note	By clicking on the associated Gantt item

Confidential

Number	2.2.1.5
Requirement	User will be able to drag experiment left and right to change its start date
Note	

Number	2.2.1.6
Requirement	User will be able to expand and contract experiment width to change its duration
Note	When changing duration in Gantt item, end date (in experiment page) will be updated accordingly

Number	2.2.1.7
Requirement	User will see signed experiments in a darker color than unsigned experiment
Note	

Number	2.2.1.8
Requirement	User will not be able to change start date and duration of signed experiments
Note	

Number	2.2.1.9
Requirement	User will see completed folder in a darker color than un-completed folders
Note	

Number	2.2.1.10
Requirement	User will see a representation of the progress of a folder as a change of color in the Gantt item
Note	Progression is measured as number of signed experiments out of number of total experiments in folder

Number	2.2.1.11
Requirement	User will see a representation of the progress of a project as a change of color in the Gantt item
Note	Progression is measured as average progression of folders in the project

Number	2.2.1.12
Requirement	User will be able to assign members to experiment/folder
Note	

2.2.2 Information columns

Number	2.2.2.1
Requirement	User will be able to view the name, owner, start date, progress and duration of each Gantt item
Note	

Number	2.2.2.2
Requirement	User will be able to hide/un-hide project/folder content
Note	

Number	2.2.2.3
Requirement	User will see under 'Progress' if an experiment is signed or not
Note	

Number	2.2.2.4
Requirement	User will see under 'Progress' a percentage representation of the number of signed experiments out total experiments in the folder
Note	

Number	2.2.2.5
Requirement	User will see under 'Progress' a percentage representation of average progress of folders in the project
Note	

2.2.3 Dependencies

Number	2.2.3.1
Requirement	User will be able to connect two experiments from the same folder
Note	

Number	2.2.3.2
Requirement	User will be able to remove a dependency link
Note	

Confidential

Number	2.2.3.3
Requirement	In end to start - If user move the end date of the parent to be after the start of the dependent, the dependent start date will move to when the parent's end date is
Note	

Number	2.2.3.4
Requirement	In end to start - User will not be able to move the dependent's start date to be before the parent's end date
Note	

Number	2.2.3.5
Requirement	In start to start - If user move the start date of the parent to be after the start of the dependent, the dependent start date will move to when the parent's start date is
Note	

Number	2.2.3.6
Requirement	In start to start - User will not be able to move the dependent's start date to be before the parent's start date
Note	

Number	2.2.3.7
Requirement	In start to end - If user move the start date of the parent to be after the end of the dependent, the dependent end date will move to when the parent's start date is
Note	

Number	2.2.3.8
Requirement	In start to end - User will not be able to move the dependent's end date to be before the parent's start date
Note	

Number	2.2.3.9
Requirement	In end to end - If user move the end date of the parent to be after the end of the dependent, the dependent end date will move to when the parent's end date is
Note	

Number	2.2.3.10
Requirement	In end to end - User will not be able to move the dependent's end date to be before the parent's end date

Confidential

Note	
------	--

Number	2.2.3.11
Requirement	User will not be able to sign a dependent experiment if the parent is not signed
Note	

Number	2.2.3.12
Requirement	User will not be able to change a dependent experiment's start date/duration from the experiment page if the parent is signed
Note	

Number	2.2.3.13
Requirement	User will not be able to move a dependent experiment to another folder if the parent is signed
Note	

Number	2.2.3.14
Requirement	User will be able to delete experiment/folder that has dependencies connection
Note	When user delete an experiment/folder it deletes all dependencies links to it

Number	2.2.3.15
Requirement	User will be able to connect two different folders, ONLY if the second (dependent) folder don't have signed experiments with in it
Note	

Number	2.2.3.16
Requirement	User will not be able to sign any experiment on the second (dependent) folder before all experiments in previous folder are signed
Note	

2.3 Labguru Collaborations module(COL)

2.3.1 Share Project

Number	2.3.1.1
Requirement	Project owner will be able to share the project with a collaborator

Note	• Project owner is required to add the collaborator name and e-mail • Collaborator will receive an email with link to the project
------	--

Number	2.3.1.2
Requirement	Project owner will be able to stop collaboration
Note	• When a collaboration is being terminated, link sent to collaborator is deactivated • If collaborator is in the project when collaboration is stopped, collaborator will be kicked out once moving to a new page in project

2.3.2 Collaborators in a project

Number	2.3.2.1
Requirement	Collaborator will not be able to edit data or perform any actions
Note	The collaboration mode is a read-only mode: • Cannot access section edit mode • Cannot add comments, tags, links, attachments, task • Start new pages • Delete pages

Number	2.3.2.2
Requirement	Collaborator will not be able to download attachments
Note	

Number	2.3.2.3
Requirement	Collaborator will be able to view: • Attachments as thumbnails • Project/folder notes - from 'Notes' tab • Project document list, without access to documents pages - from 'Documents' tab • Key papers list, without access to paper pages - from 'key papers' tab • Figures on the folder page - from 'Figures' tab • Experiment pages
Note	

Number	2.3.2.4
Requirement	Collaborator will not be able to click other Labguru page by clicking on links directing to them
Note	

2.4 Labguru Experiment module^(EXP)

2.4.1 Experiments Index Page

Number	2.4.1.1
Requirement	User will be able to add a new experiment
Note	• By clicking the ‘New Experiment’ button • User will be able to name the experiment • User will be able to select project, folder and protocol from list/search by typing the first letter • When protocol is not selected, the experiment will be empty

Number	2.4.1.2
Requirement	Requirement removed
Note	
Number	2.4.1.3
Requirement	User will be able to add a new protocol
Note	By opening the drop-down menu next to ‘New experiment’ button

Number	2.4.1.4
Requirement	User will be able to add a new document
Note	By opening the drop-down menu next to ‘New experiment’ button

Number	2.4.1.5
Requirement	User will be able to add a new project
Note	By opening the drop-down menu next to ‘New experiment’ button

Number	2.4.1.6
Requirement	User will be able to filter experiments by tags
Note	

Number	2.4.1.7
Requirement	User will be able to change the number of experiments displayed per page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 experiments per page

Confidential

Number	2.4.1.8
Requirement	User will be able to sort experiments by any column in display in ascending/descending order
Note	Filter preferences are saved

Number	2.4.1.9
Requirement	User will be able to filter experiments by any column or combination of columns in display
Note	Filter preferences are saved

Number	2.4.1.10
Requirement	User will be able to choose which columns to view on the index page
Note	• By choosing headers from list in 'Columns' option • User's preference will be saved • possible headers: Created at, Last updated, ID, Experiment name, Project, Owner, start date, End date, Status, Signed, Witnessed

Number	2.4.1.11
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	2.4.1.12
Requirement	User will be able to change the order of the columns
Note	

Number	2.4.1.13
Requirement	User will be able to enter an experiment page by clicking on its name or ID
Note	

Number	2.4.1.14
Requirement	User will see under 'Status' if the experiment was flagged as successful or failed
Note	

Number	2.4.1.15
Requirement	User will see a purple check mark if the experiment: • is signed (under 'Signed' column) • is witnessed (under 'Witness column)
Note	

Number	2.4.1.16
Requirement	User will see under 'End date' a calculation of the end date by taking the experiment's start date and duration
Note	End date for experiments without duration will be in the same day as the start date

Number	2.4.1.17
Requirement	Requirement removed
Note	

Number	2.4.1.18
Requirement	User will see list for signed experiments that need witnessing
Note	- By filtering experiments list by 'Experiments to witness' - If has 'Can witness' privileges, user will see all signed experiments from projects with default signature flow - If part of a signature group, user will see all signed experiments from projects with that signature flow

Number	2.4.1.19
Requirement	User will see a list for experiments with invalid forms
Note	- By filtering experiments list by 'Experiments with invalid form' - Only users with edit permissions to the experiment

Number	2.4.1.20
Requirement	Requirement removed
Note	

2.4.2 Experiment show page

Number	2.4.2.1
Requirement	User will be able to change experiment name
Note	

Number	2.4.2.2
Requirement	User will be able to assign experiment to another member
Note	

Confidential

Number	2.4.2.3
Requirement	User will be able to change experiment start date & time
Note	

Number	2.4.2.4
Requirement	User will be able to edit experiment duration
Note	Experiment duration can be only positive integer

Number	2.4.2.5
Requirement	User will be able to tag experiment
Note	

Number	2.4.2.6
Requirement	User will be able to sign experiment
Note	- Only user with editing permission to experiment will be able to sign it - When another user is editing experiment page, the sign function will be disabled, and a message will appear on hover indicating the page is being edited

Number	2.4.2.7
Requirement	User will be able to link experiment to other pages
Note	- By clicking the linked resources button on the 'More' menu - By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	2.4.2.8
Requirement	User will be able to view version history
Note	

Number	2.4.2.9
Requirement	User will be able to delete experiment
Note	Only user with edit permissions to the experiment

Number	2.4.2.10
Requirement	User will be able to load sections from protocols
Note	By clicking on the 'add section' link and choosing protocol to load from

Number	2.4.2.11
Requirement	User will be able to add a new section
Note	By clicking on the 'add section' link

Number	2.4.2.12
Requirement	User will be able to download all attachments added to experiment
Note	

Number	2.4.2.13
Requirement	User will be able to add a task connected to experiment
Note	Task assignee can be only a member who has view permissions to the project the experiment is part of

Number	2.4.2.14
Requirement	Experiment owner will be able to set edit permissions for other members/teams
Note	- By default, only admins, project owner and experiment owner have edit permissions - User without edit permissions will not be able to perform any of the following actions: Renaming, assigning, changing start date and duration, tagging, signing, adding links, deleting, loading/adding sections, adding tasks, downloading attachments, flagging, modifying the sections and/or modifying elements within the sections

Number	2.4.2.15
Requirement	User will be able to duplicate experiment
Note	- A pop-up notification will be displayed when duplication process has started/ completed - Notification will include experiment's location and link to the duplicate page - Usage record information will not be saved to the duplicate page

Number	2.4.2.16
Requirement	User will be able to save experiment as protocol
Note	- By clicking on the 'Save as protocol' link in the 'more' menu - All the data in the experiment will be saved in the protocol

Confidential

	Usage record information will not be saved to the protocol page
--	---

Number	2.4.2.17
Requirement	Requirement removed
Note	

Number	2.4.2.18
Requirement	User will be able to print experiment to pdf file
Note	

Number	2.4.2.19
Requirement	User will be able to add a comment to experiment
Note	

Number	2.4.2.20
Requirement	User will be able to view: - Experiment creator - Creation date - Last update date
Note	

Number	2.4.2.21
Requirement	User will be able to export experiment as a Word document
Note	By clicking on 'More' and then 'Export to Word'

Number	2.4.2.22
Requirement	User will be able to move experiment to another folder
Note	Only users with edit permissions to the experiment

Number	2.4.2.23
Requirement	User will be able to see experiment calendar
Note	Experiment calendar will show events created from date range elements that were added to the experiment

2.4.3 Flagging experiments

Number	2.4.3.1
Requirement	User will be able to flag an experiment as successful/failed/custom flag

Note	
------	--

Number	2.4.3.2
Requirement	User will be able to un-flag a flagged experiment
Note	

Number	2.4.3.3
Requirement	User will not be able to flag/un-flag a signed experiment
Note	

Number	2.4.3.4
Requirement	User will be able to see experiment flag in following lists: • Project show page • Experiments index page • Folder page
Note	

Number	2.4.3.5
Requirement	User will be able to see list of the account active experiment flags
Note	

2.4.4 Actions in experiment sections

Number	2.4.4.1
Requirement	User will be able to change experiment section name
Note	

Number	2.4.4.2
Requirement	User will be able to add to experiment section the following elements: • Text • Spreadsheet • Steps • Samples • Plate • Equipment • Sketch • Compound • Reaction • Date range • Attachments

Confidential

Note	Date range cannot be added to ‘Description’ and ‘Results’ sections
------	--

Number	2.4.4.3
Requirement	User will be able to delete an element from experiment section
Note	

Number	2.4.4.4
Requirement	User will be able to copy an experiment section to report
Note	User can search for report by typing the first 3 letters of its name in the search bar or creating a new report

Number	2.4.4.5
Requirement	Requirement removed
Note	

Number	2.4.4.6
Requirement	User will be able to move experiment section up/down
Note	

Number	2.4.4.7
Requirement	User will be able to delete an experiment section
Note	

Number	2.4.4.8
Requirement	User will be able to save an experiment section
Note	Saving a section will also close section edit mode

Number	2.4.4.9
Requirement	User will be able to add page break between elements in a section
Note	- The page break mark will be visible on screen - When experiment printed to pdf, the element after the page break mark will be printed on top a new page

Number	2.4.4.10
Requirement	User will be able to copy an experiment section to experiment

Note	User can search for experiment by typing the first 3 letters of its name in the search bar or by adding to the same experiment
------	--

Number	2.4.4.11
Requirement	User will be able to move elements in experiment section up/down
Note	By dragging the element to the desired position

2.4.5 Experiment with locked format

Number	2.4.5.1
Requirement	User will not be able to add/rename/delete/move sections in experiment with locked format
Note	

Number	2.4.5.2
Requirement	User will not be able to add/rename/delete/move element in experiment with locked format
Note	

Number	2.4.5.3
Requirement	User will not be able to add inline comments
Note	

Number	2.4.5.4
Requirement	User will not be able to unlock/lock format from experiment page
Note	

Number	2.4.5.5
Requirement	User will see in experiment page if it's format is locked
Note	

Number	2.4.5.6
Requirement	User will be able to change experiment format when experiment was created by duplicating experiment with locked format
Note	

Number	2.4.5.7
Requirement	User will be able to remove ‘Add to page view’ in experiment with locked format
Note	

2.4.6 Experiment History Report

2.4.6.1 General

Number	2.4.6.1.1
Requirement	User will be able to print experiment history report
Note	- By clicking on ‘More’, ‘Print’ and then check the ‘History report’ check box - Popup message to download PDF file will be displayed

Number	2.4.6.1.2
Requirement	User will be able to view the required entries information: - Timestamp - Creation date - User - Action - Event - Diff - Object
Note	Entry made by auto-save action will be displayed in ‘Action’ column

Number	2.4.6.1.3
Requirement	User will be able to see headers and footers in history report
Note	

Number	2.4.6.1.4
Requirement	After witnessing, user will be able to see history report as pdf attachment
Note	

2.5 Labguru Protocols module (PRT)

2.5.1 Protocols index page

Number	2.5.1.1
Requirement	User will be able to add a new protocol
Note	

Number	2.5.1.2
Requirement	User will be able to filter protocols by tags
Note	

Number	2.5.1.3
Requirement	User will be able to search for a protocol by name
Note	

Number	2.5.1.4
Requirement	User will be able to change the number of protocols he sees on single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 documents per page

Number	2.5.1.5
Requirement	User will be able to sort protocols by any column in display in ascending/descending order
Note	

Number	2.5.1.6
Requirement	User will be able to filter protocols by any combination of columns in display
Note	

Number	2.5.1.7
Requirement	User will be able to choose which columns to view on the index page
Note	- By choosing headers from list in 'Columns' option - User's preference will be saved

Number	2.5.1.8
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	2.5.1.9
Requirement	User filter/sort selection will be saved when user move to other pages in the collection
Note	

Number	2.5.1.10
Requirement	User will be able to see list of active protocols that user has access to view
Note	

2.5.2 Protocol show page

Number	2.5.2.1
Requirement	User will be able to change protocol name
Note	

Number	2.5.2.2
Requirement	User will be able to edit a protocol source
Note	Protocol source can be link to page when has value in 'url' field. If 'url' field is left empty, source will be a text field.

Number	2.5.2.3
Requirement	User will be able to edit duration of a protocol
Note	Protocol duration is measured in days (integers only)

Number	2.5.2.4
Requirement	User will be able to tag a protocol
Note	

Number	2.5.2.5
Requirement	User will be able to sign protocol
Note	- Protocol owner and admins are authorized to sign the protocol - When another user is editing protocol page, the sign function will be disabled, and a message will appear on hover indicating the page is being edited

Number	2.5.2.6
Requirement	User will be able to duplicate protocol

Confidential

Note	
------	--

Number	2.5.2.7
Requirement	User will be able to start a new experiment from protocol
Note	• By clicking the 'Start Experiment' button • User will be able to name the experiment • User will be able to select project and folder • Protocol field is not editable • User will be notified when the experiment creation is completed with a link to the experiment show page

Number	2.5.2.8
Requirement	User will be able to link protocol to other pages
Note	• By clicking the linked resources button on the 'More' menu • By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	2.5.2.9
Requirement	Protocol owner will be able to set edit permissions for other members/teams
Note	By default, only admins and protocol owner have edit permissions

Number	2.5.2.10
Requirement	Protocol owner will be able to set view permissions for other members/teams
Note	By default, protocol is visible to all members

Number	2.5.2.11
Requirement	User will be able to view version history of a protocol
Note	

Number	2.5.2.12
Requirement	User will be able to print protocol to pdf file
Note	Protocol's name, owner, duration, creation date, tags, signed by (name and id) and signed date time will appear in pdf

Number	2.5.2.13
Requirement	User will be able to delete protocol
Note	

Confidential

Number	2.5.2.14
Requirement	User will be able to load sections from protocols to protocol
Note	By clicking on the 'add section' link and choosing protocol to load from

Number	2.5.2.15
Requirement	User will be able to add a new section to protocol
Note	By clicking on the 'add section' link

Number	2.5.2.16
Requirement	User will be able to download all attachments added to protocol
Note	

Number	2.5.2.17
Requirement	User will be able to add a comment to a protocol
Note	

Number	2.5.2.18
Requirement	User will be able to view: • Protocol creator • Creation date • Last update date
Note	

Number	2.5.2.19
Requirement	User will be able to export protocol as a Word document
Note	By clicking 'More' and then 'Export to Word'

Number	2.5.2.20
Requirement	User will be able to lock/unlock protocol format
Note	• User will see indication for locked protocol in the protocol • Experiments already created from the protocol will not be affected when protocol format is locked

Number	2.5.2.21
Requirement	User will be warned when starting a new experiment and some of the samples do not have stocks

Note	Warning can be removed by clicking 'Dismiss' or ignored by staring an experiment
------	--

Number	2.5.2.22
Requirement	User will be able to assign protocol to another member
Note	- Only admins and protocol owner - Assignee can be only a member who has viewing permissions to the protocol

2.5.3 Archiving protocols

Number	2.5.3.1
Requirement	User will be able to archive protocol
Note	By clicking on 'Archive' link via 'More' menu

Number	2.5.3.2
Requirement	User will not be able to edit archived protocol
Note	

Number	2.5.3.3
Requirement	User will not be able to start an experiment from an archived protocol
Note	

Number	2.5.3.4
Requirement	User will not be able to add archived protocols to existing experiment
Note	

Number	2.5.3.5
Requirement	User will not be able to delete archived protocol
Note	

Number	2.5.3.6
Requirement	User will be able to restore archived protocol
Note	By clicking on 'Restore' link via 'More' menu Only for the owner of the protocol and user with admin privileges

Number	2.5.3.7
Requirement	User will be able view list of archived protocols
Note	By clicking on 'View archived protocols' from Protocols index page

2.5.4 Actions in protocol sections

Number	2.5.4.1
Requirement	User will be able to add to protocol section the following elements: • Text • Spreadsheet • Steps • Samples • Plate • Equipment • Form • Sketch • Compound • Reaction • Date range • Attachments
Note	With the exception of date range elements that cannot be added to 'Description' sections

Number	2.5.4.2
Requirement	User will be able to delete an element from protocol section
Note	

Number	2.5.4.3
Requirement	User will be able to copy a protocol section to report
Note	

Number	2.5.4.4
Requirement	User will be able to add an attachment to protocol section
Note	

Number	2.5.4.5
Requirement	User will be able to move protocol section up/down
Note	

Number	2.5.4.6
--------	---------

Requirement	User will be able to delete a protocol section
Note	

Number	2.5.4.7
Requirement	User will be able to save a protocol section
Note	Saving a section will also close section edit mode

Number	2.5.4.8
Requirement	User will be able to move elements in protocol section up/down
Note	By dragging the element to the desired position

2.5.5 External protocols

Number	2.5.5.1
Requirement	User will be able to copy protocol from 'Protocol directory' to account protocols list
Note	By clicking 'Add from protocols directory'

Number	2.5.5.2
Requirement	User will be able to search protocols in directory by name
Note	

Number	2.5.5.3
Requirement	User will be able to share protocol to repository
Note	- Sharing must be first activated by account owner in account settings - Samples element cannot be copied to the repository

2.6 Data elements^(DTL)

2.6.1 General

Number	2.6.1.1
Requirement	User will be able to move data element and change their location in section.
Note	By clicking on the 'Move element' icon

Number	2.6.1.2
Requirement	User will be able to remove data element from section
Note	By clicking on the 'trash' icon

2.6.2 Compound element

Number	2.6.2.1
Requirement	User will be able to add new compound by: • Drawing • Uploading structure file
Note	• In Small molecule drawer pop up when adding compound element • When structure is duplicated, compound will be linked to the original compound in compound collection • Supported file types: ○ DL Molfile (v2000 and v3000) ○ RXN file format ○ InChI String ○ ChemAxon Extended SMILES ○ ChemAxon Extended CML file formats

Number	2.6.2.2
Requirement	User will be able to edit compound
Note	By clicking on compound drawing to open Small molecule drawer

2.6.3 Reaction element

Number	2.6.3.1
Requirement	User will be able to add new reaction by: • Drawing • Uploading file
Note	• In Small molecule drawer pop up when adding reaction element • When structure is duplicated, compound will be linked to the original compound in compound collection • Supported file types: ○ DL Molfile (v2000 and v3000) ○ RXN file format ○ InChI String ○ ChemAxon Extended SMILES ○ ChemAxon Extended CML file formats • When adding new reaction by drawing, user will be able to add: ○ Reaction plus tool ○ Arrow open angle tool

Number	2.6.3.2
Requirement	User will be able to edit a reaction

Note	By clicking on reaction drawing to open Small molecule drawer
------	---

Number	2.6.3.3
Requirement	User will be able to see stoichiometric table of a reaction
Note	Stoichiometric table includes: • Component • Role • Name • Formula • MW (gr/mol) • Density (gr/mol) • Mole • Mole ratio • Volume • Weight • %yield

Number	2.6.3.4
Requirement	User will be able to see the structure of a component in the stoichiometric table
Note	By hovering component column

Number	2.6.3.5
Requirement	User will be able to edit component's role in the stoichiometric table
Note	From a pre-defined list: • Reactant • Reagent • Solvent • Catalyst

Number	2.6.3.6
Requirement	User will be able to edit component's name in the stoichiometric table
Note	

Number	2.6.3.7
Requirement	User will be able to add the following values to a reaction component: • Density • Mole (mmol/mol) • Mole ratio • Volume (µl/ml/L) • Weight

Note	• If entering 'Density', 'Mole ratio' and 'Mole' of a component, 'Volume' and 'Weight' fields will be <u>automatically</u> calculate • If entering 'Density', 'Mole ratio' and 'Volume' of a component, 'Mole' and 'Weight' fields will be <u>automatically</u> calculate • If entering 'Density', 'Mole ratio' and 'Volume' of a component, 'Mole' and 'Volume' fields will be <u>automatically</u> calculate
------	--

Number	2.6.3.8
Requirement	User will be able to edit the following values to a reaction component: • Density • Mole (mmol/mol) • Mole ratio • Volume (µl/ml/L • Weight
Note	

Number	2.6.3.9
Requirement	User will not be able to edit the following values to a reaction component: • Formula • MW • %Yield
Note	

Number	2.6.3.10
Requirement	User will be able to calculate %yield for reaction product
Note	Mole ratio and volume/weight fields are required for calculation

Number	2.6.3.11
Requirement	User will be able to map identical atoms in reagents and products
Note	By clicking the Reaction Mapping Tools options: • Reaction Auto-Mapping Tool • Reaction Mapping Tool • Reaction Unmapping Tool

2.6.4 Text element

Number	2.6.4.1
Requirement	User will be able to add a text element to an experiment or protocol
Note	By clicking the ‘Add text’ button in section toolbar

Number	2.6.4.2
Requirement	User will be able to paste a copied image into text element
Note	

Number	2.6.4.3
Requirement	User will be able to restore text element to previous version
Note	

Number	2.6.4.4
Requirement	User will be able to insert elements directly to the text section: • Image • Video • Table • Unordered/ordered list • Horizontal line • Time stamp
Note	When the cursor is placed on the beginning of a new row, by clicking on ‘Quick insert +’ button on the left side of the section

2.6.5 Spreadsheet element

Number	2.6.5.1
Requirement	User will be able to add spreadsheet element to a • Projects • Experiment • Protocol • Recipe • Document • Reports • Request template • Request • Folder • SOP • Schedule page
Note	By clicking on the ‘Add spreadsheet’ button in section toolbar

Number	2.6.5.2
Requirement	User will be able to convert the spreadsheet to PDF
Note	By clicking on the 'PDF' link in the spreadsheet element

Number	2.6.5.3
Requirement	User will be able to export spreadsheet to CSV
Note	By clicking on the 'Export to CSV' link in the spreadsheet element

Number	2.6.5.4
Requirement	User will be able to adjust the spreadsheet size.
Note	By dragging the arrows below the spreadsheet up and down.

Number	2.6.5.5
Requirement	User will be able to protect spreadsheet with password
Note	• By clicking 'Format' menu via 'Home' tab • User can select which cells will be editable for all users

Number	2.6.5.6
Requirement	User will be able to open the spreadsheet on full screen
Note	• By clicking 'Format' menu via 'Home' tab • User can select which cells will be editable for all users

2.6.6 Steps element

Number	2.6.6.1
Requirement	User will be able to add steps element to • Projects • Experiment • Protocol • Recipe • Reports
Note	By clicking on the 'Add steps' link via the 'More' menu in any section.

Number	2.6.6.2
Requirement	User will be able to add step

Note	By clicking on the ‘+’ icon a new step template will be displayed
------	---

Number	2.6.6.3
Requirement	User will be able to set/stop timer for each step.
Note	By setting time and clicking on the ‘play’/‘stop’ button.

Number	2.6.6.4
Requirement	User will be able to mark step as completed
Note	By checking the ‘completed’ checkbox

Number	2.6.6.5
Requirement	User will be able to mark all steps as completed
Note	By checking the ‘All steps completed’ check box.

Number	2.6.6.6
Requirement	User will be able to reorder steps
Note	By clicking up & down on the arrows next to the step number.

Number	2.6.6.7
Requirement	User will be able to remove steps
Note	By clicking on the ‘x’ icon

Number	2.6.6.8
Requirement	User will be able to undo deletion of last step
Note	

Number	2.6.6.9
Requirement	User will be able to restore steps element to previous version
Note	

2.6.7 Plate element

2.6.7.1 ‘Add Plate’ popup

Number	2.6.7.1.1
Requirement	User will be able to add plate element to • Experiment • Protocol • Reports
Note	By clicking on the ‘Add plate’ link via the ‘More’ menu in any section.

	A plate customization popup will be opened. User will be required to choose the size of the plate. By clicking on the 'Apply changes' button plate will be added to the experiment/protocol/report section. The name of the plate will be contained from date & time.
--	---

Number	2.6.7.1.2
Requirement	After choosing plate size, user will be able to add samples to the plate's wells.
Note	By choosing well /multiple wells and checking samples from 'Samples' checkboxes list. User can add only samples added to sample element in same section.

Number	2.6.7.1.3
Requirement	User will be able to hide samples from 'Samples' list which are not used in the plate
Note	By clicking on the 'Hide x unused' Samples list.

Number	2.6.7.1.4
Requirement	User will be able to select/unselect all wells.
Note	By clicking on the 'Select all'/'clear selection' link

Number	2.6.7.1.5
Requirement	User will be able to set concentration in a single well
Note	By clicking on a well with samples a 'set concentration' section will be displayed.

Number	2.6.7.1.6
Requirement	User will be able to setup standards for more than one well, using the following parameters: • arrange replicates in rows/columns • set number of points • set number of replicates • set the starting quantity • set the serial factor
Note	By selecting more than one well a 'setup standards' section will be displayed.

Number	2.6.7.1.7
Requirement	User will be able to reset data in selected wells
Note	By clicking on the 'Reset data in selected wells' link. The link is displayed only when selecting wells with samples or text annotation

Number	2.6.7.1.8
Requirement	User will be able to resize plate when the plate is empty
Note	By clicking on the ‘Resize plate’ link

Number	2.6.7.1.9
Requirement	User will be able to add well annotations
Note	By choosing color and adding text.

Number	2.6.7.1.10
Requirement	User will be able to show/hide samples from plate
Note	By clicking on the eye icon in ‘Sample Annotations’ section.

Number	2.6.7.1.11
Requirement	User will be able to mark well as control
Note	By checking the ‘mark as control’ checkbox. well color will be changed to green once user marks it as a control.

Number	2.6.7.1.12
Requirement	User will be able to add description to sample annotation
Note	By default, description will be the sample name

Number	2.6.7.1.13
Requirement	User will be able to add symbol to sample annotation
Note	

Number	2.6.7.1.14
Requirement	User will be able to select color for sample annotation
Note	- Once color has been selected in the associated ‘Samples’ element, it will be set as default sample annotation color - Changing sample annotation color, does not affect the sample color selected in ‘Samples’ element

2.6.7.2 Plate element in section

Number	2.6.7.2.1
Requirement	User will be able to edit plate’s name
Note	By clicking on the ‘Edit’ pen next to plate title.

Number	2.6.7.2.2
Requirement	User will be able to print label of a plate
Note	By clicking on the ‘Print label’ link.

Number	2.6.7.2.3
Requirement	Requirement removed
Note	

Number	2.6.7.2.4
Requirement	User will be able view info of the wells content
Note	By hovering wells and clicking on the info tip of each well.

2.6.7.3 Import\export plate

Number	2.6.7.3.1
Requirement	User will be able to download an empty xlsx plate template
Note	

Number	2.6.7.3.2
Requirement	User will be able to export plate to xlsx template
Note	

Number	2.6.7.3.3
Requirement	User will be able to import plate from an xlsx template
Note	

Number	2.6.7.3.4
Requirement	User will be notified when trying to import plate from an invalid file
Note	

Number	2.6.7.3.5
Requirement	User will be required to add valid values under 'Row' and 'Column' columns
Note	- By adding letters to 'Rows' - By adding numbers to 'Columns'

Number	2.6.7.3.6
Requirement	User will be able to import plate with only annotations or plate properties
Note	Values to columns 'Row' and 'Columns' still must be added

Number	2.6.7.3.7
Requirement	User will be able to set plate size in xlsx template
Note	Default size is 96-wells

Number	2.6.7.3.8
Requirement	User will be able to mark wells as control in xlsx template
Note	

Number	2.6.7.3.9
Requirement	User will be able to add well annotation in xlsx template
Note	- Limited to 255 characters - Limited to one annotation per well

Number	2.6.7.3.10
Requirement	User will be notified when trying to import a plate with well annotation with more than 255 characters
Note	

Number	2.6.7.3.11
Requirement	User will be able to add existing inventory items to wells
Note	- By entering name or SysID of item that exists in collection - In case there are multiple items with the same name, the first will be added to plate

Number	2.6.7.3.12
Requirement	User will be able to create new inventory items
Note	By adding item name that does not exist in collection

Number	2.6.7.3.13
Requirement	User will be able to add existing stocks to wells
Note	- By entering name or id of existing stock - In case there are multiple stocks with the same name, the first will be added to plate

Number	2.6.7.3.14
Requirement	User will be able to create new stocks
Note	

Number	2.6.7.3.15
Requirement	User will be able to add concentration to wells
Note	

Number	2.6.7.3.16
Requirement	User will be able to add/edit plate properties
Note	

Number	2.6.7.3.17
Requirement	User will be notified when: • trying to import a plate with missing values under 'Rows' or 'Column' • trying to import a plate with invalid values under 'Rows' or 'Column' • trying to import a plate with values under 'Rows' or 'Column' that do not fit the plate size • trying to import a plate with non-numeric values under 'Concentration'
Note	

2.6.8 Labguru Form element

2.6.8.1 Form element in protocols

Number	2.6.8.1.1
Requirement	User will be able to add: • text • images • links • videos • tables • special characters • horizontal lines to form element
Note	

Number	2.6.8.1.2
Requirement	User will be able to add the following form input fields: • Text • Numbers only • Paragraph • Collection • Checkbox • Radio buttons • Select box • Date picker • Time picker • Members • Equipment
Note	

Number	2.6.8.1.3
Requirement	User will be able to add a hint to: • Text • Numbers only • paragraph form input fields

Confidential

Note	
------	--

Number	2.6.8.1.4
Requirement	User will be able to name input fields
Note	- Required when connecting form to dataset - Except for collection input field

Number	2.6.8.1.5
Requirement	User will be able to mark input fields as required
Note	Except Checkbox and Radio button

Number	2.6.8.1.6
Requirement	User will be able to group multiple checkboxes or radio buttons by name
Note	Required when connecting form to dataset

Number	2.6.8.1.7
Requirement	User will be able to label checkboxes and radio buttons by name
Note	Label will be visible next to the checkboxes/radio button after input field was added

Number	2.6.8.1.8
Requirement	User will be able to set form behavior when out-of-range values will be added to the form in experiment page: - Allow and highlight out-of-range values in red - Allow out-of-range value only if user entered reason of deviation - Don't allow out-of-range values
Note	

Number	2.6.8.1.9
Requirement	User will not be able to add values to form input fields in protocol
Note	

Number	2.6.8.1.10
Requirement	User will be able to allow row duplication in tables in form, for experiments started from protocol
Note	By checking 'Allow row duplication' in form settings

Number	2.6.8.1.11
Requirement	User will be able to set character number limit to text input fields
Note	

Number	2.6.8.1.12
Requirement	User will be able to give a name to form
Note	

2.6.8.2 Form element in experiments

Number	2.6.8.2.1
Requirement	User will not be able to edit form element expect interaction with form input fields
Note	

Number	2.6.8.2.2
Requirement	User will not be able to add text to numbers only input fields
Note	

Number	2.6.8.2.3
Requirement	User will be warned when entering an out of range value in a numbers-only object with min/max values
Note	The object will be highlighted in red and when hovering there will be a bubble with explanation of the error

Number	2.6.8.2.4
Requirement	User will be able to check multiple checkboxes from the same group
Note	

Number	2.6.8.2.5
Requirement	If form has more than one group of radio buttons, user will be able to choose one from each group
Note	

Number	2.6.8.2.6
Requirement	User will be able to select date in date picker input field
Note	

Number	2.6.8.2.7
Requirement	User will be able to enter time in time picker input field
Note	Time format is hh-mm-AM/PM

Number	2.6.8.2.8
Requirement	User will be able to clear data in date/time picker input fields

Confidential

Note	
------	--

Number	2.6.8.2.9
Requirement	User will be able to delete form element
Note	

Number	2.6.8.2.10
Requirement	User will be able to duplicate the last row in table form
Note	Option needs to be enabled in the protocol in advance

Number	2.6.8.2.11
Requirement	User will be able to remove a row from table
Note	- Only when row duplication was enabled in protocol - Double confirmation will be required

Number	2.6.8.2.12
Requirement	User will not be able to sign experiment if there is an empty required input filed
Note	

Number	2.6.8.2.13
Requirement	User will not be able to select archived equipment in equipment input field
Note	

Number	2.6.8.2.14
Requirement	User will be required to select reason for change from a pre-defined list when editing form input field
Note	- Only if enabled in account settings - In the 'Select reason for change' pop up - When any value is changed, reason for change pop up will be displayed only after the first entry is saved

Number	2.6.8.2.15
Requirement	User will be able to see change description/full location of the changed value (if in table - column, row)
Note	In the 'Select reason for change' pop up

Number	2.6.8.2.16
Requirement	User will be able to enter free text in the 'Select reason for change' pop up
Note	Only if enabled in account settings

Number	2.6.8.2.17
Requirement	User will be able to add value to a target form input field based on the field's name, and on the experiment ID, using the Measurements API request
Note	- In request body: - Value is 'measurement' - Input field name is 'input_name' - Experiment ID is 'experiment_id' - After the data is populated into the target field, the field should be highlighted to indicate a change was made and the user needs to save a section - If experiment has multiple fields with the same name, only the first (that is in the form that was created first) field will be updated

Number	2.6.8.2.18
Requirement	User will not be able to edit collection auto-populated field
Note	

2.6.8.3 Forms to datasets

Number	2.6.8.3.1
Requirement	User will be able to connect form element to a newly created dataset
Note	

Number	2.6.8.3.2
Requirement	User will not be able to create a dataset if form element do not have any objects with names or any other "mandatory for dataset" fields
Note	

Number	2.6.8.3.3
Requirement	For each form object with name, a column will be added to the dataset with the value in the name field as the column header
Note	

Number	2.6.8.3.4
Requirement	Every experiment started from a protocol with form element that is connected to a serial dataset will create a new row in the dataset
Note	

Confidential

Number	2.6.8.3.5
Requirement	When user edit a form element that is connected to a dataset, the dataset will be updated accordingly
Note	

Number	2.6.8.3.6
Requirement	User will be able to click on a link to get to the connected dataset
Note	

Number	2.6.8.3.7
Requirement	User will be able to go to the experiment page of which the data came from for each row in the dataset
Note	From datasets show page

Number	2.6.8.3.8
Requirement	User will be able to go to the protocol connected to the dataset
Note	From datasets show page

Number	2.6.8.3.9
Requirement	When user delete rows in form element linked to a dataset, the dataset will be updated accordingly
Note	

Number	2.6.8.3.10
Requirement	User will be able to select dataset type: - Table dataset - Serial dataset
Note	- On 'Form settings' pop up - Only when forms contain an inserted table - In table dataset type, user can select from 2 options: - Create a new dataset for each new experiment initiated from the protocol - Consolidate data into a single dataset across all experiments initiated from the same protocol

Number	2.6.8.3.11
Requirement	User will be required to name the dataset
Note	

Number	2.6.8.3.12
Requirement	User will not be able to duplicate row in forms with table linked to a serial dataset type
Note	

Number	2.6.8.3.13
Requirement	Every experiment started from a protocol with form element that is connected to a table dataset will create a new dataset
Note	• Link to the new dataset will be created • When duplicating an experiment/ copying a section containing a form linked to a table dataset, a new dataset will be created

Number	2.6.8.3.14
Requirement	Every experiment started from a protocol with form element that is connected to an aggregational table dataset will create a new row in the dataset for each row created in table within the for element
Note	When duplicating an experiment/ copying a section with a form element linked to an aggregated table dataset, a link to the original aggregative dataset will be presented

Number	2.6.8.3.15
Requirement	User will be able to remove link to dataset
Note	• By clicking 'Unlink' • In newly created dataset with table inserted

2.6.8.4 Forms to datasets reason TBD

2.6.9 Samples and reagents element

2.6.9.1 General

Number	2.6.9.1.1
Requirement	User will be able to add inventory items to Samples and reagents element
Note	• By choosing collection and then typing the name of item or any other field from the collection • User will not be able to add archived items to element

Number	2.6.9.1.2
Requirement	User will be able to create a new inventory item from Samples and reagents element
Note	• By selecting collection and entering the item name • User will be required to add unique name when unique name field is enabled in the collection

Confidential

	User will be notified when trying to create an item with an existing name
--	---

Number	2.6.9.1.3
Requirement	User will be required to add required fields when creating a new item in the Samples and reagents element (if exist in collection)
Note	Pop up will open to fill the required field

Number	2.6.9.1.4
Requirement	User will be able to create a new inventory item when he searched for an item that do not exists in the inventory
Note	

Number	2.6.9.1.5
Requirement	User will be able to go to a item's show page in the inventory by clicking on its name When item is selected to a sample element
Note	

Number	2.6.9.1.6
Requirement	User will see samples grouped by collection, each collection with its own default columns
Note	

Number	2.6.9.1.7
Requirement	User will be able to remove one sample
Note	

Number	2.6.9.1.8
Requirement	User will be able to remove all samples of one collection
Note	

Number	2.6.9.1.9
Requirement	User will be able to select stock/s to add for a sample
Note	- Only on experiment page - User will not be able to add stocks to sample of an archived item

Confidential

Number	2.6.9.1.10
Requirement	User will be able to go to the stocks show page from the 'Select stocks' pop-up
Note	Only on experiment page

Number	2.6.9.1.11
Requirement	User will be able to go to the box show page from the 'Select stocks' pop-up
Note	Only on experiment page

Number	2.6.9.1.12
Requirement	User will be able to mark stock as consumed
Note	Only on experiment page

Number	2.6.9.1.13
Requirement	User will be able to add a remark to stock
Note	Only on experiment page

Number	2.6.9.1.14
Requirement	User will see time stamp for when Samples and reagents element was added
Note	Only on experiment page

Number	2.6.9.1.15
Requirement	User will be able to remove a stock from Samples and reagents element
Note	Only on experiment page

Number	2.6.9.1.16
Requirement	User will be able to select color for sample
Note	Selected color will be set as default for plate elements added to the same section

Number	2.6.9.1.17
Requirement	User will be able to add a stock to Samples and reagents element by scanning stock's barcode or QR code
Note	Stock with matching values in Id or barcode fields will be added to element

Number	2.6.9.1.18
Requirement	User will be able to select which columns to present for each collection
Note	- Limited for max 7 column - User preference is saved is all samples element for that collection - User will be able to add parent data to sample element

Number	2.6.9.1.19
Requirement	User will be able to re-order columns' order for each collection
Note	

Number	2.6.9.1.20
Requirement	User will be able to create new stocks directly from the sample element
Note	

2.6.9.2 Protocols and samples

Number	2.6.9.2.1
Requirement	User will be able to lock Samples and reagents element
Note	In experiments, user will be able to only add/remove or update stocks for samples added at protocol

Number	2.6.9.2.2
Requirement	When item was added to Samples and reagents element, user will see in inventory item show page, under 'Experiments' or 'Protocols' tabs
Note	

2.6.9.3 Search results in sample element

Number	2.6.9.3.1
Requirement	User will be able to remove/add columns from search results table
Note	

Number	2.6.9.3.2
Requirement	User will be able to change the order of the columns
Note	

Number	2.6.9.3.3
Requirement	User will be able to sort columns in ascending or descending order
Note	

Number	2.6.9.3.4
Requirement	User will be able add all samples in view in a single button
Note	

Number	2.6.9.3.5
Requirement	User will be able to see the same preference setting the next time he will search for items in that collection
Note	

2.6.9.4 Update stock directly from sample element

Number	2.6.9.4.1
Requirement	User will be able to mark the option 'Automatically mark stocks as consumed', in the account settings
Note	

Number	2.6.9.4.2
Requirement	User will be able to enter the amount used for a selected stock
Note	By updating 'Volume used' and 'Weight used' columns

Number	2.6.9.4.3
Requirement	User will be able to specify units used in 'Volume used' and 'Weight used' columns
Note	

Number	2.6.9.4.4
Requirement	User will be able to select from pre-defined list of volume/weight units
Note	

Number	2.6.9.4.5
Requirement	User will be able to select other types of volume/weight units if they were added to the account in advance
Note	

Confidential

Number	2.6.9.4.6
Requirement	Requirement removed
Note	

Number	2.6.9.4.7
Requirement	User will not be able to update the volume/weight left in a stock from the sample element without entering volume/weight used first
Note	

Number	2.6.9.4.8
Requirement	User will be able to revert the volume/weight used after he clicked on "subtract", as long as page was not refreshed
Note	By clicking on 'Revert' button

Number	2.6.9.4.9
Requirement	User will see indication for volume/weight used on the stock page
Note	In 'History' tab, together with experiment name in which the volume/weight was reduced

Number	2.6.9.4.10
Requirement	User will not be able to bring stock volume/weight to below zero
Note	Will get an error message when the final stock volume/weight is supposed to be below zero

Number	2.6.9.4.11
Requirement	Requirement removed
Note	

Number	2.6.9.4.12
Requirement	Requirement removed
Note	

Number	2.6.9.4.13
Requirement	User will not be able to remove stock that was used (volume or weight) from sample element
Note	

2.6.9.5 Output Stocks2.6.9.5.1 Mark stock as output

Number	2.6.9.5.1.1
Requirement	User will be able to mark stock as output
Note	

Number	2.6.9.5.1.2
Requirement	User will not be able to remove stocks that were marked as output stock from sample element
Note	

Number	2.6.9.5.1.3
Requirement	User will not be able to mark a consumed stock as an output stock
Note	

Number	2.6.9.5.1.4
Requirement	User will not be able to mark a stock as an output stock, if it was already defined as an output stock on a different experiment
Note	

Number	2.6.9.5.1.5
Requirement	Requirement removed
Note	

Number	2.6.9.5.1.6
Requirement	User will be able to unmark output stock
Note	

2.6.9.6 Sample properties

Number	2.6.9.6.1
Requirement	User will be able to add samples properties to sample element
Note	• By clicking on 'Properties' button • Only in protocols • Samples properties added in protocol will be visible in experiments that will start from that protocol

Number	2.6.9.6.2
Requirement	User will be able to insert a value into a sample property
Note	In experiments and protocols

Number	2.6.9.6.3
Requirement	User will be required to select name for samples property
Note	- If name field is left empty, property will not be created - Default properties names will be "Property X"

Number	2.6.9.6.4
Requirement	User will be able to select type for each samples property
Note	From text, numbers or date

Number	2.6.9.6.5
Requirement	User will be able to add more than one property to one 'Samples element'
Note	- By default, one property will be available - User will be able to add up to 3 properties

Number	2.6.9.6.6
Requirement	User will be able to delete property
Note	By clicking on the bin icon next to the property

Number	2.6.9.6.7
Requirement	User will be able to edit property name and type
Note	By clicking on the field and change the value

2.6.9.7 Check for updates

Number	2.6.9.7.1
Requirement	User will be able to check if changes were made to items and stocks that were added to the samples element
Note	

Number	2.6.9.7.2
Requirement	User will be able to accept or reject a change
Note	

Number	2.6.9.7.3
Requirement	User will be able to accept or reject all changes
Note	

2.6.10 Attachments element

Number	2.6.10.1
Requirement	User will be able to add attachment by uploading files to Attachments element
Note	• By clicking on the attachments box and browsing for the files • By dragging and dropping the files onto the attachments box

Number	2.6.10.2
Requirement	User will be able to add unattached attachment to the element
Note	By clicking on 'Add unattached files'

Number	2.6.10.3
Requirement	User will be able to download all attachments in an Attachments element
Note	• By clicking on 'Download all attachments' • Download will be as zip file

Number	2.6.10.4
Requirement	User will be able to add to page images and XLSX, PDF, PPTX files to page
Note	By clicking on arrow down button next to attachment and then 'Add to page'

Number	2.6.10.5
Requirement	User will see in the added element to which attachment it is synced with
Note	

Number	2.6.10.6
Requirement	User will be able to edit the attachment by editing the added element
Note	Images and Excel attachments only

Number	2.6.10.7
Requirement	User will be able to open attachment using Local Edit
Note	See section 6.9 in FRS

Number	2.6.10.8
Requirement	User will be able to download attachment
Note	By clicking on arrow down button next to attachment and then 'Download'

Number	2.6.10.9
Requirement	User will be able to delete attachment
Note	• By clicking on arrow down button next to attachment and then 'Delete' • User will not be able to delete an attachment that was added to page • If attachment was added to existing dataset, the associated vectors will be removed from dataset

Number	2.6.10.10
Requirement	User will be able to save XLSX, TXT and CSV files as a dataset
Note	• By clicking on arrow down button next to attachment and then 'Add as dataset' • The dataset will be created only for file with a valid structure for dataset • Only in experiments • User will be able to select to save as new dataset or add to existing • Only datasets with matching headers will be available to select from

Number	2.6.10.11
Requirement	User will be able to rename element's title
Note	

2.6.11 Code element

Number	2.6.11.1
Requirement	User will be able to type code into Code element
Note	Syntax will be shown based on the language selected

Number	2.6.11.2
--------	----------

Requirement	User will be able to name the Code element
Note	By clicking on the pencil icon next to the name

Number	2.6.11.3
Requirement	User will be able to select the code language
Note	Selection from: GO, Javascript, Python, R, Ruby, SQL

Number	2.6.11.4
Requirement	User will be able to select Code element theme
Note	Selection from: Monokai, Elegant, XQ Dark, XQ light, Solarzied

Number	2.6.11.5
Requirement	User will be able to restore code element to previous version
Note	

2.6.12 Equipment element

2.6.12.1 General

Number	2.6.12.1.1
Requirement	User will be able to add an equipment
Note	- By typing the name of equipment or scanning their barcode - User will not be able to add archived equipment to element

Number	2.6.12.1.2
Requirement	User will be able to create a new equipment
Note	- By clicking 'Create new equipment' and entering it's name in the search field - When user search for equipment that do not exists in the inventory, he will be able to create a new one - User will be required to add unique name

Number	2.6.12.1.3
Requirement	User will be able to go to equipment's show page in the by clicking on its name
Note	

Number	2.6.12.1.4
Requirement	User will be able to remove an equipment/all element from the list

Confidential

Note	User will not be able to delete equipment if it was added to experiment/protocol
------	--

Number	2.6.12.1.5
Requirement	User will be able to view equipment's calendar
Note	

Number	2.6.12.1.6
Requirement	User will be able to add a remark to instrument
Note	

Number	2.6.12.1.7
Requirement	User will be able to select which columns to present for each instrument
Note	Limited for max 7 column

Number	2.6.12.1.8
Requirement	User will be able to re-order columns' order for each instrument
Note	

Number	2.6.12.1.9
Requirement	User will be able to record equipment usage during an experiment
Note	• By clicking 'Add usage record' • Multiple usage lines can be added

Number	2.6.12.1.10
Requirement	User will be able to edit equipment title
Note	

2.6.12.2 Usage Record

Number	2.6.12.2.1
Requirement	User will be able to edit operator name
Note	Automatically populated with the user who enters the record

Number	2.6.12.2.2
Requirement	User will be able to edit operating date and time

Note	Automatically populated with the date and time when the record is created
------	---

Number	2.6.12.2.3
Requirement	User will be able to add operating time and units
Note	- Positive numeric value - Unit options from dropdown list: seconds, hours, minutes, days

Number	2.6.12.2.4
Requirement	User will be able to add usage details
Note	

Number	2.6.12.2.5
Requirement	User will be able to delete an equipment usage record
Note	Warning message will be displayed for confirmation

2.6.12.3 Search results in equipment element

Number	2.6.12.3.1
Requirement	User will be able to remove/add columns from search results table
Note	

Number	2.6.12.3.2
Requirement	User will be able to change the order of the columns
Note	

Number	2.6.12.3.3
Requirement	User will be able to sort columns in ascending or descending order
Note	

Number	2.6.12.3.4
Requirement	User will be able add all equipment's in view in a single button
Note	

2.6.12.4 Check for updates

Number	2.6.12.4.1
Requirement	User will be able to check if changes were made to equipment that was added
Note	

Number	2.6.12.4.2
Requirement	User will be able to accept or reject a change
Note	

Number	2.6.12.4.3
Requirement	User will be able to accept or reject all changes
Note	

2.6.13 Date Range element2.6.13.1 General

Number	2.6.13.1.1
Requirement	User will be able to add date range element
Note	• In Experiments: User can select between relative or absolute date range • In Protocols: User can add only relative date range • User will be able to add only one date range element to a section • User will see indication of start and end times in the top of the section • Date range event will be created

Number	2.6.13.1.2
Requirement	User will be able to edit date range element
Note	• In Experiments: by clicking: 'Edit Date Range' in the 'More' menu • In Protocols: by clicking: 'Edit relative date Range' in the 'More' menu • Date range event will be updated

2.6.13.2 Relative

Number	2.6.13.2.1
Requirement	User will be able to select start point for relative date range elements
Note	• Selection between start date of experiment or previous section end date

	Default value is start date of experiment
--	---

Number	2.6.13.2.2
Requirement	User will be able to set how long after selected start point the date range should start (months/days/hours)
Note	Experiment/protocol duration will be updated according to time set

Number	2.6.13.2.3
Requirement	User will be required to set duration for date range
Note	Experiment/protocol duration will be updated according to time set

Number	2.6.13.2.4
Requirement	User will be able to remove date range
Note	By changing all elements to zero, when editing an existing date range element

2.6.13.3 Absolute

Number	2.6.13.3.1
Requirement	User will be required to select start and end date and time
Note	

Number	2.6.13.3.2
Requirement	User will be able to choose if date range duration will affect experiment duration
Note	

Number	2.6.13.3.3
Requirement	User will be able to remove date range
Note	By clicking on 'X' button in the Start date or End fate box.

2.6.14 Sketch element

TBA

2.6.15 Restore element to previous version

Number	2.6.15.1
Requirement	User will be able to select a version to view from list of all versions
Note	By default, current version will be displayed

Number	2.6.15.2
Requirement	User will be able to filter list of version by creating member
Note	

Number	2.6.15.3
Requirement	User will be able to filter list of version by creation date
Note	

Number	2.6.15.4
Requirement	User will be able to jump to an older or a newer version
Note	

Number	2.6.15.5
Requirement	User will be able to see number of versions available
Note	

2.7 Labguru Samples basket module(SMB)

Number	2.7.1
Requirement	User will be able to open the samples basket
Note	By clicking the samples basket icon on the purple toolbar

Number	2.7.2
Requirement	User will be able to close the samples basket
Note	- By clicking the samples basket icon on the purple toolbar - By clicking the x button on the top-right corner of the sample basket

Number	2.7.3
Requirement	User will be able to add an inventory item to the samples basket by searching for it in selection menu
Note	- User must first select sample collection from collections dropdown - User can then search for an inventory item by typing the first 3 letters of the item's name in the search bar - When typing the item's name, a table with matching results appears (only items from the selected collection will be displayed) - By choosing collection and then typing the name of item or any other field from the collection - User will not be able to add archived items to basket

Number	2.7.4
Requirement	User will be able to add an inventory item to the samples basket by scanning a barcode
Note	Stock with matching values in stock id or barcode fields will be added to the basket

Number	2.7.5
Requirement	User will be able to remove/add columns from search results table when searching for an item in samples basket
Note	

Number	2.7.6
Requirement	User will be able to change the order of the columns when searching for an item in samples basket
Note	

Number	2.7.7
Requirement	User will be able to sort columns in ascending or descending order
Note	

Number	2.7.8
Requirement	User will be able add all samples in view in a single button
Note	

Number	2.7.9
Requirement	User will be able to see the same preference setting the next time he will search for items in that collection
Note	

Number	2.7.10
Requirement	User will be able to add stocks to the sample basket
Note	- By clicking the 'Select stocks' button and choosing the desired stocks - Selected stocks will be automatically added to the sample basket - User will not be able to add stocks to item that was archived after added to basket

Number	2.7.11
Requirement	User will be able to navigate to inventory items that were added to the sample basket
Note	Clicking on the name of an item will direct User to the item show page

Confidential

Number	2.7.12
Requirement	User will be able to navigate to stocks that were added to the sample basket
Note	Clicking on the name of a stock will direct User to the stock show page

Number	2.7.13
Requirement	User will be able to create new stocks directly from the sample basket
Note	By clicking the 'Create new stock/s' button on the 'Select stock' popup

Number	2.7.14
Requirement	User will be able to add remarks to a sample
Note	By adding text to the remarks section next to each sample

Number	2.7.15
Requirement	Requirement removed
Note	

Number	2.7.16
Requirement	User will be able to remove a sample
Note	By clicking the trash icon next to the sample

Number	2.7.17
Requirement	User will be able to remove a stock
Note	By clicking the trash icon next to the stock

Number	2.7.18
Requirement	User will be able to start a new experiment that includes the inventory items and stocks that were added to the sample basket
Note	- By clicking the 'New experiment' button on the bottom of the sample basket - Clicking the 'New experiment' button opens a popup where User can choose the project, folder, and the experiment name - The new experiment will include a samples element with all the items and stocks that were added to the sample basket - Consumed stocks will not be added to new experiments

Confidential

Number	2.7.19
Requirement	User will be able to scan storage locations into samples basket
Note	- By scanning storage location custom barcode, ID or selecting manually from list - By entering its auto-generated barcode with prefix STR+ ID or a custom barcode value into the search bar

Number	2.7.20
Requirement	User will be able to scan multiple locations
Note	Scanned storage locations will be displayed below scanned stocks

Number	2.7.21
Requirement	User will be able to remove a storage location from the samples basket
Note	by clicking on 'X'

Number	2.7.22
Requirement	User will be able to add storage location by manual searching
Note	when storage location selected already (either manually or by scanning) the option for manual search will not be visible

Number	2.7.23
Requirement	User will be able to move stocks to a new target storage location
Note	By clicking on 'Move stocks'

Number	2.7.24
Requirement	User will be able to clear basket content
Note	

Number	2.7.25
Requirement	User will be able to create new inventory item from the sample basket
Note	- By clicking 'Create new sample' button or via search field for non-existing item - User will be required to add unique name when unique name field is enabled in the collection

Number	2.7.26
Requirement	User will be required to add mandatory fields when creating a new item in the Sample basket (if exist in collection)
Note	Pop up will open to fill the required field

Number	2.7.27
Requirement	User will be able to add basket to existing experiment
Note	- Only users with editing permissions to experiment 'Add to current experiment' will be available only when the user is in an experiment page

Number	2.7.28
Requirement	User will be able to select which experiment section to add basket to
Note	In the 'Select Section' pop up

Number	2.7.29
Requirement	User will be able to create new experiment section
Note	Basket will be added in section called 'Custom' at the bottom of the experiment

Number	2.7.30
Requirement	User will not be able to add basket to signed experiment
Note	

Number	2.7.31
Requirement	User will not be able to add basket to experiment from archived project
Note	

Number	2.7.32
Requirement	User will not be able to add basket to experiment from protocol with locked format
Note	

Number	2.7.33
Requirement	User will be able to pool stocks
Note	More than one stock must be selected to pool a new stock

2.8 Folders module (FLD)

2.8.1 Folder show page

Number	2.8.1.1
Requirement	User will be able to add a new experiment
Note	By clicking the 'New Experiment' button

Confidential

	• User will be able to name the experiment • User will be able to select project, folder and protocol from list/search by typing the first letter • When protocol is not selected, the experiment will be empty
--	---

Number	2.8.1.2
Requirement	User will be able to tag folder
Note	

Number	2.8.1.3
Requirement	User will be able to print folder to PDF file
Note	Printout will include list of experiment that are directly in the printed folder

Number	2.8.1.4
Requirement	User will be able to set default edit permissions for experiments to selected members/teams
Note	By clicking on 'Permissions' in the 'More' menu • By default, only admins, project owner and experiment owner have edit permissions • User will be able to select 'All' or 'None' from the offered members list • Only admins, project owner and folder owner can edit experiment default edit permissions at folder level

Number	2.8.1.5
Requirement	User will be able to move folder to another project/folder
Note	• By selecting project/folder from the list • Only admins, project owner and folder owner can move folder

Number	2.8.1.6
Requirement	User will be able to delete folder
Note	• By clicking on 'Delete folder' in the 'More' menu • Only admins, project owner and folder owner can delete folder • Folder cannot be deleted if it contains a signed experiment

Number	2.8.1.7
Requirement	User will be able to add comment to folder
Note	

Number	2.8.1.8
Requirement	User will be able to link folder to other pages
Note	• By clicking the linked resources button on the 'More' menu • By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	2.8.1.9
Requirement	User will be able to add a task connected to folder
Note	Task assignee can be only a member who has view permissions to the project the folder is part of

Number	2.8.1.10
Requirement	User will be able to edit folder name
Note	

Number	2.8.1.11
Requirement	User will be able to assign another member to be the folder owner
Note	Can assign only member that have access to project the folder is part of

Number	2.8.1.12
Requirement	User will be able to set folder due date
Note	User can set who receives notification on due date and when notification should be sent

Number	2.8.1.13
Requirement	User will be able to see folder start date, duration and progress
Note	• Start date is calculated by the earliest experiment start date • Duration is calculated by time from start date of first experiment to end date of last experiment • Folder progress is calculated by number of signed experiments from total number of experiments • For calculation of the end date refer 2.4.1.16

Number	2.8.1.14
Requirement	User will be able to add folder description
Note	Available elements: Text box, spreadsheet, steps, sketch, attachments and code

Number	2.8.1.15
Requirement	User will be able to expand/collapse folders
Note	• User preferences will be saved • New folder will collapse by default

2.8.2 Experiments tab

Number	2.8.2.1
Requirement	User will be able to view list of experiments in the folder
Note	User will be able to enter an experiment page by clicking on its name or ID

Number	2.8.2.2
Requirement	User will be able to change the number of items displayed per page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 experiments per page

Number	2.8.2.3
Requirement	User will be able to sort experiments by any column in display in ascending/descending order
Note	Filter preferences are saved

Number	2.8.2.4
Requirement	User will be able to filter experiments by any combination of columns in display
Note	Filter preferences are saved

Number	2.8.2.5
Requirement	User will be able to choose which columns to view on the index page
Note	• By choosing headers from list in 'Columns' option • User's preference will be saved • possible headers: Created at, Last updated, ID, Experiment name, Project, Owner, start date, End date, Status, Signed, Witnessed

Number	2.8.2.6
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	2.8.2.7
Requirement	User will be able to change the order of the columns
Note	

Number	2.8.2.8
Requirement	User will be able to enter an experiment page by clicking on its name or ID
Note	

2.8.3 Folders tab

Number	2.8.3.1
Requirement	User will be able to see list of subfolders and experiments that are under the folder
Note	

Number	2.8.3.2
Requirement	User will be able to add new sub-folder
Note	By clicking 'Add new folder'

Number	2.8.3.3
Requirement	User will be able to enter folder page
Note	By clicking folder name

Number	2.8.3.4
Requirement	User will be able to enter experiment show page
Note	By clicking experiment name

Number	2.8.3.5
Requirement	User will be able to re-arrange subfolders
Note	- By dragging and dropping the folder - User preferences will be saved

2.8.4 Figures tab

Number	2.8.4.1
Requirement	User will be able to add figures
Note	- By clicking 'Add new figure' - User will be able to add figure title/ description - User will be able to select existing figures or upload a new one

3. Knowledgebase

3.1 Labguru Datasets module (DTS)

3.1.1 'Datasets' index page

Number	3.1.1.1
Requirement	User will be able to add a new dataset
Note	- By clicking the 'New Dataset' button on datasets index page - User will be required to add dataset name and upload a data file (xlsx or csv) - User will be notified when the dataset creation was completed with a link to the dataset show page

Number	3.1.1.2
Requirement	User will be able to search datasets by name
Note	By using the search bar on the index page of the datasets

Number	3.1.1.3
Requirement	User will be able to filter datasets by tags
Note	

Number	3.1.1.4
Requirement	User will be able to change the number of datasets he sees on single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 documents per page

Number	3.1.1.5
Requirement	User will be able to sort datasets by any column in display in ascending/descending order
Note	

Number	3.1.1.6
Requirement	User will be able to filter datasets by any combination of columns in display
Note	

Number	3.1.1.7
Requirement	User will be able to choose which columns to view on the index page
Note	• By choosing headers from list in 'Columns' option • User's preference will be saved

Number	3.1.1.8
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	3.1.1.9
Requirement	User filter/sort selection will be saved when user move to other pages in the collection
Note	

Number	3.1.1.10
Requirement	User will be able to tag/bulk tag datasets
Note	

Number	3.1.1.11
Requirement	User will be able to delete/bulk delete datasets
Note	Double confirmation is required

3.1.2 Dataset show page

Number	3.1.2.1
Requirement	User will be able to access all data points dataset
Note	• By clicking on the datapoint • Each data point is linked to the datapoint page

Number	3.1.2.2
Requirement	User will be able to edit dataset name and description
Note	• By clicking the 'Edit' button • Clicking the 'Edit' button directs User to the 'Edit dataset' page where he can edit the name and description fields

Number	3.1.2.3
Requirement	User will be able to delete dataset

Confidential

Note	- By clicking the 'Delete' button - When clicking the 'Delete' button a pop-up appears with a warning 'Please note that all stocks, attachments and data associated with this item will also be deleted'
------	---

Number	3.1.2.4
Requirement	User will be able to export the dataset to csv
Note	- By clicking the 'Export' button - Clicking the 'Export' button automatically downloads a csv file that includes the dataset

Number	3.1.2.5
Requirement	If dataset was created from a Labguru page, user will be able to see a link to the page
Note	- A page in which a data file was attached and then a dataset was created accordingly - The page will appear as a link on the 'Data' tab

Number	3.1.2.6
Requirement	User will be able to add a comment to dataset
Note	

Number	3.1.2.7
Requirement	User will be able to tag dataset
Note	By clicking the 'Add tag' button on the dataset show page

Number	3.1.2.8
Requirement	User will be able to link dataset to other Labguru pages
Note	By clicking the 'Add link' button on the dataset show page

Number	3.1.2.9
Requirement	User will be able to see inventory items related data points on inventory items show pages
Note	On the 'Datasets' tab

Number	3.1.2.10
Requirement	Requirement removed
Note	

Number	3.1.2.11
Requirement	User will be able to create a dataset from a plate show page
Note	By clicking the 'Add new dataset' button

Number	3.1.2.12
Requirement	Requirement removed
Note	

Number	3.1.2.13
Requirement	User will be able to limit dataset visibility to selected members or teams
Note	• By clicking on 'Visible to...' • Only available to dataset owner and admins • User without view permission will not see the dataset in the datasets index and will not be able to view it's page • User with role with "can administer" privilege will be able to view dataset irrespective of the permissions assigned

Number	3.1.2.14
Requirement	User will be able to add Excel file to dataset
Note	• By clicking on 'Add to existing' • The vectors associated with the file will show the file's name in the 'Origin' column • The file will be added to the 'Attachments' box in dataset page • Only files with matching headers will be added to dataset

Number	3.1.2.15
Requirement	User will be able to remove vectors that were added to dataset
Note	By removing the Excel attachment from 'Attachments' box

Number	3.1.2.16
Requirement	User will be able to choose which columns to view on the show page
Note	• By choosing headers from list in 'Columns' option • User's preference will be saved

Confidential

Number	3.1.2.17
Requirement	User will be able to add attachments to datasets
Note	

Number	3.1.2.18
Requirement	User will be able to change the number of vectors displayed on single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 documents per page

Number	3.1.2.19
Requirement	User will be able to re-order columns
Note	

Number	3.1.2.20
Requirement	User will be able to sort datasets by any column in display in ascending/descending order
Note	

Number	3.1.2.21
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	3.1.2.22
Requirement	User filter/sort selection will be saved when user move to other pages
Note	

Number	3.1.2.23
Requirement	User will be able to adjust the column data type: string, number, date or time
Note	- By choosing a value from the 'Type' option - When creating a new dataset or when saving an excel attachment as dataset in experiments page

Number	3.1.2.24
Requirement	User will be able to view all dataset rows/vectors in a list view
Note	

Number	3.1.2.25
Requirement	User will be able to see vector numbering
Note	

Number	3.1.2.26
Requirement	User will be able to redirect to vector origin experiment page
Note	

3.2 Labguru SOPs module^(SOP)

3.2.1 SOP index page

Number	3.2.1.1
Requirement	User will be able to create a new SOP
Note	- By clicking on the 'New SOP' button - Clicking on 'New SOP' directs User into a new SOP page

Number	3.2.1.2
Requirement	User will be able to filter SOPs by tags
Note	

Number	3.2.1.3
Requirement	User will be able to search SOPs by name
Note	

3.2.2 SOP show page

Number	3.2.2.1
Requirement	User will be able to add a comment to SOP
Note	

Number	3.2.2.2
Requirement	User will be able to tag an SOP
Note	By clicking the 'Tag this item' button

Number	3.2.2.3
Requirement	User will be able to link SOP to other Labguru pages
Note	- By clicking on 'Linked resources' button in the 'More' drop down menu - By clicking on 'Add link' button on the 'Linked resources' box

Confidential

Number	3.2.2.4
Requirement	User will be able to add: • A text element • An attachments element • A spreadsheet element To section in SOP
Note	

Number	3.2.2.5
Requirement	User will be able to add protocols to the SOP
Note	

Number	3.2.2.6
Requirement	User will be able to delete a protocol from the SOP
Note	By clicking on the delete icon next to the protocol's name

Number	3.2.2.7
Requirement	User will be able to create a PDF copy of the SOP
Note	By clicking the 'Print' button

Number	3.2.2.8
Requirement	Authorized user will be able to sign and witness the SOP
Note	• User is authorized if he has 'Can witness' privilege and he didn't create the SOP • When another user is editing SOP page, the sign function will be disabled, and a message will appear on hover indicating the page is being edited

Number	3.2.2.9
Requirement	User will be able to change the title of the SOP
Note	

Number	3.2.2.10
Requirement	User will be able to assign the SOP to another member
Note	

Number	3.2.2.11
Requirement	User will be able to start a new folder from SOP
Note	• By clicking on the 'Start experiments' button • User can select destination project for the new folder, and the folder's name • In the new folder, an experiment will be created for each protocol added to the SOP

Number	3.2.2.12
Requirement	User will be able to set visibility permissions for SOP
Note	By clicking on 'Visible to...' in the 'More' menu

Number	3.2.2.13
Requirement	User will be able to view previous versions of the SOP in PDF format
Note	By clicking on 'Version history' in the 'More' menu

Number	3.2.2.14
Requirement	User will be able to delete SOP
Note	

Number	3.2.2.15
Requirement	User will be able to access folders started from the SOP
Note	Folder started from SOP will be created with a link to SOP

Number	3.2.2.16
Requirement	User will be able to export SOP as a Word document
Note	By clicking on the Export to Word button

Number	3.2.2.17
Requirement	User will be able to move elements in SOP section up/down
Note	By dragging the element to the desired position

3.2.3 Archiving SOP

Number	3.2.3.1
Requirement	User will be able to archive SOP
Note	

Number	3.2.3.2
Requirement	User will not be able to edit archived SOP
Note	

Number	3.2.3.3
Requirement	User will not be able to start a folder from an archived SOP
Note	

Number	3.2.3.4
Requirement	User will not be able to link archived SOP to other pages
Note	

Number	3.2.3.5
Requirement	User will not be able to delete archived SOP
Note	

Number	3.2.3.6
Requirement	User with 'Can administer' privilege will be able to re-activate archived SOP
Note	

Number	3.2.3.7
Requirement	User will be able view list of archived SOPs
Note	By clicking on 'View archived SOPs' from SOPs index page

Number	3.2.3.8
Requirement	User will be able to restore archived SOP
Note	- By clicking on the 'Restore' link via 'More' menu - Only SOP owner or admin member

3.2.4 Gantt chart in SOP

Number	3.2.4.1
Requirement	User will be able to view Gantt chart representation of protocols in SOP
Note	- Under 'Gantt' chart in SOP page - The Gantt will show representation of the SOP, with protocols ordered the same way as in SOP

Number	3.2.4.2
Requirement	User will be able to drag protocol to change its start day
Note	- By default, all protocols added to the SOP will have day 1 as start day - SOP duration will be updated when protocols are moved

Number	3.2.4.3
Requirement	User will not be able to drag a protocol to start before day 1
Note	

Number	3.2.4.4
Requirement	User will be able to expand and contract experiment width to change its duration
Note	

Number	3.2.4.5
Requirement	User will be able to add dependencies
Note	- As describe in 'Dependencies' section 2.2.3 - Dependencies will be kept when experiments are created from the SOP

Number	3.2.4.6
Requirement	User will be able to Delete dependencies
Note	

Number	3.2.4.7
Requirement	User will be able to go to the protocol page by clicking on its block
Note	

Number	3.2.4.8
Requirement	User will see for each protocol its name, start day and duration (in days)
Note	In the information columns on the left

3.3 Labguru Recipes module (RCP)

3.3.1 Recipes index page

Number	3.3.1.1
Requirement	User will be able to add a recipe
Note	

Number	3.3.1.2
Requirement	User will be able to filter recipes by tags
Note	

Number	3.3.1.3
Requirement	User will be able to search recipes by name
Note	

3.3.2 Recipe show page

Number	3.3.2.1
Requirement	User will be able to create a PDF file of the recipe
Note	By clicking the 'Print' button

Number	3.3.2.2
Requirement	User will be able to edit: 1. Source of a recipe 2. Ph field 3. Concentration field 4. Final volume field 5. Storage conditions field
Note	

Number	3.3.2.3
Requirement	User will be able to duplicate a recipe
Note	

Number	3.3.2.4
Requirement	User will be able to link recipe to other pages
Note	• By clicking the linked resources button on the 'More' menu • By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	3.3.2.5
Requirement	User will be able to tag a recipe
Note	

Number	3.3.2.6
Requirement	User will be able to add an ingredient to the recipe
Note	

Number	3.3.2.7
Requirement	User will be able to set the name of a new ingredient

Confidential

Note	In the 'Add recipe ingredient' pop up
------	---------------------------------------

Number	3.3.2.8
Requirement	User will be able to add a consumable as an ingredient to the recipe
Note	By typing the consumable name, or part of its name (3 letters minimum) and then selecting it from the drop-down list in 'Name' field

Number	3.3.2.9
Requirement	User will be able to set the MW/concentration of the new ingredient
Note	- In the 'Add recipe ingredient' pop up - Molecular weight should be added by using the unit gr/mol - Concentration should be added by using the unit mM

Number	3.3.2.10
Requirement	User will be able to set the final concentration of the new ingredient
Note	- In the 'Add recipe ingredient' pop up - Concentration should be added by using the unit mM

Number	3.3.2.11
Requirement	User will be able to set the weight/volume of the new ingredient
Note	- In the 'Add recipe ingredient' pop up - Weight should be added by using the unit gr - Volume should be added by using the unit ml

Number	3.3.2.12
Requirement	User will be able to edit an ingredient's name, MW/concentration, Final concentration and Weight/volume
Note	By clicking on the pencil icon in the ingredient line

Number	3.3.2.13
Requirement	User will be able to remove an ingredient
Note	By clicking on the 'x' icon in the ingredient line

Number	3.3.2.14
Requirement	User will be able to edit the 'Preparation' section

Confidential

Note	- The title 'Preparation' is not editable - Text, spreadsheet, steps elements can be added to the section
------	--

Number	3.3.2.15
Requirement	Requirement removed
Note	

Number	3.3.2.16
Requirement	User will be able to add an attachment to the preparation section
Note	

Number	3.3.2.17
Requirement	User will be able to add a comment to a recipe
Note	

Number	3.3.2.18
Requirement	User can add an in-line comment to the preparation section
Note	

Number	3.3.2.19
Requirement	User will be able to recalculate recipe ingredients
Note	By changing 'Final volume'

Number	3.3.2.20
Requirement	User will be able to edit the recipe name
Note	By clicking on the current name

Number	3.3.2.21
Requirement	User will be able to download all the recipe's attachments
Note	By clicking on the download icon on the attachments section on the bottom of the screen

Number	3.3.2.22
Requirement	User will be able to move elements in recipe's section up/down
Note	By dragging the element to the desired position

3.4 Labguru Papers module^(PAP)

3.4.1 Papers index page

Number	3.4.1.1
Requirement	User will be able to add a paper
Note	By clicking on the ‘New paper’ button on papers index page

Number	3.4.1.2
Requirement	User will be able to sort papers alphabetically
Note	

Number	3.4.1.3
Requirement	User will be able to sort papers by update date
Note	

Number	3.4.1.4
Requirement	User will be able to search papers by title
Note	By using the search bar

Number	3.4.1.5
Requirement	User will be able to filter papers by tags
Note	

Number	3.4.1.6
Requirement	User will be able to import paper from xml file
Note	

Number	3.4.1.7
Requirement	User will be able to search and view paper from PubMed
Note	By clicking on the ‘PubMed’ button in papers index page

3.4.2 Paper show page

Number	3.4.2.1
Requirement	User will be able to mark paper as favorite
Note	By clicking on the ‘Add to favorite papers’ link in paper show page

Confidential

Number	3.4.2.2
Requirement	User will be able to un-mark paper as favorite
Note	By clicking on the 'Add to favorite papers' link in paper show page

Number	3.4.2.3
Requirement	User will be able to delete a paper
Note	

Number	3.4.2.4
Requirement	User will be able to add/edit: • Authors • Publication date • Journal • Volume • Pages • URL • DOI • Abstract And upload a PDF for a paper
Note	

Number	3.4.2.5
Requirement	User will be able to upload a PDF for a paper
Note	By clicking on 'Upload a PDF'

Number	3.4.2.6
Requirement	Requirement removed
Note	

Number	3.4.2.7
Requirement	User will be able to add a comment to a paper in the discussion section
Note	

Number	3.4.2.8
Requirement	User will be able to link a paper to other pages
Note	• By clicking the 'Linked resources' button on the 'More' menu • By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	3.4.2.9
Requirement	User will be able to tag a paper
Note	By clicking on 'Tag this item'

Number	3.4.2.10
Requirement	User will be able to change the title of a paper
Note	By clicking on the existing title

3.5 Labguru Documents module(DOC)

3.5.1 Documents index page

Number	3.5.1.1
Requirement	User will be able to add a new document
Note	

Number	3.5.1.2
Requirement	User will be able to filter documents by tags
Note	

Number	3.5.1.3
Requirement	User will be able to search documents by name
Note	

Number	3.5.1.4
Requirement	User will be able to change the number of documents he sees on single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 documents per page

Number	3.5.1.5
Requirement	User will be able to sort documents by any column in display in ascending/descending order
Note	

Number	3.5.1.6
Requirement	User will be able to filter documents by any combination of columns in display
Note	

Number	3.5.1.7
Requirement	User will be able to choose which columns to view on the index page
Note	- By choosing headers from list in 'Columns' option - User's preference will be saved

Number	3.5.1.8
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	3.5.1.9
Requirement	User filter/sort selection will be saved when user move to other pages in the collection
Note	

3.5.2 Document show page

Number	3.5.2.1
Requirement	User will be able to edit the title of the document
Note	

Number	3.5.2.2
Requirement	User will be able to assign the document to another member
Note	

Number	3.5.2.3
Requirement	User will be able to tag the document
Note	

Number	3.5.2.4
Requirement	User will be able to print the document to a PDF file
Note	

Number	3.5.2.5
Requirement	User will be able to sign the document
Note	- By clicking the 'Sign' button - Clicking the 'Sign' button opens a 'sign document' popup - Only document owner and admin users can sign the document

	When another user is editing document page, the sign function will be disabled, and a message will appear on hover indicating the page is being edited
--	--

Number	3.5.2.6
Requirement	User will be able to duplicate the document
Note	

Number	3.5.2.7
Requirement	User will be able to link the document to other pages
Note	- By clicking the linked resources button on the 'More' menu - By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	3.5.2.8
Requirement	The document owner will be able to set visibility permissions
Note	

Number	3.5.2.9
Requirement	User will be able to view the version history of the document
Note	

Number	3.5.2.10
Requirement	User will be able to delete a document
Note	

Number	3.5.2.11
Requirement	User will be able to add a comment on the discussion box of the document page
Note	

Number	3.5.2.12
Requirement	User will be able to add: - A text element - An attachment element - A spreadsheet element - An In-line comment To a document
Note	

Number	3.5.2.13
Requirement	User will be able to download all the document's attachments
Note	By clicking on the download icon on the attachments section on the bottom of the screen

Number	3.5.2.14
Requirement	User will be able to add a task from the document page
Note	

Number	3.5.2.15
Requirement	User will be able to export document as a Word document
Note	By clicking on the Export to Word button

Number	3.5.2.16
Requirement	User will be able to move elements in Documents section up/down
Note	By dragging the element to the desired position

3.5.3 Archiving documents

Number	3.5.3.1
Requirement	User will be able to archive the document
Note	

Number	3.5.3.2
Requirement	User will not be able to edit archived document
Note	

Number	3.5.3.3
Requirement	User will not be able to link archived document to other pages
Note	

Number	3.5.3.4
Requirement	User will not be able to delete archived documents
Note	

Number	3.5.3.5
Requirement	User with 'Can administer' privilege will be able to re-activate archived document
Note	

Number	3.5.3.6
Requirement	User will be able to view list of archived documents
Note	By clicking on 'View archived documents' from Documents index page

Number	3.5.3.7
Requirement	User will be able to restore archived documents
Note	• By clicking on 'Restore' link via 'More' menu • Only document owner or admin member

3.6 Labguru Reports module (REP)

3.6.1 Reports index page

Number	3.6.1.1
Requirement	User will be able to create a new report
Note	• By clicking on the 'New report' button • Clicking on 'New report' directs User into a new report page

Number	3.6.1.2
Requirement	User will be able to filter reports by tags
Note	

Number	3.6.1.3
Requirement	User will be able to search reports by name
Note	

Number	3.6.1.4
Requirement	Removed requirement
Note	

Number	3.6.1.5
Requirement	User will be able to change the number of reports he sees on single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 documents per page

Number	3.6.1.6
Requirement	User will be able to sort reports by any column in display in ascending/descending order

Note	
------	--

Number	3.6.1.7
Requirement	User will be able to filter reports by any combination of columns in display
Note	

Number	3.6.1.8
Requirement	User will be able to choose which columns to view on the index page
Note	• By choosing headers from list in 'Columns' option • User's preference will be saved

Number	3.6.1.9
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	3.6.1.10
Requirement	User filter/sort selection will be saved when user move to other pages in the collection
Note	

3.6.2 Report show page

Number	3.6.2.1
Requirement	User will be able to edit the title of the report
Note	

Number	3.6.2.2
Requirement	User will be able to assign the report to another member
Note	

Number	3.6.2.3
Requirement	User can add cover page to the report
Note	

Number	3.6.2.4
Requirement	User will be able to print the report to a PDF file
Note	

Confidential

Number	3.6.2.5
Requirement	User will be able to duplicate the report
Note	

Number	3.6.2.6
Requirement	User will be able to link the report to other pages
Note	• By clicking the linked resources button on the 'More' menu • By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	3.6.2.7
Requirement	The report owner will be able to set visibility permissions for the report
Note	

Number	3.6.2.8
Requirement	User can see the version history of the report
Note	By clicking on 'Version history' in the 'More' button and from the bottom of the page

Number	3.6.2.9
Requirement	User will be able to view new sections that were added to the report as sections on the left-hand side of the page
Note	

Number	3.6.2.10
Requirement	User will be able to add a comment to the discussion box
Note	From the bottom of the page

Number	3.6.2.11
Requirement	User will be able to collapse and open sections on the report
Note	By clicking on the black arrow on the top right corner of the section

Number	3.6.2.12
Requirement	User will be able to add a new section to the report
Note	By clicking on 'New section' adds a custom section to the report

Confidential

Number	3.6.2.13
Requirement	User will be able to re-order sections
Note	• By clicking on ‘Move up’/’Move down’ in the section ‘More’ menu • By dragging and dropping the sections thumbnails to the desired position

Number	3.6.2.14
Requirement	User will be able to add: • A text element • A spreadsheet element • A sketch element • An attachments element • An In-line comment • A reaction element • A compound element to sections in report
Note	From section tool bar and ‘More’ menu

Number	3.6.2.15
Requirement	User will be able to delete a section
Note	By selecting ‘Delete section’ from the section’s ‘More’ menu

Number	3.6.2.16
Requirement	User will be able to sign a report
Note	When another user is editing the report, the sign function will be disabled, and a message will appear on hover indicating the page is being edited

Number	3.6.2.17
Requirement	User will be able to revert a signed report
Note	

Number	3.6.2.18
Requirement	User will be able to reject a signed report
Note	Only users with “can witness” privilege

Number	3.6.2.19
Requirement	User will be able to witness a report
Note	Only users with “can witness” privilege

Number	3.6.2.20
Requirement	User will be able to move elements in report sections
Note	By dragging the element to the desired position

3.6.3 Archiving reports

Number	3.6.3.1
Requirement	User will be able to archive the report
Note	

Number	3.6.3.2
Requirement	User will not be able to edit archived report
Note	

Number	3.6.3.3
Requirement	User will not be able to link archived report to other pages
Note	

Number	3.6.3.4
Requirement	User will not be able to delete archived report
Note	

Number	3.6.3.5
Requirement	User with 'Can administer' privilege will be able to re-activate archived report
Note	

Number	3.6.3.6
Requirement	User will be able view list of archived reports
Note	By clicking on 'View archived reports' from Reports index page

Number	3.6.3.7
Requirement	User will be able to restore the report
Note	• By clicking on 'Restore' link via 'More' menu • Only report owner and admin member

3.7 Labguru Image Bank module (IMG)

3.7.1 Images index page

Number	3.7.1.1
Requirement	User will be able to filter images by: • Owner • State (attached/unattached)
Note	

Number	3.7.1.2
Requirement	User will be able to filter images by tags
Note	

Number	3.7.1.3
Requirement	User will be able to search images by name
Note	

3.7.2 Image show page

Number	3.7.2.1
Requirement	User will be able to download image as it was when it was initially attached
Note	By clicking on the 'download original' button

Number	3.7.2.2
Requirement	User will be able to download annotated image
Note	By clicking on the 'download annotated' button

Number	3.7.2.3
Requirement	User will be able to re-attach image to inventory items
Note	

Number	3.7.2.4
Requirement	User will be able to add a comment
Note	By adding text to the comments box

Number	3.7.2.5
Requirement	User will be able to link an image to other pages
Note	• By clicking the linked resources button on the 'More' menu

	By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page
--	---

Number	3.7.2.6
Requirement	User will be able to delete an image
Note	By clicking on the 'Delete item' in the 'more' menu

4. Inventory

4.1 Biocollections^(BIO)

4.1.1 Index page

4.1.1.1 General

Applicable to all the following biocollection: Antibodies, Bacteria, Cell lines, Flies, Fungi, Genes, Lipids, Plasmids, Primers, Proteins, Sequences, Tissues, Viruses, Worms, Yeasts, Zebrafishes, Compounds, Botany, Rodents and custom collections

Number	4.1.1.1.1
Requirement	User will be able to add a new inventory item
Note	- User will not be able to add a new item without adding its name - Each item has a unique auto-generated Sysid - User will be required to add unique name when unique name field is enabled in the collection - SysID format: {collection letters} - {2 digit year} . {serial number} (e.g., PL-25.0010)

Number	4.1.1.1.2
Requirement	User will be able to filter collection by tags
Note	

Number	4.1.1.1.3
Requirement	User will be able to search for items by name and Sysid
Note	

Number	4.1.1.1.4
Requirement	User will be able to change the number of items in single page
Note	User will see 20 per page by default and can change to 5, 50, 100, and 200 items per page

Number	4.1.1.1.5
Requirement	User will be able to sort the collection by any column in display in ascending/descending order

Note	
------	--

Number	4.1.1.1.6
Requirement	User will be able to filter the collection by any combination of columns in display
Note	

Number	4.1.1.1.7
Requirement	User will be able to choose which columns to view on the index page
Note	- By choosing headers from list in 'Columns' option - User's preference will be saved

Number	4.1.1.1.8
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Number	4.1.1.1.9
Requirement	User will be able to re-order columns
Note	By dragging and dropping the columns

Number	4.1.1.1.10
Requirement	User will be able to delete selected items
Note	- Authorized users only - By selecting items and clicking on 'delete' - User will not be able to delete the item if it was added to experiment/protocol

Number	4.1.1.1.11
Requirement	User will be able to bulk edit selected items
Note	- By marking up to 100 items and clicking on 'Edit' - Purchasable attribute should not be part of the 'Edit selected' option

Number	4.1.1.1.12
Requirement	User will be able to tag items
Note	By selecting items and clicking on 'Tag'

Number	4.1.1.1.13
Requirement	User will be able to view shared repositories for the collection

Confidential

Note	• By clicking on ‘Repositories’ • Not applicable for custom collections
------	--

Number	4.1.1.1.14
Requirement	User will be able to import/update items
Note	By clicking on the cog menu and ‘import/update’

Number	4.1.1.1.15
Requirement	User will be able to export items
Note	• By clicking on the cog menu and ‘Export’ • Export all items by not selecting any stock • Export only selected items by first checking their boxes • Export all filtered items by clicking on the header check box and ‘Select all XXX items’ • Export file will be sent in email • User can select whether to include stocks data or not • User can export sequences from collections that include sequence

Number	4.1.1.1.16
Requirement	User will be able view previous imports/updates
Note	By clicking on the cog menu and ‘Previous imports/updates’

Number	4.1.1.1.17
Requirement	User will be clear all filters/sorts
Note	By clicking on the cog menu and ‘Clear all filters/sorts’

Number	4.1.1.1.18
Requirement	User will be able to edit items values from index page
Note	• By clicking on ‘Edit’ in the left of each row • Applicable to editable fields of type text, numbers, date and pre-defined list • Item’s name must be unique when unique name field is enabled

Number	4.1.1.1.19
Requirement	User will be able to archive items from index page
Note	• By selecting multiple items -> ‘Archive’ in collection index page • User will not be able to archive items that has active stocks

Number	4.1.1.1.20
Requirement	User will be able to see how many items in collection
Note	

4.1.1.2 Section removed

4.1.1.3 'Antibodies page

Number	4.1.1.3.1
Requirement	Requirement removed
Note	

4.1.1.4 'Plasmids page

Number	4.1.1.4.1
Requirement	User will be able to add plasmid by uploading gb file
Note	By clicking on 'GeneBank'

Number	4.1.1.4.2
Requirement	User will be able to add clone plasmids
Note	By clicking on the cog menu and 'Batch Cloning'

Number	4.1.1.4.3
Requirement	User will be able to export plasmids with all sequences
Note	By clicking on the cog menu and 'Export sequences'

4.1.1.5 'Proteins page

Number	4.1.1.5.1
Requirement	User will be able to export proteins with all sequences
Note	By clicking on the cog menu and 'Export sequences'

4.1.1.6 'Compounds page

Number	4.1.1.6.1
Requirement	User will be able to search for compounds by Full structure search
Note	By clicking on 'Structure search'

Number	4.1.1.6.2
Requirement	User will be able to search for compounds by Sub structure search
Note	By clicking on 'Structure search'

Number	4.1.1.6.3
Requirement	User will be able to view structures in compounds index page
Note	

Number	4.1.1.6.4
Requirement	User will be able to add structure with SMILE or InChI notation
Note	- By adding a value in 'Structure (SMILES)' or 'Structure (InChI)' fields in an import template - When adding value in both fields, the structure will be according to InChI notation

Number	4.1.1.6.5
Requirement	User will be able to export structures in SMILES and InChI notation
Note	

4.1.1.7 Botany's 'Seeds' page

Number	4.1.1.7.1
Requirement	User will be able to add new plants from selected seeds
Note	By clicking on 'Saw seeds'

Number	4.1.1.7.2
Requirement	User will be able to add new crosses from selected seeds
Note	By clicking on 'Cross seeds'

4.1.1.8 Botany's 'Plants' page

Number	4.1.1.8.1
Requirement	User will be able to add new seeds from selected plants
Note	By clicking on 'Harvest plants'

Number	4.1.1.8.2
Requirement	User will be able to add new crosses from selected plants

Note	By clicking on 'Cross plants
------	------------------------------

Number	4.1.1.8.3
Requirement	User will not be able to delete plant/s used for crosses
Note	

4.1.1.9 Rodents' 'Specimens' page

Number	4.1.1.9.1
Requirement	User will be able to reserve selected specimens
Note	By clicking on 'Make reservation'

Number	4.1.1.9.2
Requirement	User will be able to add treatment to selected specimens
Note	By clicking on 'Add treatment'

4.1.2 Biocollection's show page

Applicable to all of the following biocollection: Antibodies, Bacteria, Cell lines, Flies, Fungi, Genes, Lipids, Plasmids, Primers, Proteins, Sequences, Tissues, Viruses, Worms, Yeasts, Zebrafishes, Compounds, Botany, Rodents and Custom collections.

Number	4.1.2.1
Requirement	User will be able to edit item
Note	- By clicking on 'Edit' - Item's name must be unique when unique name field is enabled

Number	4.1.2.2
Requirement	User will be able to duplicate item
Note	- By clicking on 'Duplicate', purchasable attributes are not duplicated to new item - Item's name must be unique when unique name field is enabled - When an item contains linked resources, they should be associated with the duplicated item

Number	4.1.2.3
Requirement	User will be able to delete item
Note	- By clicking on 'delete' - Item cannot be deleted if it has aliquoted stocks - Item cannot be deleted if it was added to experiments/protocols

Number	4.1.2.4
Requirement	User will be able to print pdf version of item show page
Note	By clicking on 'print'

Number	4.1.2.5
Requirement	User will be able to add item to the shopping list
Note	By clicking on 'Add to shopping list'

Number	4.1.2.6
Requirement	User will be able to share item to repository
Note	- By clicking on 'Share to repository' - Not applicable to custom collections

Number	4.1.2.7
Requirement	User will be able to link item to other pages
Note	By clicking on 'add link' in 'Linked resources' box

Number	4.1.2.8
Requirement	User will be able to tag item
Note	By clicking on 'add tag' in 'Tags' box

Number	4.1.2.9
Requirement	User will be able to add tasks to item
Note	By clicking on 'add task' in 'Tasks' box

Number	4.1.2.10
Requirement	User will be able to attach files to item
Note	By clicking on 'Upload a file' box

Number	4.1.2.11
Requirement	User will be able to add comments to item
Note	By typing in 'Discussion box' and clicking 'Post'

Number	4.1.2.12
Requirement	User will be able to archive item
Note	- By clicking on 'Archive' - Only user with 'Inventory Edit' permission will be able to archive item - Archived items will keep their links, tags, tasks, comments, and attachments that were added to them before the archive

	• User will be able to archive item that was added to the samples element • User will not be able to archive item that has active stocks • User will not be able to archive item that was added to the shopping list
--	--

4.1.3 Biocollection's show pages tabs

4.1.3.1 'Info' tab

Number	4.1.3.1.1
Requirement	User will be view item information
Note	Appears in all biocollection show pages

Number	4.1.3.1.2
Requirement	User will be able to see values entered to fields customized for collection
Note	For parents fields- user will be redirected to each parent item show page, by clicking on the item's name

4.1.3.2 'Sequence & Map' tab

Number	4.1.3.2.1
Requirement	User will be able to view sequences that were added to item and their maps
Note	Appears in items from Plasmids, Sequences and custom collections

Number	4.1.3.2.2
Requirement	User will see DNA viewer for DNA sequences
Note	When multiple sequences are added to an item, DNA viewer will show map of the first DNA sequence added

4.1.3.3 'Order History' tab

Number	4.1.3.3.1
Requirement	User will be able to view the order history of an item
Note	Appears in items from collections that have the purchasable attributes selected

Number	4.1.3.3.2
Requirement	User will be able to add item to shopping list

Note	
------	--

4.1.3.4 'Stocks' tab

Number	4.1.3.4.1
Requirement	User will be able to view stocks that belong to item
Note	- Appears in all biocollection show pages -except Rodent Strains, Rodent Specimens) - Only active stocks will be shown

Number	4.1.3.4.2
Requirement	User will be able to view consumed stocks that belonged to item
Note	- By clicking on 'Show consumed stocks' - User will be able to go back to view active stocks

Number	4.1.3.4.3
Requirement	User will be able to add new stocks to item
Note	- By clicking on 'Add Stock' - Limited to 5000 stocks

Number	4.1.3.4.4
Requirement	Requirement removed
Note	

Number	4.1.3.4.5
Requirement	Requirement removed

Number	4.1.3.4.6
Requirement	User will be able to edit selected stocks
Note	- By marking up to 100 stocks and clicking on 'Edit selected' - By clicking on the pencil icon to edit one stock

Number	4.1.3.4.7
Requirement	User will be able to tag stock/s
Note	By selecting stocks and clicking on 'Tag'

Number	4.1.3.4.8
Requirement	User will be able to mark stock/s as consumed
Note	- By selecting stocks and clicking on 'Mark as consumed' - By clicking on the stamp icon to consume one stock

Number	4.1.3.4.9
Requirement	User will be able to dispatch stock/s
Note	By selecting stocks and clicking on 'Dispatch'

Number	4.1.3.4.10
Requirement	User will be able to print labels for stocks/s
Note	By selecting stocks and clicking on 'Print labels'

Number	4.1.3.4.11
Requirement	User will be able to add stock/s to samples basket
Note	By selecting stocks and clicking on 'Add to basket'

Number	4.1.3.4.12
Requirement	User will be able to view inventory item total amount of material by weight\volume
Note	- Total material will be presented with the largest units of measurement in stocks - Calculation of total material is done on stocks with default units only

4.1.3.5 'Datasets' tab

Number	4.1.3.5.1
Requirement	User will be able to view the datasets vectors that item was added to
Note	- From 'Datasets' tab in item page - Appears in all biocollection show pages - except Rodent Strains

Number	4.1.3.5.2
Requirement	User will be able to export the list of dataset vectors to csv
Note	

Number	4.1.3.5.3
Requirement	User will be able to export the list of dataset vectors to xlsx

Note	In the xlsx file, groups of similar vectors will be separated to sheets
------	---

4.1.3.6 ‘Experiments’ tab

Number	4.1.3.6.1
Requirement	User will be able to see experiment that item was added to
Note	- From ‘Experiments’ tab in item page - See requirements :2.4.1.7 - 2.4.1.16 User will be able to view the experiment in which the compound was created as part of a reaction/compound element

Number	4.1.3.6.2
Requirement	User will be able to view the number of stocks that were used in the experiment
Note	

Number	4.1.3.6.3
Requirement	User will be able to view the number of stocks that were marked as output in that experiment
Note	

4.1.3.7 ‘Protocols’ tab

Number	4.1.3.7.1
Requirement	User will be able to view protocols that the item was added to
Note	- From ‘Protocols’ tab - See requirements: 2.5.1.1-2.5.1.8

4.1.3.8 ‘Specimens’ tab

Number	4.1.3.8.1
Requirement	User will be able to view specimens’ belong to item
Note	- From ‘Specimens’ tab - Appears in Rodent Strains collection

Number	4.1.3.8.2
Requirement	User will be able to create new specimens to item
Note	- By clicking on ‘Create specimens’ - Available fields in ‘Add specimens’ page: name, DoB, status and cage

	• Required fields” Number of specimens
--	--

Number	4.1.3.8.3
Requirement	User will be able to mark specimen/s as reserved
Note	• By selecting specimens and clicking on ‘Make reservation’ • Possible fields in ‘Reserve specimen’ page: Select more specimens, Description and notification settings. • Mandatory fields: Required at date
Number	4.1.3.8.4
Requirement	User will be able to edit specimen/s
Note	By selecting specimens and clicking on ‘Edit selected

4.1.3.9 ‘Tissue & Samples’ tab

Number	4.1.3.9.1
Requirement	User will be able to view tissues derived from the item
Note	• From ‘Tissues & Samples’ tab • Appears in items from Rodent Specimens collection

Number	4.1.3.9.2
Requirement	User will be able to a new derived tissue
Note	By clicking on ‘Add new tissue’

4.1.3.10 ‘Treatments’ tab

Number	4.1.3.10.1
Requirement	User will be able to view treatments specimen was added to
Note	• From ‘Treatments’ tab • Appears in items from Rodent Specimens collection

4.1.3.11 ‘Sequences’ tab

Number	4.1.3.11.1
Requirement	User will be able to view sequences added to item
Note	• From ‘Sequences’ tab • Appears by default in items in collections: Genes, Primers, Proteins

4.1.3.12 'Derived item' tab

Number	4.1.3.12.1
Requirement	User will be able to view items that were derived from item
Note	- From 'Derived item' tabs - Appears in items from collection that have derived collections set up in the collection customization page

4.1.4 Inventory item 'Add new/Edit item' page

Number	4.1.4.1
Requirement	User will be able to populate fields customized for collection
Note	Fields types can be from the following: - selected pre-defined fields - custom fields - Parent items

Number	4.1.4.2
Requirement	User will not see blacklisted pre-defined fields
Note	

Number	4.1.4.3
Requirement	User will not be able to save if required fields left empty
Note	

Number	4.1.4.4
Requirement	User will be able to cancel without adding new/editing item
Note	
Number	4.1.4.5
Requirement	User will not see read only fields in page
Note	

Number	4.1.4.6
Requirement	User will be able to add parent items
Note	- Item can have multiple parent from the same collection and/or from multiple collections - Limited to 10 parents selection

Number	4.1.4.7
Requirement	User will be able to remove chosen parent item/s
Note	

Number	4.1.4.8
Requirement	User will not be able to add an item as a parent more than once
Note	

Number	4.1.4.9
Requirement	User will be able to mark item as non-commercial
Note	- Only in collection with purchasable attributes - User will be required to fill 'Produced by', 'Manufacturer' and 'Catalog number' will be greyed out - By default, Non-commercial will be turned off

Number	4.1.4.10
Requirement	User will be required to add Manufacturer and Catalog Number
Note	- For catalog number- maximum character length is 50 - Only in commercial item 'Produced by' will be greyed out - Catalog number must be unique - If the manufacturer exists, it will be displayed by typing the first letter

Number	4.1.4.11
Requirement	User will be required to add a unique name to item
Note	When unique name field is enabled in the collection

Number	4.1.4.12
Requirement	User will be able to add size and size units to item
Note	- Only in collection with purchasable attributes - User will be able to select units from drop down list - By default, user will see kg unit

Number	4.1.4.13
Requirement	User will be able to add units remarks to item
Note	

4.1.5 Labguru antibody module (ANT)

Number	4.1.5.1
Requirement	User will be able to add applications to an antibody
Note	By selecting an application from the 'Application' drop down menu

	• By clicking on the ‘Add another application’ link • in the ‘Add new antibody’ page
--	---

Number	4.1.5.2
Requirement	Removed requirement
Note	

Number	4.1.5.3
Requirement	User will be able to add dilution
Note	By filling the field in the ‘Add new antibody’ page

Number	4.1.5.4
Requirement	User will be able to delete an application
Note	By clicking on the ‘x’ icon in the ‘Add new antibody’ page

Number	4.1.5.5
Requirement	User will be able to define the organism of the antibody
Note	By clicking on the ‘Raised in’ and selecting an organism in the ‘Add new antibody’ page

Number	4.1.5.6
Requirement	User will be able to define what the antibody reacts with
Note	By clicking on the ‘Reacts with’ in the ‘Add new antibody’ page

Number	4.1.5.7
Requirement	User will be able to set clonality
Note	By clicking on ‘Clonality’ drop down menu in the ‘Add new antibody’ page

4.1.6 Section removed

4.1.7 Labguru Bacteria module

Number	4.1.7.1
Requirement	User will be able to add strain to bacterium
Note	By filling the ‘Strain’ field in ‘Add new bacterium’ page

Number	4.1.7.2
Requirement	User will be able to add Sensitivity to a bacterium
Note	By filling the ‘Sensitivity’ field in the ‘Add new bacterium’ page

4.1.8 Labguru Botany module(BTN)

4.1.8.1 Seeds

Number	4.1.8.1.1
Requirement	User will be able to set: • Genotype • Phenotype • Male/Female parent For a seed
Note	In the 'Add new seed' page

Number	4.1.8.1.2
Requirement	User will be able to saw seeds
Note	By selecting seeds and clicking on 'Saw seeds' button in seeds index page

Number	4.1.8.1.3
Requirement	User will be able to cross seeds
Note	By selecting seeds and clicking on 'Cross seeds' button in seeds index page

4.1.8.2 Plants

Number	4.1.8.2.1
Requirement	User will be able to set: • Genotype • Phenotype • Generation • Harvest date For a plant
Note	In 'Add new plant' page

Number	4.1.8.2.2
Requirement	User will be able to harvest plants
Note	By selecting plants and clicking on the 'Harvest plants' button In plants index page

Number	4.1.8.2.3
Requirement	User will be able to cross plants
Note	By selecting seeds and clicking on 'Cross plants' button in plants index page

Number	4.1.8.2.4
Requirement	User will not be able to delete plant if it was used for crosses
Note	

4.1.9 Labguru Cell lines module [\(CEL\)](#)

Number	4.1.9.1
Requirement	User will be able to set: • Organism • Medium & serum • Tissue For a cell line
Note	In 'Add new cell line' page

4.1.10 Labguru Compounds Module [\(COM\)](#)

Number	4.1.10.1
Requirement	User will be able to set: • Cas • Molar mass • Density • Boiling point • Melting point • Hazards • Structure For a compound
Note	

Number	4.1.10.2
Requirement	User will be able to view structures in compounds show page
Note	

Number	4.1.10.3
Requirement	User will be able to view auto-populated chemical properties: • Formula • MW (g/mol) • cLogP • TPSA • Elemental analysis • SMILES • Canonical SMILES • InChI • HBD • HBA • Non-hydrogen atom count • Rotatable bond count

	• FractionCSP3 • Aromatic ring count • Aliphatic ring count • Lipinski Rule of Five For a compound
Note	• Admin users can select in the collection customization page which properties to show • Fields are uneditable

Number	4.1.10.4
Requirement	User will not be able to edit compound that was created in an experiment
Note	

Number	4.1.10.5
Requirement	User will be able to import/bulk import compounds and their stock using SDF file
Note	• By clicking on the import option • The stocks will be stored in a temporary box named after the SDF file • User receives an email notification when the import is completed • Existing compounds- ○ New item will not be created ○ A new stock will be added to the original item

Number	4.1.10.6
Requirement	User will be able to create a compound with an existing structure
Note	• Duplicates will be linked to their original compound

Number	4.1.10.7
Requirement	User will be able to select compound attributes for import from SDF file
Note	• By selecting from a list of existing compound/stock attribute fields • User will be able to add new compound/stock field • User will be able to ignore attribute field • Attribute cannot be selected more than once

Number	4.1.10.8
Requirement	User will be able to view Import compound summary
Note	• Including entries found, new compounds created, duplicated entries

	• Linked to compounds index page
--	--

Number	4.1.10.9
Requirement	User will be able to view legacy data
Note	• Legacy data includes previously selected chemical properties and manually entered Formula values • Will be available only for compounds with legacy data presentation enabled • The data will be presented on compound's show page • It will not be presented in the edit page • The legacy data will be removed when the structure is edited

Number	4.1.10.10
Requirement	User will be able to add/remove stereochemistry representations for a chemical structure
Note	Each representation will be registered with a unique structure ID

4.1.11 Labguru Flies module (FLY)

Number	4.1.11.1
Requirement	User will be able to set: • Genotype • Phenotype For a fly
Note	In 'Add new Fly' page

4.1.12 Labguru Fungi module(FNG)

Number	4.1.12.1
Requirement	User will be able to set: • Species • Phenotype • Host • Virulent • Sporulate • Mycelia • Sporocarps (fruiting bodies) For a fungus
Note	In 'Add new Fungus' page

4.1.13 Labguru Genes module([GEN](#))

Number	4.1.13.1
Requirement	User will be able to set: • Expression location • Metabolic pathway For a gene
Note	In 'Add new gene' page

Number	4.1.13.2
Requirement	User will be able to add new associated primer to gene
Note	In 'Primers' tab in gene show page

Number	4.1.13.3
Requirement	User will be able to see list of primers associated to gene
Note	In 'Primers' tab in gene show page

4.1.14 Labguru Lipids module([LIP](#))

Number	4.1.14.1
Requirement	User will be able to set: • Molecular formula/weight • Cas number • Stock solution prep • Solubility • Conditions for use • Conditions for storage • Safety information • Impurities For a lipid
Note	In 'Add new lipid' page

4.1.15 Labguru Plasmids module([PLA](#))

Number	4.1.15.1
Requirement	User will be able to set: • Base vector • Length • Host • Usage • Resistance • Insert • Clone number • Sequence For a plasmid
Note	In 'Add new plasmid' page

Number	4.1.15.2
Requirement	User will be able to add a new plasmid with sequence
Note	By filling the sequence field with sequence in the plasmids show page

Number	4.1.15.3
Requirement	User will be able to add sequence from GenBank to plasmid
Note	

Number	4.1.15.4
Requirement	User will be able to update a GenBank file
Note	By clicking on the 'Upload GenBank file' link in plasmids index page

Number	4.1.15.5
Requirement	Requirement removed
Note	

4.1.16 Labguru Primer module (PRM)

Number	4.1.16.1
Requirement	User will be able to set: • Label/tag • Orientation • Tm • Bind position • Fragment size • Organism • Application (pcr / rt / sdm / etc) • Restriction site • Sequence For a primer
Note	In 'Add new primer' page

4.1.17 Labguru Protein module (PRO)

Number	4.1.17.1
Requirement	User will be able to set: • Gene • Species • Mutations • Chemical modifications • Tag • Purification method • Molecular weight • Extinction coefficient at 280 nm Storage buffer

	- Storage buffer - Storage temperature For a protein
Note	In 'Add new protein' page

4.1.18 Labguru Sequences module [\(SEQ\)](#)

Number	4.1.18.1
Requirement	User will be able to set: - Type - Sequence - Accession - Organism For a sequence
Note	In 'Add new sequence' page

Number	4.1.18.2
Requirement	Requirement removed
Note	

4.1.19 Labguru Tissue module [\(TIS\)](#)

Number	4.1.19.1
Requirement	User will be able to set: - Species - Genotype/phenotype - Animal details - Tissue type - Harvest date - Fixation/embedding procedure - Applications - Storage conditions For a tissue
Note	In 'Add new tissue' page

Number	4.1.19.2
Requirement	User will be able to add a tissue via rodent specimen's show page
Note	By clicking on the plus icon in the 'Tissue & Samples' tab

Number	4.1.19.3
Requirement	User will be able to see the specimen that the tissue was derived from in the tissue's show page
Note	In case a tissue was originated from specimens via rodent specimen's show page

4.1.20 Labguru Worms module (WRM)

Number	4.1.20.1
Requirement	User will be able to set: • Made by • Genotype/phenotype • Gene • Mutagen • Growth conditions • Outcrossed • Dauer defective For a worm
Note	In 'Add new worm' page

4.1.21 Labguru yeasts module (YST)

Number	4.1.21.1
Requirement	User will be able to set: • Reproduction • Genetic background • Transgenic features • Phenotype For a yeast
Note	In 'Add new yeast' page

4.1.22 Labguru Viruses module (VIR)

Number	4.1.22.1
Requirement	User will be able to set: • Gene/insert • Virus type • Plasmid • Serotype/strain • Mutations/deletions • Tag • Selectable marker • Producer cell type • Viral coat • Tropism • Species • Storage conditions For a virus
Note	In 'Add new virus' page

4.1.23 Labguru Zebrafishes module(ZEB)

Number	4.1.23.1
Requirement	User will be able to set: • Genotype • Gene • Mutant allele • Transgene • Phenotype For a zebrafish
Note	In 'Add new zebrafish' page

4.1.24 Labguru Rodents module4.1.24.1 Rodent strains(RDS)

Number	4.1.24.1.1
Requirement	User will be able to set: • Transgene • Genotype/phenotype For a rodent strain
Note	In 'Add new rodent strain' page

4.1.24.2 Rodent cages(RDC)

Number	4.1.24.2.1
Requirement	User will be able to set: • Breeding • Storage location For a rodent cage
Note	In 'Add new rodent cage' page

Number	4.1.24.2.2
Requirement	Cage storage location & specimens will be shown in cages index page by default
Note	In rodent cage index page

Number	4.1.24.2.3
Requirement	User will be able to add a litter to a mating cage
Note	By clicking on the 'Create litter specimen' link in the cage show page on 'Specimens' tab

Number	4.1.24.2.4
Requirement	User will be able to add a specimen to a cage
Note	By clicking on the 'Add specimens' link in the cage show page on 'Specimens' tab

Number	4.1.24.2.5
Requirement	User will be able to reserve rodent specimens
Note	By clicking on the 'Make reservation' link in the cage show page on 'Specimens' tab

Number	4.1.24.2.6
Requirement	User will be able to add treatments to specimens in cage
Note	By clicking on the 'Add Treatment' link in the cage show page on 'Treatment' tab

4.1.24.3 Rodent specimen_(RDP)

Number	4.1.24.3.1
Requirement	User will be able to set: • Ethics approval id • Experiment • Strain • Sex • Female/male parent • Date of birth/death • Ear/coat tag • status For a rodent specimen
Note	In 'Add new rodent specimen' page

Number	4.1.24.3.2
Requirement	User will be able to reserve specimens
Note	By clicking on the 'Make reservation' link in rodent specimen index page

Number	4.1.24.3.3
Requirement	User will be able to add treatments to specimens
Note	By clicking on the 'Add Treatment' link in rodent specimen index page

4.1.24.4 Rodent treatments(RDT)4.1.24.4.1 Treatments index page

Number	4.1.24.4.1.1
Requirement	User will be able to filter treatments by member
Note	By selecting a member via the 'Fitter by user' drop down menu

Number	4.1.24.4.1.2
Requirement	User will be able to print treatments displayed on screen
Note	By clicking the 'Print' link

Number	4.1.24.4.1.3
Requirement	User will be able to search for a treatment by name
Note	By entering the treatment name in search box

Number	4.1.24.4.1.4
Requirement	User will be able to view Today's treatments details
Note	The table will include the following data: • Is the treatment done or not • Treatment data • Administered by • Specimen cage

Number	4.1.24.4.1.5
Requirement	User will be able to view Future treatments details
Note	The table will include the following data: • User will be able to filter by date • Treatment data • Administered by • Specimen cage

Number	4.1.24.4.1.6
Requirement	User will be able to view treatments history details
Note	The table will include the following data: • User will be able to filter by date • Is the treatment done or not • Treatment data • Administered by • Applied by • Specimen • Cage

Number	4.1.24.4.1.7
Requirement	User will be able to add new treatment
Note	By clicking on the 'Add treatment' button

4.1.24.4.2 'Add Treatment' page

Number	4.1.24.4.2.1
Requirement	User will be able to connect treatment to a protocol
Note	By clicking on the 'Protocol' field a list of protocols will be displayed and the option to load steps to it.

Number	4.1.24.4.2.2
Requirement	User will be able to load a section from a connected protocol into treatment
Note	By clicking on the 'Step' field

Number	4.1.24.4.2.3
Requirement	User will be able to manage protocols that will be used as rodent protocols
Note	- By clicking on the 'manage protocols' link - By dragging and dropping protocols between the two lists: List no1 - Rodent protocols List no2 - All protocols

Number	4.1.24.4.2.4
Requirement	User will be required to set name of treatment
Note	By filling the 'Name' field

Number	4.1.24.4.2.5
Requirement	User will be required to select scheduled date of treatment
Note	- By clicking on the 'Scheduled at' field - A date picker will be opened

Number	4.1.24.4.2.6
Requirement	User will be able to select a cage for treatment
Note	- By clicking on the 'applied at' field - A list of cages will be displayed - User will be able to search for a cage in the search box at the top of the list

Number	4.1.24.4.2.7
Requirement	User will be able to add more than one cage for the treatment

Note	By clicking on the ‘add more cages’
------	-------------------------------------

Number	4.1.24.4.2.8
Requirement	User will be able to add remarks to the treatment
Note	By filling the ‘Remarks’ field

4.1.25 Archived item

Number	4.1.25.1
Requirement	User will be able to see all archived items of a collection in a list
Note	

Number	4.1.25.2
Requirement	User will not be able to: • Edit • Duplicate • Archive • Add to shopping list • Share to repository • Tag • Add sequence to • Link to other pages • Add tasks to • Add comments to • Add attachments to archived item
Note	

Number	4.1.25.3
Requirement	User will be able to restore archived item
Note	Item cannot be restored when unique field contains a value already exists

Number	4.1.25.4
Requirement	User will be able to delete archived item
Note	

Number	4.1.25.5
Requirement	User will be able to print archived item to pdf
Note	

4.2 Labguru Consumables module (CONS)

4.2.1 Consumables index page

Number	4.2.1.1
Requirement	User will be able to add new consumable
Note	- By clicking on the 'New Material' - User will not be able to add a new consumable without adding its name

Number	4.2.1.2
Requirement	User will be able to add consumables to the shopping list
Note	- By clicking on the cart icon - User will be required to enter quantity & price - By clicking on the 'Add to shopping list' button in the 'Add to shopping list' pop up

Number	4.2.1.3
Requirement	User will be able to see how many purchase requests are pending in the shopping list for a consumable
Note	- By hovering the cart icon - The number is displayed under the cart

Number	4.2.1.4
Requirement	User will be able to see how many stocks of a consumable exists in storage
Note	By clicking the 'Tube' link in the 'Location' column

Number	4.2.1.5
Requirement	User will be able to see a consumable manufacturer information
Note	By clicking on the manufacturer name in 'Manufacturer' column

Number	4.2.1.6
Requirement	User will be able to export items
Note	- By clicking on the cog menu and 'Export' - Export all items by not selecting any stock - Export only selected items by first checking their boxes - Export all filtered items by clicking on the header check box and 'Select all XXX items' - Export file will be sent in email - User will be able to export items without stocks data

4.2.2 Consumable show page

Number	4.2.2.1
Requirement	User will be able to duplicate consumable
Note	

Number	4.2.2.2
Requirement	User will be able to move consumable to different collections
Note	By using the 'Move to' link in the consumable show page

Number	4.2.2.3
Requirement	User will be able to delete consumable
Note	

Number	4.2.2.4
Requirement	User will be able to print consumable page to pdf
Note	

4.2.3 'Add new Material' page

Number	4.2.3.1
Requirement	Requirement removed
Note	

Number	4.2.3.2
Requirement	Requirement removed
Note	

4.3 Labguru Shopping List module (SHP)4.3.1 'Add item to shopping list' popup -from inventory item show page

Number	4.3.1.1
Requirement	User will be notified when an item already exists on the shopping list when he tries to add it again
Note	A notification will appear in the popup that opens while adding an item to the list

Number	4.3.1.2
Requirement	User will be able to set: • Quantity • Price • Currency • Remarks

	• Description • Attachments • Pre-defined custom fields For the ordered item
Note	• By default, quantity will be set to 1 • Max quantity is 100

4.3.2 'Shopping list' index page

Number	4.3.2.1
Requirement	User will be able to add items to shopping list
Note	• By using 'Select material to add to shopping list' • Only for items from collections with 'Purchasable attributes' selected

Number	4.3.2.2
Requirement	User will be able to view all items added to shopping list sorted by request date in descending order (newest is first)
Note	

Number	4.3.2.3
Requirement	Only users with 'Can approve order' privilege will be able to approve an order
Note	• Under 'Status' column • When order is approved, the button is replaced with a 'Submit order' button

Number	4.3.2.4
Requirement	Only users with 'Can approve order' privilege will be able to submit an order
Note	• Under 'Status' column • When order is submitted, the button is replaced with a 'Mark as arrived' button

Number	4.3.2.5
Requirement	Only users with 'Can approve order' privilege will be able to approve and submit an order in one action
Note	• Under 'Status' column • When order is approved and submitted, the button is replaced with a 'Mark as arrived' button

Number	4.3.2.6
Requirement	User will be able to cancel an order

Note	• Users without 'Can approve orders' will be able to cancel only orders added by them • Users with 'Can approve orders' will be able to cancel all orders
------	--

Number	4.3.2.7
Requirement	User will be able to mark an order as arrived
Note	• By clicking the 'Mark as arrived' button • Clicking the 'Mark as arrived' button opens a popup where User can add a location for the arrived stock • When stocks are added to box, user is required to select box position for the stocks • When stocks are added to general storage location or gridless box, user is required to select how many stocks arrived • When stocks are created, user will see notification with link to item • When user mark order as partially arrived, the order is not removed from shopping list and quantity is updated • User cannot add more stocks than ordered

Number	4.3.2.8
Requirement	User will be able to edit the order
Note	• By clicking the pencil icon on the first column • Clicking the pencil icon opens a popup where User can edit the price, quantity, budget, remarks, attachments, custom fields data

Number	4.3.2.9
Requirement	User will be able to filter the orders
Note	• By clicking the 'Filter items' button • Clicking the 'Filter items' button opens a popup where User can choose how he wishes to filter the shopping list • Can filter by: Requested by, Cat. Number, Manufacturer, Representative, Price range, Budget, Status, Request date, Order date and custom fields created in account

Number	4.3.2.10
Requirement	User will be able to export the shopping list to a CSV file
Note	• By clicking the 'Export to CSV' button • Clicking the 'Export to CSV' button automatically download the CSV file

Number	4.3.2.11
Requirement	User will be able to create a purchase request form
Note	• By clicking the 'Print' button • Clicking the 'Print' button opens a popup where User can choose the vendors that will appear on the file • The purchase request form is a PDF file • Only orders that were not submitted will be added to the pdf print out

Number	4.3.2.12
Requirement	User will be able to search for items by item name
Note	By using the search bar

Number	4.3.2.13
Requirement	User will be able to clear filters after searching the shopping list or filtering it
Note	By clicking the 'Clear filters' button

4.3.3 Section Removed

4.3.4 Order History page

Number	4.3.4.1
Requirement	User will be able to view all items that were ordered and marked as arrived or partially arrived
Note	• By clicking the 'Order history' button • The item's name (linked to the item's show page), catalog number, and manufacturer (linked to the company page) appear on the 'Material' column

Number	4.3.4.2
Requirement	User will be able to edit orders in 'Order history'
Note	• By clicking the pencil icon on the first column • Clicking the pencil icon opens a popup where User can edit the price, quantity, budget, remarks, attachments, custom fields data

Number	4.3.4.3
Requirement	User will be able to filter the order history
Note	• By clicking the 'Filter items' button

	Clicking the 'Filter items' button opens a popup where User can choose how he wishes to filter the order history
--	--

Number	4.3.4.4
Requirement	User will be able to export the order history to a CSV file
Note	- By clicking the 'Export to CSV' button - Clicking the 'Export to CSV' button automatically download the CSV file

Number	4.3.4.5
Requirement	User will be able to search for items on the order history by item name
Note	By using the search bar

4.3.5 Service orders page

Number	4.3.5.1
Requirement	User will be able to add a new service order
Note	- By clicking the 'New service request' button - Clicking the 'New service request' button opens a popup where User selects a member (from a dropdown), company, adds a description, quantity, price, and budget

Number	4.3.5.2
Requirement	User will be able to export the service orders list to a CSV file
Note	- By clicking the 'Export to CSV' button - Clicking the 'Export to CSV' button automatically download the CSV file

Number	4.3.5.3
Requirement	User will be able to edit a service order
Note	- By clicking the pencil icon on the first column - Clicking the pencil icon opens a popup where User can edit the service order details

Number	4.3.5.4
Requirement	User will be able to see who added the service order to the shopping list

Note	- User who added the item appears on the 'User' column with his private name and image - Clicking on User's image directs User to User's profile page - Image of past users is greyed out
------	---

Number	4.3.5.5
Requirement	User will be able to search for service orders by description
Note	By using the search bar

Number	4.3.5.6
Requirement	User will be able to mark a service order as arrived
Note	- By clicking the 'Mark as arrived' button - When a purchase order is being added, the 'Mark as arrive' button appears on the status column - When clicking the 'Mark as arrived' button, the service order is being removed from the list

Number	4.3.5.7
Requirement	User will be able to cancel a service order
Note	- Users without 'Can approve orders' will be able to cancel only service orders added by them - Users with 'Can approve orders' will be able to cancel all service orders

4.4 Labguru Companies module(CMP)

Number	4.4.1
Requirement	User will be able to view all companies that were added to the account
Note	- By clicking on 'Manage vendors' in Cog menu in consumables index page - Including the company's name, the representative (if added), URL (if added), phone number (if added) - The name of the company is linked to the company's page

Number	4.4.2
Requirement	User will be able to add a new company
Note	- By clicking the 'New company' button - The name field is mandatory

	• When a company is being created a green 'Company was successfully created' message appears on the top of the page • When a company is being created User is being directed to the company show page • The URL is linked to the company's website
--	--

Number	4.4.3
Requirement	User will be able to define a company as a distributor or a representative
Note	• By marking this option while creating or editing a company • Companies which are defined as 'Distributor/representative', will appear on the 'Distributed/Represented by' dropdown that appears on the 'New company' and 'Edit company' pages • When marking a company as 'Distributor/representative', the 'Distributed/Represented by' dropdown becomes disabled • If company X was defined as a representative and was selected on the 'Distributed/Represented by' dropdown in company Y show page, company Y will appear on the 'Represents' field of company X details

Number	4.4.4
Requirement	User will be able to edit a company
Note	• By clicking the 'Edit company' button on the company show page • When a company is being created a green 'Company was successfully updated' message appears on the top of the page • When a company is being updated User is being directed to the company show page

Number	4.4.5
Requirement	User will be able to search for a company by name
Note	• By using the search bar under the 'Companies' title • The search locates text within all fields

Number	4.4.6
Requirement	User will be able to add comments to companies
Note	

Number	4.4.7
Requirement	User will be able to link companies to other pages
Note	

Number	4.4.8
Requirement	User will be able to delete company
Note	• Company which contains materials cannot be deleted • Double confirmation is required

4.5 Customize collection page (CMZ)

4.5.1 Pre-defined list

Number	4.5.1.1
Requirement	User will be able to select from list of pre-defined fields
Note	• By checking/unchecking fields from the ‘customize fields’ list • Checked fields will be seen in the index table of the collection

Number	4.5.1.2
Requirement	User will be able to select chemical properties: • Formula • MW (g/mol) • cLogP • TPSA • Elemental analysis • SMILES • Canonical SMILES • InChI • HBD • HBA • Non-hydrogen atom count • Rotatable bond count • FractionCSP3 • Aromatic ring count • Aliphatic ring count • Lipinski Rule of Five that will be auto-calculated and displayed in compound page
Note	• Relevant only for ‘Compounds’ collection • By checking/unchecking fields from the ‘customize fields’ list

Number	4.5.1.3
Requirement	User will be able to mark field name of an item as unique/non-unique
Note	As default, 'Allow Non-Unique Names' will be marked

Number	4.5.1.4
Requirement	User will be able to choose to display legacy data on compound page
Note	• Relevant only for 'Compounds' collection • It will be hidden by default • Legacy data includes previously selected chemical properties and manually entered formula values

4.5.2 Custom fields

Number	4.5.2.1
Requirement	User will be able to add custom field to collections, stocks, equipment and shopping list
Note	• By clicking on 'add custom field' • User will be required to enter field's name • User will not be able to save an item with a non unique value when unique name field is enabled

Number	4.5.2.2
Requirement	User will be able to set a custom field as: • Regular field (text or numbers) • Only numbers • URL • Date • Pre-defined list • Text area • Pre-defined list with multiple select • Simple long text area • Date and time • Account LOV
Note	

Number	4.5.2.3
Requirement	User will be able to add a hint to a custom field
Note	Adding a hint adds it on the 'Add new' form and the 'edit item' page

Number	4.5.2.4
Requirement	Requirement removed
Note	

Number	4.5.2.5
Requirement	User will be able to mark a custom field as required
Note	

Number	4.5.2.6
Requirement	User will be able to mark a custom field as read-only
Note	Read-only fields can be edited only from the API

Number	4.5.2.7
Requirement	User will not be able to mark a custom field as required and read-only
Note	Checkboxes need to be mutually exclusive

Number	4.5.2.8
Requirement	User will be able to edit custom field's properties
Note	The field's type will not be editable

Number	4.5.2.9
Requirement	User will be able to delete custom field
Note	Deletion will not be allowed for custom fields that are: • Not empty for all items in the collection • Currently used to sort or filter the index page

4.5.3 Derived Collections

Number	4.5.3.1
Requirement	User will be able to add derived collections to collection
Note	• This method is used in all inventory modules • Only active collections can be added as derived

Number	4.5.3.2
Requirement	User will be able to mark a collection as derived of itself
Note	This method is used in all inventory modules

Number	4.5.3.3
Requirement	User will be able to mark a derived collection as required field
Note	This method is used in all inventory modules

Number	4.5.3.4
Requirement	User will be able to unmark a derived collection as required field
Note	This method is used in all inventory modules

Number	4.5.3.5
Requirement	User will be able to remove a derived collection
Note	This method is used in all inventory modules

5. Storage & Equipment module

5.1 Labguru Storages module [\(STG\)](#)

5.1.1 Storages page

Number	5.1.1.1
Requirement	User will be able to view all storage locations that were created in the account
Note	

Number	5.1.1.2
Requirement	User will be able to see the details of a storage location
Note	- By clicking on the storage location name - Including Id, type, description, temperature and barcode

Number	5.1.1.3
Requirement	User will be able to delete a storage location
Note	

Number	5.1.1.4
Requirement	User will be able to edit the storage location
Note	- By clicking the 'Edit' button - Clicking the 'Edit' button enables User to change the name, temperature, description, barcode and storage location (choose from a dropdown) of the storage location - The type field is not editable if user has defined it as a rack or if boxes were already added to the storage - User cannot edit an existing storage type to a rack

Number	5.1.1.5
Requirement	User will be able to see stocks located in a storage location
Note	

Number	5.1.1.6
Requirement	User will be able to access stocks located in a storage location
Note	Stocks names are linked to the stock's pages

Number	5.1.1.7
Requirement	User will be able to view equipment located in a storage location

Note	Stocks names are linked to the stock's pages
------	--

Number	5.1.1.8
Requirement	User will be able to access a page of equipment located in a storage location
Note	By clicking on the equipment's name

Number	5.1.1.9
Requirement	User will be able to view boxes located in a storage location
Note	

Number	5.1.1.10
Requirement	User will be able to access a box located in a storage location
Note	By clicking on the box name

Number	5.1.1.11
Requirement	User will be able to add/change the image of the storage page
Note	By clicking on the 'Update image' link

Number	5.1.1.12
Requirement	User will be able to add a new storage location
Note	• By clicking on the 'add new storage ' button • Storages can be added inside other storages • Clicking on the 'add new storage location' button opens an add new location screen on the right side of the screen • All users can add storage locations (except read-only users) • Storage location type is mandatory • When selecting rack, rows and columns are mandatory

Number	5.1.1.13
Requirement	User will be able to hide/un hide the content of a storage location
Note	By clicking on the black arrow next to the storage location name

Number	5.1.1.14
Requirement	User will be able to add a comment on a storage location
Note	

Number	5.1.1.15
Requirement	User will be able to add new stocks to storage location
Note	

Number	5.1.1.16
Requirement	User will be able to rearrange storage location
Note	- User will get notification when storage location changed successfully or failed - Storage cannot be moved to rack type storage

Number	5.1.1.17
Requirement	User will be able to rearrange order of sub-storage locations located under the same storage
Note	

5.1.2 Section removed

Number	5.1.2.1
Requirement	Requirement removed
Note	

Number	5.1.2.2
Requirement	Requirement removed
Note	

Number	5.1.2.3
Requirement	Requirement removed
Note	

Number	5.1.2.4
Requirement	Requirement removed
Note	

Number	5.1.2.5
Requirement	Requirement removed
Note	

Number	5.1.2.6
Requirement	Requirement removed
Note	

5.2 Labguru Boxes module (BOX)

5.2.1 Boxes index page

Number	5.2.1.1
Requirement	User will be able to add a new box

Confidential

Note	- By clicking on the ‘New Box’ button - The max box size limit is 50 * 50
------	--

Number	5.2.1.2
Requirement	User will be able to delete/bulk delete a box
Note	By selecting box/s and clicking on ‘Delete’

Number	5.2.1.3
Requirement	User will be able to print labels for a box
Note	By selecting box/s and clicking on ‘Print labels’

Number	5.2.1.4
Requirement	User will be able to tag a box
Note	By selecting box/s and clicking on ‘Tag’

Number	5.2.1.5
Requirement	Requirement removed
Note	

Number	5.2.1.6
Requirement	User will be able to sort boxes in ascending or descending order by any column header
Note	- By clicking on a header name - By clicking on the kendo menu and select one of the sorting options

Number	5.2.1.7
Requirement	User will be able to search for boxes by name, ID or scan a barcode
Note	

Number	5.2.1.8
Requirement	User will be able to filter boxes by any combination of multiple columns
Note	

Number	5.2.1.9
Requirement	User will be able to filter boxes by to tags
Note	By clicking on tag’s name

Confidential

	- User can cancel selection by clicking on same tag again - Clicking on multiple tags will filter for boxes that are tags with all selected tags
--	---

Number	5.2.1.10
Requirement	User will be able to add and remove columns from the display
Note	- By clicking on the kendo menu -> 'Columns' and checking the wanted columns - User preferences will be saved

Number	5.2.1.11
Requirement	User will be able to set the number of boxes displayed per page
Note	Default number is 20 boxes per page. Can be changed to 50,100 or 200

5.2.2 Box show page

Number	5.2.2.1
Requirement	User will be able to edit box
Note	- By clicking the 'Edit' button - Editable fields are box name, barcode, content type, box dimensions, color, owner, description and storage location

Number	5.2.2.2
Requirement	User can export box as an xlsx file
Note	By clicking on 'Export to Excel'

Number	5.2.2.3
Requirement	User can print box as a PDF file
Note	- By clicking on the 'Print' button - Box name, owner, color, type, size, storage location and content type

Number	5.2.2.4
Requirement	User can print a label for a box
Note	By clicking on 'Print label' in box toolbar

Number	5.2.2.5
Requirement	User can view a box details: • Owner • Content-type • Personal/shared • Color • Storage • Barcode
Note	• The owner's name is linked to the member show page • The storage is links to all locations in the storage path

Number	5.2.2.6
Requirement	User will be able to see all stocks from the box ordered by position
Note	An illustration of the stock type and color appears in the box position

Number	5.2.2.7
Requirement	User will be able to see stock information from box page
Note	When hovering stock icon

Number	5.2.2.8
Requirement	User will be able to mark stock/s from the box as consumed
Note	• For a single stock, by clicking 'Consume' in stock position • For multiple stocks, by selecting their positions and then clicking on 'Consume' in 'Bulk actions' menu

Number	5.2.2.9
Requirement	User will be able to edit stock/s from the box
Note	• For a single stock, by clicking 'Edit' in stock position • For multiple stocks, by selecting their positions and then clicking on 'edit' in the 'Bulk actions' menu

Number	5.2.2.10
Requirement	User will be able to dispatch stock/s from the box
Note	• For a single stock, by clicking 'Dispatch' in stock position • For multiple stocks, by selecting their positions and then clicking on 'Dispatch' in the 'Bulk actions' menu

Confidential

Number	5.2.2.11
Requirement	User will be able to print labels for stock/s from the box
Note	- For a single stock, by clicking the 'Print Label' button in stock position - For multiple stocks, by selecting their positions and then clicking on 'Print labels' in the 'Bulk actions' menu

Number	5.2.2.12
Requirement	User will be able to add stock to selected position in the box
Note	By clicking the '+' icon inside the box

Number	5.2.2.13
Requirement	User will be able to tag box
Note	By clicking on 'Add tag'

Number	5.2.2.14
Requirement	User will be able to link the box to other Labguru pages
Note	By clicking on 'add link'

Number	5.2.2.15
Requirement	User will be able to add an attachment to box
Note	- By clicking on 'Upload a file' - By dragging file and drop it into 'Upload a file'

Number	5.2.2.16
Requirement	User will be able to delete box
Note	By clicking on 'Delete'

Number	5.2.2.17
Requirement	Box owner will be required to confirm box deletion action
Note	By checking to check-boxes and clicking on 'I understand - delete this box' in a pop-up that will appear after clicking on 'Delete'

Number	5.2.2.18
Requirement	User will be able to see the name of the member that added the box, date of creation and date of last update
Note	On the bottom of the page

	The member's name is linked to the member's page
--	--

Number	5.2.2.19
Requirement	User will be able to add stock/s to samples basket
Note	- By selecting stock positions and clicking on 'Add to basket' in the 'Bulk actions' menu - In gridless boxes - by selection stocks and clicking on 'Add to basket'

Number	5.2.2.20
Requirement	User will be able to see list of members with permission to edit box
Note	

Number	5.2.2.21
Requirement	Authorized user will be able to edit box permissions
Note	- By default, box permissions will be given to all users - User will be able to select members and/or teams from the offered members list - User will be able to select 'All' or 'None' from the offered members list - Box owner and admin users' permissions are uneditable - Description of the actions a user with box permission is able to perform is described in Appendix A - Stocks located in a box 'inherent' box permissions - Stock owner will hold their permission to use stock, even if located in a box he has no permission to

Number	5.2.2.22
Requirement	User will be able to change stock position from the box to selected one
Note	By clicking 'Reposition' in stock position

Number	5.2.2.23
Requirement	User will be able to duplicate stock/s from the box
Note	- By clicking 'Duplicate' in stock position - For multiple stocks, by dragging the mouse over their positions - Stocks created via pooling cannot be duplicated

Number	5.2.2.24
Requirement	User will not be able to delete box with aliquoted stocks
Note	

5.3 Labguru Stocks module (STK)

5.3.1 Stocks index page

Number	5.3.1.1
Requirement	User will be able to sort stocks in ascending or descending order by any column header
Note	• By clicking on a header name • By clicking on the kendo menu and select one of the sorting options

Number	5.3.1.2
Requirement	User will be able to search for stocks by name, ID or scan a barcode
Note	

Number	5.3.1.3
Requirement	User will be able to filter stocks by any combination of multiple columns
Note	

Number	5.3.1.4
Requirement	User will be able to filter stocks by to tags
Note	• User can cancel selection by clicking on same tag again • Clicking on multiple tags will filter for stocks that are tags with all selected tags

Number	5.3.1.5
Requirement	User will be able to add and remove columns from the display
Note	• By clicking on the kendo menu -> 'Columns' and checking the wanted columns • User preferences will be saved

Number	5.3.1.6
Requirement	User will be able to edit selected stocks
Note	By marking up to 100 stocks and clicking on 'Edit selected'

Number	5.3.1.7
Requirement	User will be able to tag stock(s)

Confidential

Note	By selecting stocks and clicking on ‘Tag’
------	---

Number	5.3.1.8
Requirement	User will be able to mark stock(s) as consumed
Note	By selecting stocks and clicking on ‘Mark as consumed’

Number	5.3.1.9
Requirement	User will be able to dispatch stock(s)
Note	By selecting stocks and clicking on ‘Dispatch’

Number	5.3.1.10
Requirement	User will be able to print labels for stock(s)
Note	By selecting stocks and clicking on ‘Print labels’

Number	5.3.1.11
Requirement	User will be able to import/update stocks
Note	By clicking on the cog icon and ‘import/update’

Number	5.3.1.12
Requirement	User will be able to view list of consumed stocks
Note	By clicking on ‘View consumed stocks’

Number	5.3.1.13
Requirement	Requirement removed
Note	

Number	5.3.1.14
Requirement	User will be able to set the number of stocks displayed per page
Note	Default number is 20 stocks per page. Can be changed to 50,100 or 200

Number	5.3.1.15
Requirement	User will be able to adjust the width of the columns
Note	- By dragging and dropping the columns' borders - Auto-adjust by double clicking the border

Number	5.3.1.16
Requirement	User will be able to add stock/s to sample basket
Note	By selection stocks and clicking on ‘Add to basket’

Confidential

Number	5.3.1.17
Requirement	User will be able to export stocks
Note	• By clicking on cog menu and 'Export' • Export all stocks by not selecting any stock • Export only selected stock by first checking their boxes • Export all filtered stocks by clicking on the header check box and 'Select all XXX stocks'

Number	5.3.1.18
Requirement	User will be able to view previous import/update events
Note	By clicking on cog menu and 'Previous Imports/Updates'

Number	5.3.1.19
Requirement	User will be clear all filters/sorts
Note	By clicking on the cog menu and 'Clear all filters/sorts'

Number	5.3.1.20
Requirement	User will be able to select all filtered stocks
Note	By clicking on 'Select all' after selecting all stocks in page

5.3.2 Consumed stocks index page

Number	5.3.2.1
Requirement	User will be able to search for consumed stocks by name
Note	

Number	5.3.2.2
Requirement	User will be able to filter consumed stocks using the Advanced search
Note	Searching according to the advanced search criteria, an item filters only the relevant rows

Number	5.3.2.3
Requirement	User will be able to filter consumed stocks by to tags
Note	

Number	5.3.2.4
Requirement	User will be able to customize table view
Note	The number of columns and order is being updated according to the table view preferences

Number	5.3.2.5
Requirement	User will be able to edit selected consumed stocks
Note	By marking up to 100 stocks and clicking on 'Edit selected'

Number	5.3.2.6
Requirement	User will be able to print labels for consumed stocks/s
Note	By selecting stocks and clicking on 'Print labels'

Number	5.3.2.7
Requirement	User will be able to delete consumed stocks
Note	By selecting stocks and clicking on 'Delete'

Number	5.3.2.8
Requirement	User will be able to set the number of consumed stocks displayed per page
Note	Default number is 20 stocks per page. Can be changed to 50,100 or 200

5.3.3 Stock show page

5.3.3.1 General

Number	5.3.3.1.1
Requirement	User will be able to see stock details: • Name and ID • Collection • Storage location • Owner • Stored/frozen by • Stored/ frozen on • Privacy • Concentration • Volume • Weight • Units • Stock color • Barcode • Size
Note	

Number	5.3.3.1.2
Requirement	User will be able to edit stock
Note	

Number	5.3.3.1.3
Requirement	User will be able to mark stock as consumed
Note	

Confidential

Number	5.3.3.1.4
Requirement	Requirement removed
Note	

Number	5.3.3.1.5
Requirement	User will be able to dispatch stock
Note	

Number	5.3.3.1.6
Requirement	User will be able to print label of stock
Note	By clicking on 'Print label'

Number	5.3.3.1.7
Requirement	User will be able to tag stock
Note	By clicking on 'add tag' in 'Tags' side box

Number	5.3.3.1.8
Requirement	User will be able to link stock to other pages
Note	By clicking on 'add link' in 'Linked Resources' side box

Number	5.3.3.1.9
Requirement	User will be able to add attachment to stock
Note	By clicking on 'Upload a file' in 'Attachments' side box

Number	5.3.3.1.10
Requirement	User will see list of attached images/files
Note	

Number	5.3.3.1.11
Requirement	User will be able to download a file attached to stock
Note	

Number	5.3.3.1.12
Requirement	User will be able to view an image attached to stock in full size
Note	

Number	5.3.3.1.13
Requirement	User will be able to remove a file attached to stock
Note	

Number	5.3.3.1.14
Requirement	User will be able to add a comment
Note	

Number	5.3.3.1.15
Requirement	User will see list of experiments that the stock was added to (via Samples element)
Note	- By clicking on the 'Experiment' tab - List includes: - Experiment name - Project name - Folder name - Experiment owner - If stock is output of the experiment - See requirements 2.4.1.7 - 2.4.1.16

Number	5.3.3.1.16
Requirement	User will be able to add stocks to samples basket
Note	By clicking on 'Add to basket'

Number	5.3.3.1.17
Requirement	User will be able to add a task to stock
Note	By clicking on 'Add Task'

5.3.3.2 'Stock history' tab

Number	5.3.3.2.1
Requirement	User will be able to see a stock's history of events
Note	- From 'History' tab - Events will be ordered chronologically

Number	5.3.3.2.2
Requirement	User will see the following activities: - Creation - Location change - Ownership changes - Usage in experiments - Editing (other than location, owner) - Consumption - Aliquoting - Dispatch
Note	

Number	5.3.3.2.3
Requirement	User will see for each event - name of event, event details, member name and id, date of event

Note	
------	--

Number	5.3.3.2.4
Requirement	User will be able to filter stock history by member
Note	

Number	5.3.3.2.5
Requirement	User will be able to filter stock history be date
Note	

Number	5.3.3.2.6
Requirement	User will be able to print stock history
Note	

5.3.3.3 Consumed stocks

Number	5.3.3.3.1
Requirement	User will see name and collection of the item associated with stock
Note	

Number	5.3.3.3.2
Requirement	User will be able to edit consumed stock
Note	

Number	5.3.3.3.3
Requirement	User will be able to delete consumed stock
Note	Deletion is not allowed if the consumed stock was added to experiment/ has aliquots

Number	5.3.3.3.4
Requirement	User will be able to print label of consumed stock
Note	By clicking on 'Print label'

Number	5.3.3.3.5
Requirement	User will be able to tag consumed stock
Note	By clicking on 'add tag' in 'Tags' side box

Number	5.3.3.3.6
Requirement	User will be able to link consumed stock to other pages
Note	By clicking on 'add link' in 'Linked Resources' side box

Confidential

Number	5.3.3.3.7
Requirement	User will be able to add attachment to consumed stock
Note	By clicking on 'Upload a file' in 'Attachments' side box

Number	5.3.3.3.8
Requirement	User will see list of attached images/files to consumed stock
Note	

Number	5.3.3.3.9
Requirement	User will be able to download a file attached to consumed stock
Note	

Number	5.3.3.3.10
Requirement	User will be able to view an image attached to consumed stock in full size
Note	

Number	5.3.3.3.11
Requirement	User will be able to remove a file attached to consumed stock
Note	

Number	5.3.3.3.12
Requirement	User will be able to add a comment
Note	

Number	5.3.3.3.13
Requirement	User will see list of experiments that the consumed stock was added to
Note	By clicking on the 'Experiment' tab

Number	5.3.3.3.14
Requirement	User will be able to add a task to consumed stock
Note	By clicking on 'Add Task'

Number	5.3.3.3.15
Requirement	User will be able to aliquot consumed stock
Note	- By clicking on 'Add Aliquots' - Limited to 5000 aliquots

5.3.4 Stocks audit

Number	5.3.4.1
Requirement	The action 'show' will be shown in audit
Note	

Number	5.3.4.2
Requirement	The action 'audit' will be shown in audit
Note	

Number	5.3.4.3
Requirement	Requirement removed
Note	

Number	5.3.4.4
Requirement	The creation of a new stock from item show page will be shown in audit
Note	

Number	5.3.4.5
Requirement	The creation of a new stocks via collection import will be shown in audit
Note	

Number	5.3.4.6
Requirement	The creation of a new stocks via stock import will be shown in audit
Note	

Number	5.3.4.7
Requirement	The creation of a new stocks via box will be shown in audit
Note	

Number	5.3.4.8
Requirement	Shows 'comment added' action in audit
Note	

Number	5.3.4.9
Requirement	Shows 'tag added' action in audit
Note	

Number	5.3.4.10
Requirement	Shows ‘tag removed’ action in audit
Note	

Number	5.3.4.11
Requirement	Shows ‘mark as used’ action in audit
Note	

Number	5.3.4.12
Requirement	Shows ‘label printed’ action in audit
Note	

Number	5.3.4.13
Requirement	Shows ‘link added’ action in audit
Note	

Number	5.3.4.14
Requirement	Shows ‘link removed’ action in audit when removing single link
Note	

Number	5.3.4.15
Requirement	Shows ‘link removed’ action in audit when removing all links
Note	

Number	5.3.4.16
Requirement	Shows ‘attachment added’ action in audit
Note	

Number	5.3.4.17
Requirement	Shows ‘attachment archived’ action in audit
Note	

Number	5.3.4.18
Requirement	Shows ‘added to experiment’ action in audit
Note	All elements added

Number	5.3.4.19
Requirement	Shows ‘update’ action in audit
Note	

5.3.5 Low stock alert

Number	5.3.5.1
Requirement	Requirement removed
Note	

Number	5.3.5.2
Requirement	Requirement removed
Note	

Number	5.3.5.3
Requirement	Subscribed user will be notified when low stock alert is triggered
Note	- By default, user will see notifications in Dashboard, Shopping list and Weekly Digest email - Notifications will disappear when stocks are added, and the stock count/weight/volume of stocks passes the set threshold

Number	5.3.5.4
Requirement	User will be able to know when a stock alert has already been set for an item
Note	In 'Stocks' tab in item show page

Number	5.3.5.5
Requirement	User will be able to know when the low stock alert is triggered
Note	In 'Stocks' tab in item show page

Number	5.3.5.6
Requirement	User will be able to subscribe to a low stock alert
Note	- Alert can be based on stock count, total volume or total weight - For volume/weight - up to 4 decimal numbers - For volume/weight - user need to select units

Number	5.3.5.7
Requirement	User will be able to edit low stock alert
Note	

Number	5.3.5.8
Requirement	User will be able to unsubscribe from a low stock alert
Note	By clicking on 'Stock alert' and then on 'Unsubscribe'

5.3.6 Aliquoting stocks5.3.6.1 Creating new aliquots

Number	5.3.6.1.1
Requirement	User will be able to aliquot a stock from the stock show page
Note	• By clicking on 'Add Aliquots' from stock show page • Limited to 5000 aliquots

Number	5.3.6.1.2
Requirement	User will be able to select how much to aliquot from the original stock
Note	• Default will be to aliquot all the amount of the stock • User will be able to see the amount available in the stock, and then change it • The field should be numbers only • User will not be able to select a number that is bigger than the origin stock or below 0 • Aliquots volume unit will be the same as the origin's unit

Number	5.3.6.1.3
Requirement	User will be able to aliquot stocks without volume calculation
Note	• Applicable to volume-less stocks and stocks with value in volume • If original stock has volume, by setting the amount to reduce to 0 • If original stock is volume-less, default value to reduce should be 0 • When aliquoting without volume calculation, aliquots created without volume

Number	5.3.6.1.4
Requirement	User will be required to select how many stocks to create or what should be the final volume of each stock
Note	• User will have 2 windows to select either number of stocks or final volume • Once user add something to one of the fields, the other field is automatically calculated: amount to aliquot = number of aliquots * final amount per aliquot • Both fields will be numbers only

Number	5.3.6.1.5
Requirement	User will be able to select if origin stock will be consumed after aliquoting
Note	Default will not consume the stock

Number	5.3.6.1.6
Requirement	User will be required to select location for aliquots
Note	By selecting storage location or box

Number	5.3.6.1.7
Requirement	User will be required to select positions in box in the exact number as the number of aliquots selected
Note	- User will not be able to select more positions than aliquots - As long as the number of aliquots and number of positions are not the same, the indicator's color will be red, and the 'Next' button will be unresponsive - When the correct number of positions was selected, the indicator will turn green, and the 'Next' button will be active

Number	5.3.6.1.8
Requirement	User will be able to update aliquots non-required fields
Note	- The data will be added to all created aliquots - By default, the fields 'Name', 'Type', 'Color', 'Content', 'lot' and 'Concentration' will have the origin stock's values - User will not be able to edit 'Volume', 'Content' and 'Collection' - By default, the member who did the aliquoting, will be added as the owner of the aliquots

Number	5.3.6.1.9
Requirement	User will be required to update aliquots required custom fields
Note	The data will be added to all created aliquots

Number	5.3.6.1.10
Requirement	User will see a stock that was created by aliquoting the following information: 'Aliquoted by' (member who did the aliquot) 'Aliquoted on' (Date of aliquoting) 'Aliquoted from' (origin stock)
Note	Name of origin stock will be link to origin stock page

Number	5.3.6.1.11
Requirement	User will see notification at left bottom page when aliquoting has started/completed
Note	

Number	5.3.6.1.12
Requirement	User will be able to view the completed new aliquots
Note	By clicking on the link in the pop-up message or under Aliquots tab

5.3.6.2 'Aliquots' tab in stock show page

Number	5.3.6.2.1
Requirement	User will be able to view all aliquots made from a stock
Note	• By clicking on 'Aliquots' tab in stock show page • List will show active and consumed aliquots

Number	5.3.6.2.2
Requirement	User will be able to choose which columns to view
Note	• By choosing headers from list in 'Columns' option • User's preference will be saved

Number	5.3.6.2.3
Requirement	User will be able to sort aliquots in ascending or descending order by column headers
Note	• By clicking on a header name • Applicable for '#', 'Name' columns

Number	5.3.6.2.4
Requirement	User will be able to filter aliquots by column headers
Note	• By clicking in the filter line • Applicable for '#', 'Name' columns

Number	5.3.6.2.5
Requirement	User will be able to filter aliquots by tags
Note	

Number	5.3.6.2.6
Requirement	User will be able to change the number of aliquots in single page
Note	User will see 20 per page by default and can change to 5, 50, 100, and 200 aliquots per page

Number	5.3.6.2.7
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

Confidential

Number	5.3.6.2.8
Requirement	User will be able to re-order columns
Note	By dragging and dropping the columns

Number	5.3.6.2.9
Requirement	User will be able to edit selected aliquots
Note	By marking up to 100 aliquots and clicking on 'Edit selected'

Number	5.3.6.2.10
Requirement	User will be able to tag aliquots
Note	By selecting aliquots and clicking on 'Tag'

Number	5.3.6.2.11
Requirement	User will be able to mark aliquots as consumed
Note	By selecting aliquots and clicking on 'Mark as consumed'
Number	5.3.6.2.12
Requirement	User will be able to dispatch aliquots
Note	By selecting aliquots and clicking on 'Dispatch'

Number	5.3.6.2.13
Requirement	User will be able to print labels for aliquots
Note	By selecting aliquots and clicking on 'Print labels'

5.3.7 Add/Edit stock page

5.3.7.1 General

Number	5.3.7.1.1
Requirement	User will be required to select stock type
Note	As default stock type will be 'Tube'

Number	5.3.7.1.2
Requirement	User will be able to select concentration, weight and volume units
Note	By selecting from default units or creating a new one

Number	5.3.7.1.3
Requirement	Admin user will be able to create custom unit for concentration, weight and volume
Note	

Number	5.3.7.1.4
Requirement	Admin user will be able to remove custom unit for concentration, weight and volume
Note	Unit that is used in existing stocks cannot be removed

Number	5.3.7.1.5
Requirement	User will be required to add stock owner
Note	- As default stock owner will be the member who added the stock - By clicking on 'Owner' search option or selecting from members list

Number	5.3.7.1.6
Requirement	User will be able to add stored/frozen by
Note	- As default Stored/frozen by name will be the member who added the stock - By clicking on 'Stored/frozen by' search option or selecting from active members

Number	5.3.7.1.7
Requirement	User will be able to add 'Stored/frozen on'
Note	- As default Stored/frozen on will be 'today's' date - User will be able to change the date by clicking on 'Stored/frozen on' field

5.3.7.2 From Item show page, Shopping list, Sample element, Sample basket

Number	5.3.7.2.1
Requirement	User will be required to choose location
Note	By clicking on 'Suggested boxes' search option or selecting from 'Other storage locations'

Number	5.3.7.2.2
Requirement	User will be able to select box from 'Suggested boxes'
Note	When selecting from 'Suggested boxes' only boxes with content type 'General' or same as the collection of the stock, will be available

Number	5.3.7.2.3
Requirement	User will be able to select location from 'Other storage location'
Note	User can select either a box (unstored boxes included) or a storage location

	- When selecting from ‘other storage location’ only boxes with content type ‘General’ or that is the same as the collection of the stock, will be available - In rack cell location type, user will be able to select a box within the rack cell
--	---

Number	5.3.7.2.4
Requirement	User will be able to search for location in ‘Other storage locations’
Note	User will see the path of search result

Number	5.3.7.2.5
Requirement	User will be required to name the stock
Note	As default stock name will be equal to item

Number	5.3.7.2.6
Requirement	User will be required to select privacy status from ‘Private’ or ‘Public’
Note	As default privacy will be ‘Public’

Number	5.3.7.2.7
Requirement	User will be required to select box positions
Note	- Only when selected location is grid box - When adding from Shopping list user can select positions as number of stocks in the order or less

Number	5.3.7.2.8
Requirement	User will be required to choose how many stocks to add
Note	- Only when selected location is gridless box - Limited to 5000 items - When adding from shopping list, user cannot add more stocks than ordered

5.3.7.3 Add stock from Box

Number	5.3.7.3.1
Requirement	User will be required to name the stock
Note	

Number	5.3.7.3.2
Requirement	User will be required to select item’s collection
Note	Box with content type general: by default collection will be consumable

	Box with content type other than general: collection is same as box content type and uneditable
--	---

Number	5.3.7.3.3
Requirement	User will be required to select stock's content
Note	- Selection is possible only from the collection selected in 'Collection' field - User will be able to search by Item name/ SysID/ Catalogue number

5.3.8 Stocks customization page

5.3.8.1 Custom stocks type

Number	5.3.8.1.1
Requirement	User will be able to add custom stock type
Note	

Number	5.3.8.1.2
Requirement	User will not be able to add stock type if it already exists in the list
Note	

Number	5.3.8.1.3
Requirement	User will be able to remove stock type
Note	Possible for custom and default stock types, except for 'Other' stock type

Number	5.3.8.1.4
Requirement	User will not be able to remove stock type that is in use
Note	Either for active or consumed stock

5.3.9 Stock Pooling

5.3.9.1 Creating a new pooled stock

Number	5.3.9.1.1
Requirement	User will be required to add a unique pool name
Note	Each pool has a unique auto-generated ID

Number	5.3.9.1.2
Requirement	User will be able to pool by volume/weight
Note	Pooling by volume will be displayed by default

Number	5.3.9.1.3
Requirement	User will be able to draw by defined options of pooling logics: • Draw custom amount from each stock (default) • Draw fixed amount from all stocks • Draw the entire amount and consume
Note	Draw custom amount from each stock will be displayed by default

Number	5.3.9.1.4
Requirement	User will be able to define volume/weight draw amount from each selected stock
Note	• In 'Draw custom amount from each stock' option • By selecting volume/weight unit from drop-down list • By entering numeric input for each selected stock • User will not be able to pool more than the available amount in stock/s

Number	5.3.9.1.5
Requirement	User will be able to define a fixed, identical volume/weight amount to be drawn from all selected stocks
Note	• In 'Draw fixed amount from all stocks' option • By selecting volume/weight used unit • By entering numeric input for the new pooled stock • User will not be able to select more than the available amount in stock/s

Number	5.3.9.1.6
Requirement	User will be able to draw all the amount of selected stocks into a pool stock
Note	• In 'Draw the entire amount and consume' option • Amount will be calculated by volume/ weight • Volume/weight can be equaled to Zero • All selected stocks will be consumed and the amount by defined logic will be subtracted accordingly

Number	5.3.9.1.7
Requirement	User will be able to view the following fields of selected stock for pooling sample: • ID • Name • Storage location • Volume/Weight • Lot number

Confidential

	• Expiration date • Barcode • Volume used/ Weight used
Note	

Number	5.3.9.1.8
Requirement	User will be able to see the total amount for the new pooled stock
Note	For all 3 defined options of pooling logics: • Draw custom amount from each stock (default) • Draw fixed amount from all stocks • Draw the entire amount and consume

Number	5.3.9.1.9
Requirement	User will be required to name the created pooled stock
Note	

Number	5.3.9.1.10
Requirement	User will be required to select 'Collection' and 'Content' for pooled stock
Note	

Number	5.3.9.1.11
Requirement	User will be required to select location for pooled stock
Note	

Number	5.3.9.1.12
Requirement	User will be able to update pooled stock non required fields
Note	The calculated value of the pooled stock will carry over, and are not editable

Number	5.3.9.1.13
Requirement	User will be required to update pooled stock required custom collection fields
Note	

5.3.9.2 In stock show page

Number	5.3.9.2.1
Requirement	User will be able to view information related to the pooling event and the details of the current stock
Note	By clicking 'Pool' tab -> 'Parent stocks'

Number	5.3.9.2.2
Requirement	User will be able to view all information on stocks derived from the current stock
Note	By clicking 'Pool' tab -> 'Derived stocks'

5.3.10 Dispatch

5.3.10.1 Dispatch stocks

Applicable to stock index page, stock show page, box show page and collection show page (stocks tab)

Number	5.3.10.1.1
Requirement	User will be able to dispatch stock to another account
Note	- Dispatched stocks will be marked as consumed - Private stocks can be dispatched by stock owner or admin - User receives an email notification when the dispatch was completed

Number	5.3.10.1.2
Requirement	User will be required to select workspace to transfer the stock to
Note	- User can select from workspaces they are part of - The collection, including custom fields, needs to be the same in both accounts - All mandatory fields for the collection/stock item in the receiving account must be filled out - Stocks linked to collection items that do not exist in the receiving account will automatically create new items

Number	5.3.10.1.3
Requirement	User will be able to add notes
Note	

5.3.10.2 Dispatch events show page

Number	5.3.10.2.1
Requirement	User will be able to view dispatched items sent
Note	- By clicking on 'View dispatch events' - User will be able to enter stock page by clicking on its name

Number	5.3.10.2.2
Requirement	User will be able to view dispatched items received
Note	• By clicking on 'View dispatch events' • User will be able to enter stock page by clicking on its name

Number	5.3.10.2.3
Requirement	User will be able to assign the stock to himself/another member
Note	• Assigned date will be updated once user assigns the stock to a member • Dispatch stock will be placed into new general input box • User will be able to enter stock page by clicking on its name

Number	5.3.10.2.4
Requirement	User will be able to edit the member who receives the stock
Note	Assigned date will be updated

5.4 Labguru Plates module^(PLT)

5.4.1 Plates index page

Number	5.4.1.1
Requirement	User will be able to sort plates in ascending or descending order by any column header
Note	• By clicking on a header name • Not applicable for 'Size', 'Description' and 'Storage' columns

Number	5.4.1.2
Requirement	User will be able to search for plates by name
Note	

Number	5.4.1.3
Requirement	User will be able to tag plate/s
Note	By selecting plates and clicking on 'Tag'

Number	5.4.1.4
Requirement	User will be able to print labels for plates/s
Note	By selecting plates and clicking on 'Print labels'

Number	5.4.1.5
Requirement	User will be able to set the number of plates displayed per page
Note	Default number is 20 stocks per page. Can be changed to 50,100 or 200

5.4.2 Plate show page

Number	5.4.2.1
Requirement	User will see plate Id, plate owner, created at date, name of experiment it was created in, plate description, storage location and the plate layout display
Note	

Number	5.4.2.2
Requirement	User will be able to edit the plate
Note	

Number	5.4.2.3
Requirement	User will be able to add a new dataset from the plate page
Note	By clicking the 'Add new dataset' button

Number	5.4.2.4
Requirement	User will be able to print plate to a PDF file
Note	By clicking the 'Print' button

Number	5.4.2.5
Requirement	User will be able to print a label for plate
Note	By clicking the 'Print Label' button

Number	5.4.2.6
Requirement	User will be able to access the plate owner page
Note	By clicking on the owner's name at the owner field or on the bottom of the page

Number	5.4.2.7
Requirement	User will be able to access the experiment which is linked to the plate
Note	By clicking on the experiment name

Number	5.4.2.8
Requirement	User will be able to see the content (item and stock) of each well
Note	By hovering the 'i' icon that appears when hovering over the well

Number	5.4.2.9
Requirement	User will be able to access the page of item and stock that is in a well
Note	By clicking on their names inside the popup that opens from the 'i' icon on each well

Number	5.4.2.10
Requirement	User will be able to tag plate
Note	By clicking 'Add tag'

Number	5.4.2.11
Requirement	User will be able attach a file to plate
Note	

5.5 Labguru Equipment Module (EQU)

5.5.1 Equipment index page

Number	5.5.1.1
Requirement	User will be able to create a new equipment
Note	By clicking on the 'New Equipment' button

Number	5.5.1.2
Requirement	User will be able to sort equipment in ascending or descending order by any column header
Note	- By clicking on the column's title - By clicking on the kendo menu and select one of the sorting options - Sort preferences are saved

Number	5.5.1.3
Requirement	User will be able to search for equipment by name
Note	

Number	5.5.1.4
Requirement	User will be able to filter equipment by any combination of columns

Confidential

Note	Filter preferences are saved
------	------------------------------

Number	5.5.1.5
Requirement	User will be able to filter equipment by tags
Note	• By clicking on the tag's name • User can cancel selection by clicking on same tag again • Clicking on multiple tags will filter for stocks that are tags with all selected tags

Number	5.5.1.6
Requirement	Requirement removed
Note	

Number	5.5.1.7
Requirement	User will be able to tag equipment
Note	By selecting equipment and clicking on 'Tag'

Number	5.5.1.8
Requirement	User will be able to bulk delete selected equipment
Note	• By selecting equipment and clicking on 'Delete' • User will not be able to delete equipment if it was added to experiment/protocol • Double confirmation is required

Number	5.5.1.9
Requirement	User will be able to print labels for equipment
Note	By selecting equipment and clicking on 'Print labels'

Number	5.5.1.10
Requirement	User will be able to import/update equipment
Note	By clicking on cog menu -> 'Import/Update'

Number	5.5.1.11
Requirement	User will be able to view lists of maintenance events
Note	By clicking on the cog icon and 'Maintenance log'

Number	5.5.1.12
Requirement	User will be able to view list of archived equipment
Note	By clicking on 'View archived equipment'

Number	5.5.1.13
Requirement	User will be able to set the number of equipment displayed per page

Confidential

Note	Default number is 20 equipment per page. Can be changed to 50,100 or 200. User preferences will be saved
------	--

Number	5.5.1.14
Requirement	User will be able to see the next due maintenance date
Note	- In 'Maintenance date' column - If equipment has multiple maintenances, the date will be of the next from all

Number	5.5.1.15
Requirement	User will be clear all filters/sorts
Note	By clicking on the cog menu and 'Clear all filters/sorts'

Number	5.5.1.16
Requirement	User will be able to edit/bulk edit equipment
Note	- By clicking 'Edit' icon/ in-line edit in each row - Name field must be unique

Number	5.5.1.17
Requirement	User will be able to add and remove columns from the display
Note	- By clicking on the column menu -> 'Columns' and checking the wanted columns - Unchecking a column will remove it from view - User preferences will be saved

Number	5.5.1.18
Requirement	User will be able to adjust the width of the columns
Note	- By dragging and dropping the columns' borders - Auto-adjust by double clicking the border

Number	5.5.1.19
Requirement	User will be able to re-order columns
Note	- By dragging and dropping the columns - User preferences will be saves

Number	5.5.1.20
Requirement	User will be able to view previous import/update events
Note	By clicking on cog menu -> 'Previous Imports/Updates'

Number	5.5.1.21
Requirement	User will be able to export equipment
Note	- By clicking on cog menu -> 'Export' - Export all equipment not selecting any equipment - Export only selected equipment by first checking their boxes

	Export all filtered equipment by clicking on the header check box -> 'Select XXX equipment'
--	--

5.5.2 Equipment show page

Number	5.5.2.1
Requirement	User will be able to select a storage location for equipment
Note	

Number	5.5.2.2
Requirement	User will be able to edit equipment
Note	By clicking 'Edit'

Number	5.5.2.3
Requirement	User will be able to delete equipment
Note	- By clicking 'Delete' - Deletion is not allowed if equipment was added to experiment/protocol by using the Equipment element

Number	5.5.2.4
Requirement	User will be able to archive equipment
Note	By clicking 'Archive'

Number	5.5.2.5
Requirement	User will be able to print a label for equipment
Note	By clicking on 'Print label'

Number	5.5.2.6
Requirement	User will be able to book time on an equipment
Note	By clicking 'New Event' or by clicking on cells in the equipment calendar

Number	5.5.2.7
Requirement	User will be able to filter equipment events by other equipment
Note	By clicking on 'Filter' in the equipment calendar

Number	5.5.2.8
Requirement	User will be able to set maintenance to an equipment
Note	Via 'Maintenance' tab

Number	5.5.2.9
Requirement	User will be able to view list of attachments uploaded and connect to equipment by UpFolder
Note	By clicking on 'Output files' tab

Number	5.5.2.10
Requirement	User will be able to share equipment with other labs in a department
Note	If account is associated with a department account

Number	5.5.2.11
Requirement	User will be able to see list of experiments and protocols that the equipment was added to (via Equipment element)
Note	

Number	5.5.2.12
Requirement	User will be able to see list of equipment's usage records from experiments
Note	

5.5.3 Archived equipment

Number	5.5.3.1
Requirement	User will be able to restore archived equipment
Note	By clicking on 'Restore in archived equipment show page

Number	5.5.3.2
Requirement	User will not be able to add archived equipment to equipment element in experiments and protocols
Note	

5.5.4 Labguru Maintenance Module

5.5.4.1 Maintenance log page

Number	5.5.4.1.1
Requirement	User will be able to view maintenance events filtered by: • Performed by • Maintenance type • Equipment name • Equipment type • Maintenance data • Performed by
Note	

Confidential

Number	5.5.4.1.2
Requirement	User will be able to add new maintenance type
Note	

Number	5.5.4.1.3
Requirement	User will be able to remove maintenance type
Note	

Number	5.5.4.1.4
Requirement	User will be able to bulk add maintenance events for equipment of the same type
Note	

Number	5.5.4.1.5
Requirement	User will be able to delete maintenance event
Note	

Number	5.5.4.1.6
Requirement	User will be able to sort columns in ascending or descending order
Note	

Number	5.5.4.1.7
Requirement	User will be able to filter maintenance events by any combination of columns
Note	Filter preferences are saved

Number	5.5.4.1.8
Requirement	User will be able to search for equipment/maintenance by name
Note	

Number	5.5.4.1.9
Requirement	User will be able to adjust the width of the columns
Note	- By dragging and dropping the columns' borders - Auto-adjust by double clicking the border

Number	5.5.4.1.10
Requirement	User will be able to set the number of maintenance displayed per page
Note	Default number is 10 equipment per page. Can be changed to 5,20, 50 or 100. User preferences will be saved

Number	5.5.4.1.11
Requirement	User will be able to clear all filters
Note	

Number	5.5.4.1.12
Requirement	User will be able to clear all sorting
Note	

5.5.4.2 ‘Maintenance’ tab in equipment page

Number	5.5.4.2.1
Requirement	User will see past maintenance events added to equipment
Note	

Number	5.5.4.2.2
Requirement	User will be able to edit past maintenance events
Note	

Number	5.5.4.2.3
Requirement	User will be able to delete past maintenance event
Note	

Number	5.5.4.2.4
Requirement	User will be able to add new maintenance event
Note	By clicking on ‘Perform’

Number	5.5.4.2.5
Requirement	User will be able to add maintenance type to equipment
Note	

Number	5.5.4.2.6
Requirement	User will be required to choose frequency added maintenance type
Note	

Number	5.5.4.2.7
Requirement	User will be able to add maintenance type to all equipment of the same type
Note	

Number	5.5.4.2.8
Requirement	User will be able to add custom maintenance types

Note	By clicking 'Manage maintenance types' and then 'Add maintenance type'
------	--

Number	5.5.4.2.9
Requirement	User will be able to delete maintenance types
Note	By clicking 'Manage maintenance types' and delete

5.5.4.3 Adding new maintenance event

Number	5.5.4.3.1
Requirement	User will be able to select the member who performed the maintenance event
Note	By default, the user will be selected as event owner

Number	5.5.4.3.2
Requirement	When equipment has no type, user will be required to add equipment type when trying to add new maintenance event
Note	

Number	5.5.4.3.3
Requirement	When adding maintenance event for an equipment type for the first time, the user will be required to design new maintenance template
Note	In 'Design template' pop-up

Number	5.5.4.3.4
Requirement	User will be able to add maintenance data fields to maintenance template (name and type)
Note	Type of data fields: Text, numbers and date

5.5.4.4 Maintenance notifications

Number	5.5.4.4.1
Requirement	User will see in equipment element an indication on equipment maintenance status
Note	In experiment pages

Number	5.5.4.4.2
Requirement	User will see in dashboard notifications for equipment require maintenance
Note	

Confidential

Number	5.5.4.4.3
Requirement	User will see list of equipment that require maintenance
Note	By clicking on link in dashboard user will be directed to equipment index filtered to show only equipment require maintenance

Number	5.5.4.4.4
Requirement	User will see notification on equipment require maintenance in weekly digest email
Note	

Number	5.5.4.4.5
Requirement	User will see on the equipment show page notifications for each maintenance type that require maintenance
Note	

Number	5.5.4.4.6
Requirement	User will see in equipment index page an indication on equipment maintenance status
Note	

5.5.5 Labguru Scheduler Module(SCH)

5.5.5.1 Index page

Number	5.5.5.1.1
Requirement	User will be able to add a new schedule
Note	

Number	5.5.5.1.2
Requirement	User will be able to see schedules in a table list
Note	

Number	5.5.5.1.3
Requirement	User will be able to choose which columns to view on the index page
Note	By choosing headers from list in 'Columns' option

Number	5.5.5.1.4
Requirement	User will be able to search for schedules by name
Note	

Number	5.5.5.1.5
Requirement	User will be able to filter schedules by tags
Note	

Number	5.5.5.1.6
Requirement	User will be able to change the number of schedules he sees on single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 schedules per page

Number	5.5.5.1.7
Requirement	User will be able to filter schedules by any combination of columns in display
Note	

Number	5.5.5.1.8
Requirement	User will be able to sort schedules by any column in display in ascending/descending order
Note	

Number	5.5.5.1.9
Requirement	User will be able to edit the width and order of the columns
Note	

Number	5.5.5.1.10
Requirement	User will be able to see the name of the member that added the schedule, date of creation and date of last update
Note	- On the bottom of the page - The member's name is linked to the member's page

5.5.5.2 Schedule page

Number	5.5.5.2.1
Requirement	User will be able to edit schedule name
Note	

Number	5.5.5.2.2
Requirement	User will be able to add description to schedule
Note	

Number	5.5.5.2.3
Requirement	User will be able to add comments to schedule
Note	

Number	5.5.5.2.4
Requirement	User will be able to link schedule to other Labguru pages
Note	

Confidential

Number	5.5.5.2.5
Requirement	User will be able to duplicate schedule
Note	

Number	5.5.5.2.6
Requirement	User will be able to select instruments to schedule
Note	Added instruments will be listed below in the calendar view

Number	5.5.5.2.7
Requirement	User will be able to remove instruments from schedule
Note	When removing instrument from the list, all the events that were created in this schedule will be deleted

Number	5.5.5.2.8
Requirement	User will see all events of selected instruments in calendar view
Note	

Number	5.5.5.2.9
Requirement	User will be able to create events for instruments
Note	

Number	5.5.5.2.10
Requirement	User will be required to select event title, start date, end date and time
Note	

Number	5.5.5.2.11
Requirement	User will be able to choose color for schedule events
Note	

Number	5.5.5.2.12
Requirement	User will be able to delete an event from schedule
Note	Non admin user can delete own events

Number	5.5.5.2.13
Requirement	User will be able to delete schedule
Note	

Number	5.5.5.2.14
Requirement	User will be able to tag schedule
Note	

Number	5.5.5.2.15
Requirement	User will be able to add an attachment to schedule
Note	

Number	5.5.5.2.16
Requirement	User will be able to create a full day event
Note	By checking the 'All day' checkbox

Number	5.5.5.2.17
Requirement	User will be able to select between 'Timeline Month' or 'Timeline Day' view
Note	By default, user will see 'Timeline Day'

Number	5.5.5.2.18
Requirement	User will be able to see event details when hovering over the event in the timeline
Note	Event name, owner name, start and end (date&time) will be displayed

5.5.5.3 "Find next available date" wizard

Number	5.5.5.3.1
Requirement	User will be able to search for the next available date and time for a selected instrument
Note	Only for instruments that were added to schedule

Number	5.5.5.3.2
Requirement	User will be required to enter duration for the desired event (in hours or days)
Note	- Duration field will allow only integer numbers - By default, event duration will be '1' and 'hours'

Number	5.5.5.3.3
Requirement	User will be able to approve the suggested time and create an event
Note	

Number	5.5.5.3.4
Requirement	User will be able to reject the suggested time and go back to the previous window - to choose new instrument and duration
Note	

Number	5.5.5.3.5
Requirement	User will be able to search for the next available instrument from a selected equipment type

Note	- The option to select equipment type and the option to select instrument are mutually exclusive. - User will be able to choose any equipment type, regardless if the instruments that were added to schedule
------	--

Number	5.5.5.3.6
Requirement	Requirement removed
Note	

Number	5.5.5.3.7
Requirement	User will get as a result a suggestion to create an event for a specific instrument from the selected type
Note	The suggested event will be the first available where an instrument from the selected type does not have any other events for the selected period of time

Number	5.5.5.3.8
Requirement	User will see in the "Suggested event" pop-up the name of the specific instrument
Note	

Number	5.5.5.3.9
Requirement	User will be able to choose to "Find next" available event
Note	the search will re-start for all the instruments of the type, not only for the instrument that was suggested in the first event that was found

Number	5.5.5.3.10
Requirement	User will be able to create a new event in the schedule for the suggested instrument
Note	

Number	5.5.5.3.11
Requirement	User will be able to select date range to search in
Note	- Default search window is set to 20 weeks from today - User will be notified when selected a period of over 20 weeks

Number	5.5.5.3.12
Requirement	User will be able to search for the next available time after the suggested one

Note	By clicking on 'Next' in the suggested event pop up
------	---

5.5.5.4 Scheduler in Labguru calendars

Number	5.5.5.4.1
Requirement	User will see in calendar in equipment show page all events of schedules the equipment is part of
Note	

Number	5.5.5.4.2
Requirement	User will be able to view events of schedules in dashboard calendar
Note	

Number	5.5.5.4.3
Requirement	User will be able to go to schedule page, by clicking on an event of a schedule
Note	Either from equipment show page or from Dashboard calendar

5.5.6 Creating equipment events when concurrent events are not allowed

Number	5.5.6.1
Requirement	User will not be able to create equipment event at time the equipment has another event
Note	- When selectin start/end dates- unavailable time slots (days/hours) will be unselectable - When trying to create daily event and there are hourly events in range, user will be notified that event creation is not possible

Number	5.5.6.2
Requirement	User will be able to create equipment event at time where previously created event was deleted from
Note	

Number	5.5.6.3
Requirement	User will not be able to edit equipment event to time the equipment already has another event
Note	- When selectin start/end dates- unavailable time slots (days/hours) will be unselectable - When trying to edit daily event and there are hourly events in range, user will be notified that event editing is not possible

Number	5.5.6.4
Requirement	User will not be able to create equipment event with repeats if one or more repeats are at time the equipment already has another event
Note	When trying to create daily event with repeats and there are hourly events in range of one or more repeats, user will be notified that event creation is not possible

Number	5.5.6.5
Requirement	User will be able to add equipment event with repeats at time/s where a previously created events were deleted from
Note	

Number	5.5.6.6
Requirement	User will not be able to edit equipment event with repeats if one or more repeats are moved to time the equipment already has another event
Note	When trying to edit daily event with repeats and there are hourly events in range of one or more repeats, user will be notified that event editing is not possible

Number	5.5.6.7
Requirement	User will not be able to add an equipment event via scheduler at a time the equipment has another event
Note	

Number	5.5.6.8
Requirement	User will be able to create equipment event via scheduler at time where previously created event was deleted from
Note	

Number	5.5.6.9
Requirement	User will not be able to edit equipment event to time the equipment already has another event
Note	

5.5.7 Add/Edit Equipment page

Number	5.5.7.1
Requirement	User will be required to name the equipment
Note	

Number	5.5.7.2
Requirement	User will be required to add equipment owner

Confidential

Note	- As default, equipment owner will be the member who added the equipment - By clicking on ‘Owner’ search option or selecting from members list
------	---

Number	5.5.7.3
Requirement	User will be able to add manufacturer
Note	- Creating a new manufacture will add it in the companies page - If manufacturer already exist it will add it to the company

Number	5.5.7.4
Requirement	User will be able to add serial number
Note	

Number	5.5.7.5
Requirement	User will be able to add equipment type
Note	

Number	5.5.7.6
Requirement	User will be able to add a unique barcode
Note	

Number	5.5.7.7
Requirement	User will be able to add model number
Note	

Number	5.5.7.8
Requirement	User will be able to add purchase date
Note	

Number	5.5.7.9
Requirement	User will be able to add warranty expiration date
Note	

Number	5.5.7.10
Requirement	User will be able to add instruction SOP
Note	

Number	5.5.7.11
Requirement	User will be able to add description

Note	
------	--

Number	5.5.7.12
Requirement	User will be able to add usage rate
Note	User will be required to add the following information: • Rate amount (from zero and higher) • Currency • Type of rate (rate per hour/rate per use)

6. System functions

6.1 Labguru Linked Resources module^(LRS)

Number	6.1.1
Requirement	User will be able to link to a page that was recently viewed or edited
Note	

Number	6.1.2
Requirement	User will be able to search for page to link by name or SysID (only for inventory items)
Note	• User will see creation date, type and name of each search result • Results list is limited to 80 results, user will be able to see more by clicking on “see more” button • User will be able to multiple object from the search results list

Number	6.1.3
Requirement	User will be able to view a visual diagram of network of links around the page
Note	• By clicking on the compass icon in ‘Linked resources’ box • Network will show page maximum of 2 links far from page

Number	6.1.4
Requirement	User will be able to remove a link
Note	By clicking on the ‘x’ sign

Number	6.1.5
Requirement	User will be able to delete all linked resources
Note	By clicking on the ‘Remove all links’ in the ‘Linked Resources’ box

6.2 Labguru Tasks module(TSK)

6.2.1 Tasks page

Number	6.2.1.1
Requirement	User will see open tasks assigned to him divide into 2 lists: 'Overdue tasks' and 'Open tasks'
Note	For each overdue/open task, user will see the task's name, who assigned it and when and its due date

Number	6.2.1.2
Requirement	User will see the number of overdue and open tasks assigned to him
Note	The total number of overdue and open tasks appear next to each title

Number	6.2.1.3
Requirement	User will be able to expand and collapse lists of overdue/open tasks
Note	

Number	6.2.1.4
Requirement	User will be able to go to task page by clicking on its name or its ID
Note	

Number	6.2.1.5
Requirement	User will be able to go to the page the task was created from by clicking on the value under 'Related to'
Note	

Number	6.2.1.6
Requirement	User will be able to mark task assigned to him as completed
Note	- By clicking on 'Mark as completed' in overdue/open lists - To mark multiple tasks as completed, user can select tasks' checkboxes and click on 'Mark as completed'

Number	6.2.1.7
Requirement	User will be able to tag multiple tasks
Note	By selecting tasks and clicking on the 'Tag' icon

Number	6.2.1.8
Requirement	User will be able to delete multiple tasks

Confidential

Note	• By selecting tasks and clicking on the ‘Delete’ icon • User will be able to delete only tasks owned by him
------	---

Number	6.2.1.9
Requirement	User will be able to add a new task
Note	

Number	6.2.1.10
Requirement	User will be able to filter tasks by tags
Note	

Number	6.2.1.11
Requirement	User will be able to search tasks by name or description
Note	

Number	6.2.1.12
Requirement	User will be able to change the number of tasks he sees on single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 tasks per page

Number	6.2.1.13
Requirement	User will be able to sort tasks by any column in display in ascending/descending order
Note	

Number	6.2.1.14
Requirement	User will be able to filter tasks by any combination of columns in display
Note	

Number	6.2.1.15
Requirement	User will be able to choose which tasks to view on the index page
Note	• By choosing headers from list in ‘Columns’ option • User's preference will be saved

Number	6.2.1.16
Requirement	User will be able to change the width of the columns
Note	By dragging and dropping the columns' borders

6.2.2 From inventory item show page

Number	6.2.2.1
Requirement	User will be able to create a new task
Note	

Number	6.2.2.2
Requirement	User will be able to add: - task name - due date - description
Note	
Number	6.2.2.3
Requirement	User will be able to assign the task to himself or to another member
Note	

Number	6.2.2.4
Requirement	User will be able to choose to whom the email notification will be sent to
Note	

Number	6.2.2.5
Requirement	User will be able to set the date time of the email notification to be sent
Note	

Number	6.2.2.6
Requirement	User will be able to mark task assigned to him as completed
Note	

Number	6.2.2.7
Requirement	User will be able to delete tasks assigned to him
Note	

6.2.3 From task show page

Number	6.2.3.1
Requirement	Requirement removed
Note	

Number	6.2.3.2
Requirement	User will be able to edit a task

Note	
------	--

Number	6.2.3.3
Requirement	User will be able to mark task as completed
Note	

Number	6.2.3.4
Requirement	User will be able to comment on task
Note	

Number	6.2.3.5
Requirement	User will be able to add attachments (images & files) to task
Note	

Number	6.2.3.6
Requirement	User will be able to delete a task
Note	

Number	6.2.3.7
Requirement	User will be able to tag a task
Note	

6.3 Labguru Tags module^(TAG)

6.3.1 From item show page

Number	6.3.1.1
Requirement	Requirement removed
Note	

Number	6.3.1.2
Requirement	User will be able to remove a tag from any show page
Note	

6.3.2 From item index page

Number	6.3.2.1
Requirement	User will be able to batch-tag items via item's index page
Note	Refers to inventory items, storage & equipment items (except 'Storages'), datasets

Number	6.3.2.2
Requirement	Requirement Removed

Note	
Number	6.3.2.3
Requirement	Requirement Removed
Note	

6.3.3 From tags Index page

Number	6.3.3.1
Requirement	Tags will appear sorted alphabetically in: • The tags panel of inventory item index page • The Tags index page
Note	

Number	6.3.3.2
Requirement	User will be able to • Merge • Edit • Delete Tags
Note	

6.4 Labguru Comments module (CMN)

Number	6.4.1
Requirement	Requirement removed
Note	

Number	6.4.2
Requirement	User will be able to select members that will be notified when the comment is posted
Note	By clicking on the ‘Select participants’ select box and choosing members.

Number	6.4.3
Requirement	User name, picture and his member id will be displayed once the comment is posted along with date & time
Note	Post counter will be updated with the updated number of comments in the page.

Number	6.4.4
Requirement	User will be able to access the page the comment is related to from the email notification
Note	• By clicking on ‘Reply in Labguru’

	• If user is not logged in to Labguru, user will need to login and then be directed to the correct page • If member in multiple account, user will be directed to the account the comment was sent from
--	--

6.5 Labguru Notes module^(NTE)

6.5.1 Notes in projects and folders

Number	6.5.1.1
Requirement	User will be able to add a note to project and folder
Note	• By clicking on 'Add note' in 'Notes' tab • Project/folder name will be filled automatically with the project/folder name and will be greyed out

Number	6.5.1.2
Requirement	User will be required to add note title
Note	

Number	6.5.1.3
Requirement	User will be able to view a list of notes with the following details: • Title of the note • Date of when the note was added • Member picture • Number of comments to note have • Number of attachments the note has • Note preview
Note	In the 'Notes' tab

6.5.2 Note show page

Number	6.5.2.1
Requirement	User will be able to edit note
Note	By clicking on 'Edit'

Number	6.5.2.2
Requirement	User will be able to delete a note
Note	By clicking on 'Delete'

Number	6.5.2.3
Requirement	User will be able to comment on a note
Note	

6.6 Global search^(GLB)

Number	6.6.1
Requirement	User will be able to start a search from any page
Note	• Available from top panel • User will get the search results ordered by relevance • search term will be highlighted in each result

Number	6.6.2
Requirement	User will see results in pages divided to 20 results in a page
Note	
Number	6.6.3
Requirement	User will be able to click on the name of the result and be directed to its page
Note	For comment and attachments, User will be directed to the page where they were added

Number	6.6.4
Requirement	User will see the total number of matched found
Note	

Number	6.6.5
Requirement	User will be able to click on a category and to filter the page to show only the results belong to that category
Note	

Number	6.6.6
Requirement	User will not see categories which have no search results
Note	

Number	6.6.7
Requirement	The user will see the search term he used at the top of the page
Note	Search term will be shown in bold.

Number	6.6.8
Requirement	User will be able to filter search results by update date
Note	

Number	6.6.9
Requirement	User will be able to filter search results by owner
Note	

Number	6.6.10
Requirement	User will be able to perform exact searches
Note	By using single (') or double (") quotes

6.7 'Recently viewed' list^(RCL)

Number	6.7.1
Requirement	User will be able to see his recently viewed items in the 'recently viewed' list
Note	

Number	6.7.2
Requirement	User will be able to pin items to the 'recently viewed' list
Note	

Number	6.7.3
Requirement	Pinned item name will be updated in the 'recently viewed' list
Note	

Number	6.7.4
Requirement	Requirement removed
Note	

6.8 Labguru Events module^(EVT)

6.8.1 Calendar Page

Number	6.8.1.1
Requirement	User will be able to create a new event: - By clicking on the 'new event' button via 'Calendar tab' - By clicking on a calendar date
Note	An 'Add new event' popup will be opened User will be required to enter event start date & event end date

Number	6.8.1.2
Requirement	By default, user will see only events that he is a participant or owned by him
Note	Includes general events, equipment events, folders due date, time range elements and tasks

Confidential

Number	6.8.1.3
Requirement	Requirement Removed
Note	

Number	6.8.1.4
Requirement	User will be able to include/exclude events from the calendar view by the following categories: • Labguru calendar • Equipment calendar • Google calendar
Note	User selection will be saved after leaving the page

Number	6.8.1.5
Requirement	User will be able to include/exclude events owned to former events
Note	

Number	6.8.1.6
Requirement	User will be able to include/exclude other members' calendars to the calendar view
Note	

Number	6.8.1.7
Requirement	User will be able to filter the events displayed in calendar to show only: • Experiments • Folders • Tasks • Events
Note	

Number	6.8.1.8
Requirement	User will be able to view a specific day
Note	By selecting the day from the mini calendar

Number	6.8.1.9
Requirement	User will be able to view today's calendar
Note	By clicking on 'Today'

Number	6.8.1.10
Requirement	User will be able to select between month, week or day view
Note	

Number	6.8.1.11
Requirement	User will be able custom calendar settings: • Time format - 24H/12H • Default view - Month/Week/Day • Week starts on - Sunday/Monday
Note	

Number	6.8.1.12
Requirement	User will be able to jump one day forward or backwards
Note	

Number	6.8.1.13
Requirement	User will be able to delete own events
Note	By hover over the event in the calendar and click on 'Delete'

6.8.2 Adding/editing events

Number	6.8.2.1
Requirement	User will be required to enter event title
Note	

Number	6.8.2.2
Requirement	User will be able to add event description
Note	Description will be show displayed when user hovers over the event item in the calendar

Number	6.8.2.3
Requirement	User will be able to select new start and end date and time
Note	The default will be the date selected 9 until 18

Number	6.8.2.4
Requirement	User will not be able to select event's end date that is before the event's start date
Note	

Number	6.8.2.5
Requirement	User will be able to create a full day event by checking the 'All day' checkbox
Note	For full day events, the user will not be able to pick time

Number	6.8.2.6
Requirement	User will be able to make the event a repeatable event

Confidential

Note	By clicking on 'Repeats' and selecting either daily, weekly or monthly repeats
------	--

Number	6.8.2.7
Requirement	User will be able to select members as participants in event
Note	- By clicking on 'Select members' - Multiple members and teams can be selected - By default, all members are participants - The event will appear in participants' calendar

Number	6.8.2.8
Requirement	User will be able to add all lab members as participants in the event by clicking on 'all lab members'
Note	

Number	6.8.2.9
Requirement	User will be able to select who will receive email notification
Note	- Selection is between - only the user, all participants, or no notification - The default option will be notification only to user

Number	6.8.2.10
Requirement	User will be able to select date and time for email notification
Note	Notification date must be before event start date

6.8.3 Event show page

Number	6.8.3.1
Requirement	User will be able to edit event
Note	By clicking on the 'Edit' link

Number	6.8.3.2
Requirement	User will be able to delete own event
Note	By clicking on the 'delete' link

Number	6.8.3.3
Requirement	User will be able to link event to other pages
Note	By clicking on the 'add link' link in the 'Linked resources' box

Number	6.8.3.4
Requirement	User will be able to attach files to event
Note	By clicking on the 'Upload a file' box

Number	6.8.3.5
Requirement	User will be able to add comments to event
Note	By typing in the 'Discussion box' and clicking 'Post'

6.9 Section removed

6.10 Print to PDF

Number	6.10.1
Requirement	User will be able to print page sections
Note	- By selecting 'Sections,' the user can choose either all sections or specific elements. - By default, all sections are selected

Number	6.10.2
Requirement	User will be able to print metadata
Note	- By selecting 'Metadata,' the user can choose either all metadata options or specific elements - By default, all metadata options are selected

Number	6.10.3
Requirement	User will be able to print comments
Note	

Number	6.10.4
Requirement	User will be able to add header & footer to the PDF
Note	- By marking 'Set header&footer' - User can choose elements from a predefined list and position them on the left, right, or center of the page

Number	6.10.5
Requirement	User will be able to fit content to page
Note	

Number	6.10.6
Requirement	User will be able to print a visual layout of a plate element
Note	By selecting 'Include plate table view' when a plate element is in the page

Number	6.10.7
Requirement	User will be able to view list of added to page attachments in the PDF

Note	
------	--

6.11 Label Printing^(LBP)

Applicable to the following sections: Stocks, Boxes, Equipment and plates

6.11.1 General

Number	6.11.1.1
Requirement	User will be able to add to the label: - Barcode, QR code or Datamatrix - Company logo
Note	- Barcodes will be generated automatically - When a QR code is scanned, URL will be displayed

Number	6.11.1.2
Requirement	User will be able to select which label template to use from existing templates
Note	Number available templates is displayed on the top of the window

Number	6.11.1.3
Requirement	User will be able to print label/s to a PDF file
Note	

Number	6.11.1.4
Requirement	Admin user will be able to delete label template
Note	

Number	6.11.1.5
Requirement	Admin user will be able to edit label template
Note	

Number	6.11.1.6
Requirement	Admin user will be able to add new template
Note	

6.11.2 Label templates settings

Number	6.11.2.1
Requirement	User will be able to edit/enter template name
Note	

Number	6.11.2.2
Requirement	User will be able to set the label dimensions
Note	- Width and height in millimeters - Logo and barcode dimensions change accordingly

Number	6.11.2.3
Requirement	User will be able to add new customize line/s
Note	Able to add type of information from pre-defined list and customize title

Number	6.11.2.4
Requirement	User will be able to remove information line/s
Note	

Number	6.11.2.5
Requirement	User will be able to select item attributes from pre- defined list
Note	Name field will be named same as type information by default

Number	6.11.2.6
Requirement	User will be able to edit attribute title name
Note	

Number	6.11.2.7
Requirement	User will be able to remove selected item attributes
Note	

Number	6.11.2.8
Requirement	User will be able to access ‘Advanced settings’
Note	

6.11.3 Template advanced settings

Number	6.11.3.1
Requirement	User will be able to customize general properties of the lines: - Font size - Line margins (mm) - Line height (mm)
Note	By clicking ‘Advances setup’

Number	6.11.3.2
Requirement	User will be able to customize logo layout: - Position (vertical and horizontal) from pre-defined list - Width (mm)

Note	By selecting: Top/Middle/Bottom or Custom (in mm)
------	---

Number	6.11.3.3
Requirement	User will be able to customize barcode type layout: - Position (vertical and horizontal) from pre- defined list - Width (mm)
Note	By selecting: Top/Middle/Bottom/Custom (in mm)

Number	6.11.3.4
Requirement	User will be able to exit advanced mode
Note	

6.12 System notifications

Number	6.12.1
Requirement	User will be notified when leaving an unsaved page
Note	- For the following pages: Project, Folder, Experiment, Protocol, SOP, Report, Document, Recipe - User will be able to choose ‘Leave Page’ or Stay on Page’

Number	6.12.2
Requirement	User will be notified when the internet connection goes offline/online
Note	

7. **Signing and Witnessing** (SIG)

7.1 Signed page

(applicable to experiments, protocols, SOPs, reports and documents)

Number	7.1.1
Requirement	The signer will be able to revert signature of a signed page
Note	

Number	7.1.2
Requirement	User will not be able to add/edit/delete elements in a signed page
Note	

Number	7.1.3
Requirement	User will not be able to delete a signed page
Note	

Confidential

Number	7.1.4
Requirement	Requirement removed
Note	

Number	7.1.5
Requirement	User will not be able to add attachment to signed page
Note	

Number	7.1.6
Requirement	User will be able to archive a signed page
Note	Not applicable to experiments

Number	7.1.7
Requirement	User will be able to print signed page to pdf
Note	

Number	7.1.8
Requirement	User will be able to duplicate a signed page
Note	Signature will not be duplicated

Number	7.1.9
Requirement	User will be able to view signed page version history
Note	

Number	7.1.10
Requirement	User will be able to add a comment to a signed page
Note	

Number	7.1.11
Requirement	User will see notification of signature at top of the page
Note	Notification will include signer name, id and signature time stamp

Number	7.1.12
Requirement	User will be able to add a comment when signing or reverting a page
Note	Comment will be added to discussion section at the bottom of the page with name and id of the signer, date of signing/reverting/rejecting/witnessing

Number	7.1.13
Requirement	User will be able to select a member to receive the signing/reverting comment as email
Note	The recipient will be able to click on a link in email to get to page

Number	7.1.14
Requirement	User will be able to update list of members/teams who can view signed page
Note	Applicable only for protocols

7.2 Signed experiment

Number	7.2.1
Requirement	User will be able to save signed experiment as a protocol
Note	

Number	7.2.2
Requirement	User will not be able to flag/unflag signed experiment
Note	

Number	7.2.3
Requirement	User will be able to move signed experiment to other folders
Note	

7.3 Signed protocol and SOP

Number	7.3.1
Requirement	User will be able to start an experiment from signed protocol
Note	

Number	7.3.2
Requirement	User will be able to start experiments from signed SOP
Note	

7.4 Witness a signed page

(applicable to experiments, protocols, SOPs, reports, documents and request templates)

Number	7.4.1
Requirement	Authorized user will be able to witness or reject a signed page
Note	Witness must not be the signer

Confidential

	- User with “Can witness” privilege will be able to witness experiments that are part of a project with default signature flow - User who is listed as witness in a custom signature flow will be able to witness experiments that are part of a project with that signature flow
--	--

Number	7.4.2
Requirement	User with ‘Can witness’ privilege will not be able to witness or reject an un-signed page
Note	

Number	7.4.3
Requirement	User will see notification of witnessing at top of the page
Note	Notification will include witness name, id and witnessing time stamp

Number	7.4.4
Requirement	Removed requirement
Note	

Number	7.4.5
Requirement	User ‘Can witness’ privilege will be able to add a comment when witnessing or rejecting a signed page
Note	Comment will be added to discussion section at the bottom of the page with name and id of the witness, date of witnessing/rejecting

Number	7.4.6
Requirement	User ‘Can witness’ privilege will be able to select a member to receive the witnessing/rejecting comment as email
Note	- The recipient will be able to click on a link in email to get to page - By default, witnessing/rejecting comment will be sent to signer

Number	7.4.7
Requirement	User will not be able to revert witnessed page
Note	

Number	7.4.8
Requirement	After witnessing, user should see the creation of three files: - pdf of experiment - pdf of history report

	• Zip file that includes the pdf and attachments in original formats
Note	

7.5 User authentication in signature actions

Number	7.5.1
Requirement	User will be required to enter email address and password in order to sign/revert/witness/reject a page
Note	

Number	7.5.2
Requirement	User in an environment with Single Sign On will not be required to enter email address and password in order to sign/revert/witness/reject a page
Note	

8. Account and user management

8.1 Labguru Billing module^(BIL)

8.1.1 Index page

Number	8.1.1.1
Requirement	User will be able to view its billing info
Note	• By clicking on the account drop down menu and clicking the 'Billing' link • User will be able to see the date of his next payment will be built on, date of his next charge, and his payment method

Number	8.1.1.2
Requirement	User will be able to set payment method
Note	• By clicking the 'Set Payment Method' link • 'Set Payment method' page will be opened • User will be required to fill his first name, last name, card number, expiration date & CVV number

Number	8.1.1.3
Requirement	User will be able to get a quote
Note	• By clicking the 'Get quote' button • 'Renewal quote request' pop up will be opened • User will be requested to fill the renewal period(years) & number of seats • User will be required to fill his email and send the request

Confidential

Number	8.1.1.4
Requirement	User will be able to cancel his account
Note	• By checking the 'I understand this action is irreversible' checkbox • A 'Please cancel my account' link will show up • Delete confirmation page will be open, explaining to User that this action can't be undone • A 'I understand- delete my account' button will be shown with 'cancel' button right next to it

Number	8.1.1.5
Requirement	User will be redirected to a billing page if trial period has expired
Note	

Number	8.1.1.6
Requirement	User will be charged when entering CC info to an expired trial for the current period
Note	

Number	8.1.1.7
Requirement	Account owner/Admin will receive red notification 7 days prior to upcoming expiry of subscription
Note	• User can remove the message • Notification will be removed/not enabled if payment has been established

Number	8.1.1.8
Requirement	Account owner/Admin will receive yellow notification 7-30 days prior to upcoming expiry of subscription
Note	• User can remove the message • Notification will be removed/not enabled if payment has been established

Number	8.1.1.9
Requirement	Account owner/Admin will receive red notification when starting 7 days grace period after the end of subscription
Note	• User can remove the message • Notification will be removed/not enabled if payment has been established

8.1.2 External credit card link

Number	8.1.2.1
Requirement	User will be redirected to an external CC from account controller to Billing Controller
Note	

Number	8.1.2.2
Requirement	User will be redirected back to billing page in case the credit card is rejected
Note	

8.2 Account settings module^(ACS)8.2.1 General settings

Number	8.2.1.1
Requirement	User will be able to set - Account name - Time zone
Note	

Number	8.2.1.2
Requirement	User will be able to see the amount of storage space he has in the account
Note	

Number	8.2.1.3
Requirement	User will be able to check/uncheck the 'Automatically mark stocks as consumed' check box
Note	By checking this check box User agrees that Stocks will be automatically marked as consumed when volume or weight values are zero

Number	8.2.1.4
Requirement	User will be able to turn on an option that prevents the creation of concurrent equipment events
Note	- Will affect only equipment events - User will be notified that if part of a department, this function must be activated in all department's accounts for proper functionality"

Number	8.2.1.5
Requirement	User will be able to enable support access
Note	By entering date, access to user environment is enabled

Number	8.2.1.6
Requirement	User will be able to set date format from dropdown list
Note	• The default format is: YYYY-MM-DD • Applicable only for print PDF file

Number	8.2.1.7
Requirement	User will be able to set time format from dropdown list
Note	• The default format is: HH:mm • Only for print PDF file/ print labels

8.2.2 Roles management

Number	8.2.2.1
Requirement	Requirement removed
Note	

Number	8.2.2.2
Requirement	User will be able to view list of roles, the privileges given to each role, and the number of members assigned to each role
Note	• Only users with 'Can admin system' privilege will be able to view the page • From the 'Settings' page in the account's drop-down menu

Number	8.2.2.3
Requirement	Requirement removed
Note	

Number	8.2.2.4
Requirement	User will be able to check or un-check a privilege for each role separately
Note	• For privileges: 'System Admin', 'Read Only', 'Can approve orders', 'Can witness' and 'Can Access Workflow Editor and Dashboards' • When checking System Admin - all other privileges are updated to match admin permission and become un-editable • When checking Read-only - all other privileges are updated to match Read-only permission and become un-editable • Un-checking System Admin or Read-only bring all privileges to default setting and they become editable

Number	8.2.2.5
Requirement	User will be able to set Inventory and Equipment to one of the following: None, View, Edit, Delete or Admin
Note	- Selection in one module (Inventory or Equipment) is independent on the selection made in the other module - Privileges in one module are mutual exclusive

Number	8.2.2.6
Requirement	Requirement removed
Note	

Number	8.2.2.7
Requirement	Requirement removed
Note	

Number	8.2.2.8
Requirement	User will be able to remove un-used roles from the list
Note	- User will be notified when trying to remove a used role ‘Other’ role cannot be deleted

Number	8.2.2.9
Requirement	User will be able to add custom roles
Note	Default settings for a new role are ‘System Admin’, ‘Ready-only’, ‘Can witness’, ‘Can approve orders’ and ‘Can Access Workflow Editor and Dashboards’ are turned off. Inventory and Equipment permissions are set to ‘Edit’

Number	8.2.2.10
Requirement	User will be able to edit role’s name
Note	

Number	8.2.2.11
Requirement	Requirement removed
Note	

Number	8.2.2.12
Requirement	User will be able to see list of members assigned to a role
Note	By clicking on the number of members next to role name

Number	8.2.2.13
Requirement	User will be able to assign a member to a role
Note	By hovering on role name and clicking on ‘Assign’

- | | |
|--|---|
| | <ul style="list-style-type: none"> User will be required to enter password |
|--|---|

8.2.3 New collections and collections settings

Number	8.2.3.1
Requirement	User will be able to add new collection
Note	- By filling the 'collection name' box - By clicking the 'Add new' button

Number	8.2.3.2
Requirement	User will be able to go to a collection customization page
Note	By clicking the 'customize' link next to the relevant collection

Number	8.2.3.3
Requirement	User will be able to select one or more of the pre-defined model organisms
Note	- By checking/unchecking the collections checkboxes - Unchecked collections will not be shown in the inventory collections list

Number	8.2.3.4
Requirement	User will be required to set a unique name to new collections
Note	

Number	8.2.3.5
Requirement	User will be able to delete custom collection that has no items
Note	

8.2.4 Data export

Number	8.2.4.1
Requirement	User will be able to export his account
Note	- By clicking the 'Export your account' button 'Your export is being generated' message will appear on screen - User will be notified by mail within 24 hours, when the export is done - Export files will appear in 'Data Export' tab, in ZIP format, portioned to 1gb files and will be valid up to 72 hours - User can export the account once a day

8.2.5 Account personalization

Number	8.2.5.1
Requirement	User will be able upload an account logo
Note	By clicking on the 'Click to upload' link

Number	8.2.5.2
Requirement	User will be able to go to choose color scheme
Note	

8.2.6 Signature process

Number	8.2.6.1
Requirement	User will be able to view list of custom signature processes
Note	- Only users with 'Can admin system' privilege will be able to view the page - From the 'Settings' page in the account's drop-down menu - By clicking flow name, created signature flow is displayed

Number	8.2.6.2
Requirement	User will be able to create a new signature flow
Note	User will be required to enter name

Number	8.2.6.3
Requirement	User will be able to select sign logic for each step in signature flow
Note	- Selection between pre-defined list of members/teams or all members that are part of the project - User can sign only once

Number	8.2.6.4
Requirement	User will be able to edit signature flow
Note	Signature flows already assigned to experiments are uneditable

Number	8.2.6.5
Requirement	User will be able to delete a signature flow
Note	- By clicking on the trash can icon next to title of the flow - Double confirmation will be required in order to delete a flow - Signature flows already assigned to experiments cannot be deleted

Number	8.2.6.6
Requirement	User will be able to add steps to signature flow
Note	Steps can be signed only by their set order

Number	8.2.6.7
Requirement	User will be able to re-arrange steps order
Note	

Number	8.2.6.8
Requirement	User will be able to delete steps
Note	Witnessing step cannot be deleted

Number	8.2.6.9
Requirement	User will be required to select role of signer
Note	Step name will be derived from role of signers

Number	8.2.6.10
Requirement	User will be able to create new role of signer
Note	New custom role can be used in the created flow and will not be saved in the list

Number	8.2.6.11
Requirement	User will be able to set number of signatures required to proceed to next step
Note	By default, 1 signature, maximum 3 signatures

8.2.7 Account LOVs

Number	8.2.7.1
Requirement	User will be able to see list of account LOV's
Note	

Number	8.2.7.2
Requirement	User will be able to add new account LOV
Note	- By selecting name, description and list of values for each account LOV - Name field is mandatory - User is required to add at least one value - Multiple values, will be displayed in alphabetically/ascending order

Number	8.2.7.3
Requirement	User will be able to delete account LOV, if it is not used anywhere

Note	
------	--

Number	8.2.7.4
Requirement	User will be able to edit account LOV's name and description
Note	

Number	8.2.7.5
Requirement	User will be able to add value to existing account LOV
Note	

Number	8.2.7.6
Requirement	User will be able to inactivate a value in existing account LOV
Note	- Inactive values will become unavailable to be added to new entities - Inactive values will still appear in any entity that was added before removal - User will be able to see inactive values for each LOV

Number	8.2.7.7
Requirement	User will be able to add values previously inactivated values
Note	

Number	8.2.7.8
Requirement	User will be able to create LOV with multiple select values
Note	By marking 'Allow selecting multiple values'

Number	8.2.7.9
Requirement	User will not be able to change existing LOV from single select to multiple select and vice versa
Note	

8.2.8 Forms Settings (Reason for Change)

Number	8.2.8.1
Requirement	Admin user will be able to add 'Require reason for change'
Note	By clicking on 'Advanced' tab and then on 'Form Settings'

Number	8.2.8.2
Requirement	Admin user will be able to select which reason/s will be available to the users
Note	By selecting from a pre-existing list

Number	8.2.8.3
--------	---------

Requirement	Admin user will be able to add new reason/s for change
Note	

Number	8.2.8.4
Requirement	Admin user will be able to delete reasons
Note	

Number	8.2.8.5
Requirement	Admin user will be able to enable/disable free text insertion as reason for changes
Note	

8.2.9 PDF Settings

Number	8.2.9.1
Requirement	Admin user will be able to add header/footer to PDF printouts
Note	• By clicking on 'Advanced' tab and then 'PDF Settings' • By activating pre-define option/s: experiments, protocols, SOPs, reports and/or documents • User will be able to select which header/footer to print from drop down menu: name/ creation date/ owner name/page number/confidential/signed ny and date • Header/ footer location options: left/center/right

8.2.10 Experiments Settings

Number	8.2.10.1
Requirement	Admin user will be able to see list of experiments flags
Note	By clicking on 'Advanced' tab and then 'Experiments flags'

Number	8.2.10.2
Requirement	Admin user will be able to add new experiments flags
Note	• By selecting name, description, color, icon and status • Name field is mandatory (limited to 40 characters) • Default status is 'Active'

Number	8.2.10.3
Requirement	Admin user will be able to change the flag status
Note	By default status will be active

Number	8.2.10.4
Requirement	Admin user will be able to edit experiments flags description, color and icon

Note	- User will not be able to edit flag name - For the two default flags 'successful' and 'failed' the user will be able to edit only the status field
------	--

Number	8.2.10.5
Requirement	Admin user will be able to delete experiments flags
Note	- Only if it is not used in an experiment 'Successful' and 'failed' flags cannot be deleted

8.2.11 Requests Settings

Number	8.2.11.1
Requirement	Admin user will be able to add custom request types
Note	- By clicking on 'Advanced' tab and then 'Requests settings' - Request type name must be unique

Number	8.2.11.2
Requirement	Admin user will be able to edit custom request types
Note	- Request type name must be unique - Default request types: stocks, test and shipment cannot be edited

Number	8.2.11.3
Requirement	Admin user will be able to delete custom request types
Note	Default requests types: stocks, test, shipment cannot be deleted

8.3 Privileges in Labguru^(PRV)

Number	8.3.1
Requirement	User will be able to perform the actions as specified in Users Permissions Table (Appendix A)
Note	

8.4 Labguru Login/Logout module

8.4.1 Labguru Login^(SGN)

Number	8.4.1.1
Requirement	User will be required to enter email address and password to sign in

Confidential

Note	• By adding his email address, clicking the 'Next' button, adding his password, and clicking the 'Sign in' button • Clicking the 'Sign in' button changes the button to a 'Signing in' button • After signing in, User is being directed to the dashboard page
------	--

Number	8.4.1.2
Requirement	User will be able to reset his password in case it was forgotten
Note	By clicking the 'Forgot my password' button (see separate requirements)

Number	8.4.1.3
Requirement	User will be limited to 3 failed logins before being blocked
Note	• Blocked users will receive an email with link to account unblocking page • Blocked user will be unlocked automatically after 1 hour

Number	8.4.1.4
Requirement	User will be able to switch between account
Note	• By clicking on the account name and selecting the account to switch to from the dropdown list • Relevant only for users that are part of multiple accounts • User will be notified when switching to a different account in another tab

Number	8.4.1.5
Requirement	User will be required to enter 2FA code while logging in
Note	

Number	8.4.1.6
Requirement	User will be required to enable 2FA in the first login or up to 30 days from login
Note	• By clicking on 'Enable now' or 'Remind me later' on profile setting page • 2FA can be enabled by scanning QR code with authentication app

Number	8.4.1.7
Requirement	User will be able to see Labguru version number
Note	

8.4.2 Labguru First Login(SGN)

Number	8.4.2.1
Requirement	When invited to Labguru account, user will be required to accept invitation
Note	In organizations with SSO, user does not need to accept invite or set password when invited to account

Number	8.4.2.2
Requirement	User will be required to set new password
Note	See req. 8.5.1.5 for requirements on password strength

Number	8.4.2.3
Requirement	User will be required to confirm password
Note	

8.4.3 Labguru Logout(SGN)

Number	8.4.3.1
Requirement	User will be able to logout
Note	- Logout will terminate the session - Automatic time-out will occur after pre-define number of hours - By default, number of hours for time-out is 12 hours

8.5 Password management(PAS)

8.5.1 Changing password

Number	8.5.1.1
Requirement	User will be able to set new password
Note	- From member profile setting page - By entering previous password and new password twice

Number	8.5.1.2
Requirement	User will be required set new password every 12 months
Note	- User will be able to 'Log Out To Reset Password' under the 'Current password' field - By entering previous password and new password twice

Number	8.5.1.3
Requirement	User will be able to see the new password entered
Note	By clicking on 'Show password'

Number	8.5.1.4
Requirement	User will be required to choose a new password that is not identical to previously selected passwords
Note	Limited for the last 6 passwords used

Number	8.5.1.5
Requirement	User will be required to set a new password with ‘Strong’ password strength
Note	Password must contain the following: • At least 11 characters long • Mix of lowercase, uppercase, numbers and special characters

Number	8.5.1.6
Requirement	User will receive email notification when password is updated
Note	When admin or user itself, from the member show page

Number	8.5.1.7
Requirement	Requirement removed
Note	

8.5.2 Section removed

8.6 Labguru members module^(MEM)

8.6.1 My profile

Number	8.6.1.1
Requirement	User will be able to see current logged member
Note	By clicking on profile image

Number	8.6.1.2
Requirement	User will be able view the following personal information: • Title and full name • Role • Mobile phone • Email • Privileges • Two-factor authentication
Note	Note removed

Number	8.6.1.3
Requirement	User will be able to edit the following personal information: • First name • Last name • Title • Mobile Phone • Shopping list notifications
Note	By clicking on 'Edit' in member show page

Number	8.6.1.4
Requirement	Requirement removed
Note	

Number	8.6.1.5
Requirement	User will be able set new password
Note	• By clicking on 'Set password' • User required to enter current password • User required to enter new password twice • Not applicable when Single Sign On is enabled for organization or environment

Number	8.6.1.6
Requirement	User will be able to send an email to other members of the lab
Note	By clicking on 'Messages'

Number	8.6.1.7
Requirement	Requirement removed
Note	

Number	8.6.1.8
Requirement	User will be able to get an API token
Note	By clicking on 'Labguru's 'Get API token'

Number	8.6.1.9
Requirement	User will be able to turn on/off the side-navigation menu auto-pop up when hovering
Note	

Number	8.6.1.10
Requirement	User will be able to set the delay time when auto-pop up of side-navigation it turned on
Note	From 0ms to 1000ms

Number	8.6.1.11
Requirement	User will be able to login to Labhandy with QR code
Note	By clicking ‘Sign in to Labhandy’ and scanning the QR code

Number	8.6.1.12
Requirement	User will be able to set default account
Note	When no default account is set, user will be directed to the last used account

8.6.2 Members list

Number	8.6.2.1
Requirement	User will see list of active members in account
Note	

Number	8.6.2.2
Requirement	User will be able to see list of former members
Note	

Number	8.6.2.3
Requirement	User will be able to add new members to account
Note	• Button will be available only if account did not reach seat cap • Fields First name/Last name/Email/Role are mandatory • Email will be sent to invitee • Invite email will remain active for 3 days

Number	8.6.2.4
Requirement	User will be able to go to member’s page
Note	By clicking on the member name or image

Number	8.6.2.5
Requirement	Requirement removed
Note	

Number	8.6.2.6
Requirement	User will be able to sort member according in ascending/descending order by any column added
Note	• By clicking on the column's title • By clicking on the kendo menu and select one of the sorting options • Sort preferences are saved

Number	8.6.2.7
Requirement	User will be able to filter member by any combination of columns
Note	Filter preferences are saved

Number	8.6.2.8
Requirement	User will be clear all filters/sorts
Note	By clicking on the cog menu and 'Clear all filters/sorts'

Number	8.6.2.9
Requirement	User will be able to view list of former lab members tab
Note	By clicking on 'former lab members'

Number	8.6.2.10
Requirement	User will be able to set the number of members displayed per page
Note	• Default number is 20 members per page. Can be changed to 50,100, or 200 • User preferences will be saved

Number	8.6.2.11
Requirement	User will be able to add and remove columns from the display
Note	• By clicking on the column menu -> 'Columns' and checking the wanted columns • Unchecking a column will remove it from view • User preferences will be saved

Number	8.6.2.12
Requirement	User will be able to adjust the width of the columns
Note	• By dragging and dropping the columns' borders • Auto-adjust by double clicking the border

Number	8.6.2.13
Requirement	User will be able to re-order columns
Note	• By dragging and dropping the columns • User preferences will be saved

Number	8.6.2.14
Requirement	User will be able to export members list

Note	
------	--

8.6.3 Member show page

Number	8.6.3.1
Requirement	User will be able to see the following information on member: • Name • Title • Role and privileges • Email • 2FA active/inactive • Teams • Profile image
Note	

Number	8.6.3.2
Requirement	User will be able to edit member's following information: • First name • Last name • Title • Mobile Phone • Role • Shopping list notifications
Note	User will be required to enter password to change changes

Number	8.6.3.3
Requirement	User will be able set new password
Note	• By clicking on 'Set password' • User required to enter current password • User required to enter new password twice

Number	8.6.3.4
Requirement	User will be able to lock member
Note	• Locked member will be removed from active members list • Locked member will not be able to login

Number	8.6.3.5
Requirement	User will be able to unlock former member
Note	Option to unlock former member is possible only if seat cap was not reached

Number	8.6.3.6
Requirement	User will be able to resend member invitation message
Note	

Number	8.6.3.7
Requirement	User will be able to send member API token by email
Note	

Number	8.6.3.8
Requirement	User will be able edit member's profile image
Note	

Number	8.6.3.9
Requirement	Admin user will be able to reset a member's 2FA
Note	By clicking on 'reset 2FA' on member show page

8.7 Labguru Weekly Digest Email module (WKG)

Number	8.7.1
Requirement	The weekly digest will contain the dates of the week
Note	

Number	8.7.2
Requirement	The weekly digest will contain the name of the lab
Note	

Number	8.7.3
Requirement	The weekly digest will not show experiments to member without permissions
Note	

Number	8.7.4
Requirement	The weekly digest will show experiments for witness
Note	

Number	8.7.5
Requirement	The weekly digest will show low stocks alerts
Note	

Number	8.7.6
Requirement	The weekly digest will report if there are no active projects
Note	

Number	8.7.7
Requirement	The weekly digest will report if no new key items were added in the account
Note	

Number	8.7.8
Requirement	The weekly digest will report a graph of new items added to the account
Note	

Number	8.7.9
Requirement	The weekly digest will report if account has expired stocks
Note	

Number	8.7.10
Requirement	Locked user (former) will not receive weekly digest email
Note	

Number	8.7.11
Requirement	User will be able to unsubscribe from Weekly digest
Note	

8.8 Labguru Teams module^(TMS)

8.8.1 General

Number	8.8.1.1
Requirement	User will be able to view list of all teams
Note	For each team, the following information will be available: • Team name • Team members (name, role, date added and privileges) • Team leader • Team description

Number	8.8.1.2
Requirement	Admin user will be able to create a team
Note	Team name and team members are mandatory. Team leader and description are optional

Confidential

Number	8.8.1.3
Requirement	Admin user will be able to edit all teams
Note	This includes: • Remove/change team leader • Add/remove members • Edit team's name • Edit team's description

Number	8.8.1.4
Requirement	Team leader will be able to edit team
Note	This includes: • Add/remove members • Edit team's name • Edit team's description

Number	8.8.1.5
Requirement	Admin user will be able to delete a team
Note	

8.8.2 Using teams

Number	8.8.2.1
Requirement	User will be able to filter projects by teams
Note	

Number	8.8.2.2
Requirement	User will be able to send email notification to teams when posting a comment
Note	

Number	8.8.2.3
Requirement	User will be able to give viewing permissions to a team
Note	Apply to projects, protocols, reports and documents

Number	8.8.2.4
Requirement	User will be able to remove viewing permissions to a team
Note	Apply to projects, protocols, reports and documents

Number	8.8.2.5
Requirement	User will be able to give editing permissions to a team
Note	Apply to experiments and protocols

Number	8.8.2.6
Requirement	User will be able to remove editing permissions to a team
Note	Apply to projects, experiments and protocols

Number	8.8.2.7
Requirement	Requirement removed
Note	

Number	8.8.2.8
Requirement	User that was removed from a team will lose all permissions that were given to him by being part of the team
Note	

Number	8.8.2.9
Requirement	User that was added to a team will gain all permissions that were given to his team
Note	

8.9 Section removed (SRV)

8.10 Section removed (LGT)

8.11 Section removed

8.12 Support menu

Number	8.12.1
Requirement	User will be able to download EasyEdit2 app for mac or for windows
Note	By clicking 'Download EasyEdit2'

Number	8.12.2
Requirement	User will be able to contact support team
Note	By clicking 'Contact Us' use will be able to send a message/ book a meeting/ search for help in articles

Number	8.12.3
Requirement	User will be able to search for articles by subject
Note	By clicking 'Help Desk'

Number	8.12.4
Requirement	User will be able see Biodata privacy policy
Note	By clicking 'Privacy'

Number	8.12.5
Requirement	User will be able see change log updates of Labguru versions
Note	

9. Department account^(DEP)

9.1 General

Number	9.1.1
Requirement	User in an account connected to a department account will be able to access the department account
Note	By clicking on 'Switch to...' in account dropdown menu

Number	9.1.2
Requirement	User will see in department account only the Equipment module
Note	Other modules will not be visible in department account

Number	9.1.3
Requirement	User will be able to add new equipment in department account
Note	

Number	9.1.4
Requirement	User will be able to view list of archived equipment
Note	By clicking on 'View archived equipment'

Number	9.1.5
Requirement	User will be able to set the number of equipment displayed per page
Note	Default number is 20 equipment per page. Can be changed to 50,100 or 200. User preferences will be saved

Number	9.1.6
Requirement	User will be able to sort equipment in ascending or descending order by any column header
Note	- By clicking on the column's title - By clicking on the kendo menu and select one of the sorting options - Sort preferences are saved

Number	9.1.7
Requirement	User will be able to search for equipment by name
Note	

Number	9.1.8
Requirement	User will be able to filter equipment by any combination of columns
Note	Filter preferences are saved

Number	9.1.9
Requirement	User will be able to filter equipment by tags
Note	• By clicking on the tag's name • User can cancel selection by clicking on same tag again • Clicking on multiple tags will filter for stocks that are tags with all selected tags

Number	9.1.10
Requirement	User will be able to clear all filters/sorts
Note	By clicking on the cog menu and 'Clear all filters/sorts'

Number	9.1.11
Requirement	User will be able to add and remove columns from the display
Note	• By clicking on the column menu -> 'Columns' and checking the wanted columns • Unchecking a column will remove it from view • User preferences will be saved

Number	9.1.12
Requirement	User will be able to adjust the width of the columns
Note	• By dragging and dropping the columns' borders • Auto-adjust by double clicking the border

Number	9.1.13
Requirement	User will be able to reorder columns
Note	• By dragging and dropping the columns • User preferences will be saved

9.2 Sharing equipment

9.2.1 On Sharing Lab account

Number	9.2.1.1
Requirement	Admin user will be able to share equipment to department account
Note	By clicking on 'Share' in equipment page

Number	9.2.1.2
Requirement	Admin user will be able to unshare equipment from department account

Confidential

Note	• When clicking on 'Share', the 'Share' button is being replaced with an 'Unshare' button and vice versa • Equipment will be removed from department and subscribed accounts • Events belongs to equipment will be deleted from department and subscribed accounts • Double confirmation is required
------	---

Number	9.2.1.3
Requirement	User will see which equipment was shared to the department account
Note	By viewing the 'Shared' column on the equipment index table, on lab's workspace

Number	9.2.1.4
Requirement	User will be able to edit equipment that was shared to department
Note	

Number	9.2.1.5
Requirement	User will be able to upload attachment to equipment shared to department
Note	Attachment will not appear in department nor in subscribed account

Number	9.2.1.6
Requirement	User will be able to tag equipment shared to department
Note	Tags will not appear in department nor in subscribed account

Number	9.2.1.7
Requirement	User will be able to delete shared equipment
Note	

Number	9.2.1.8
Requirement	User will be able to archive shared equipment
Note	Archived equipment will remain available in department and subscribed accounts

Number	9.2.1.9
Requirement	User will be able to schedule event on the calendar of equipment shared to department
Note	Event will be visible in department and in subscribed account

Confidential

Number	9.2.1.10
Requirement	User will be able to add new maintenance to shared equipment
Note	New maintenance will be visible in department and subscribed account

Number	9.2.1.11
Requirement	User will be able to edit maintenance type to shared equipment
Note	Changes in maintenance type will be visible in department and subscribed account

Number	9.2.1.12
Requirement	User will be able to remove maintenance from shared equipment
Note	Removed maintenance will be removed from department and subscribed account

Number	9.2.1.13
Requirement	User will be able to add maintenance event to equipment shared to department
Note	Maintenance event will be visible in department and in subscribed account

Number	9.2.1.14
Requirement	User will be able to edit maintenance event to equipment shared to department
Note	Changes in maintenance event will be visible in department and subscribed account

Number	9.2.1.15
Requirement	User will be able to add shared equipment to schedule
Note	The equipment event created from the schedule will be visible in department and subscribed accounts

9.2.2 On Department account

Number	9.2.2.1
Requirement	User will see equipment that was created in department account or shared to department account
Note	Shared equipment will be visible in equipment index with indication where it is shared from, under 'Shared' column

Confidential

Number	9.2.2.2
Requirement	Admin user will be able to subscribe/unsubscribe to equipment that was shared from another account or created in department account
Note	• By clicking on the 'Subscribe' button on the equipment show page and selecting account to subscribe the equipment to • Subscribed equipment will be copied to user's account(s) • Subscribed equipment will be visible in equipment index with indication it was "shared with" under 'Sharing' column

Number	9.2.2.3
Requirement	User will be able to edit only equipment created in department account
Note	Changes in equipment information will be visible in subscribed account

Number	9.2.2.4
Requirement	User will be able to delete only equipment created in department account
Note	• Equipment will not be deleted from subscribed account • Admin and equipment owner only

Number	9.2.2.5
Requirement	User will be able to add attachment to equipment page
Note	Only attachments created before equipment was shared to other account will be visible in subscribed account

Number	9.2.2.6
Requirement	User will be able to import/update equipment in department account
Note	• By clicking on cog menu -> 'Import/Update' • Changes made to equipment via update will be possible only for equipment created in department account

Number	9.2.2.7
Requirement	User will be able to see previous imports/updates made in department account
Note	

Number	9.2.2.8
Requirement	User will be able to export equipment

Confidential

Note	• By clicking on cog menu -> 'Export' • Export will include both equipment created in and shared to department account
------	--

Number	9.2.2.9
Requirement	User will be able to link equipment in department account
Note	Link will not be visible in sharing or subscribed accounts

Number	9.2.2.10
Requirement	User will be able to tag equipment shared to department
Note	Tag will not be visible in sharing or in subscribed accounts

Number	9.2.2.11
Requirement	User will be able to archive equipment created in department account
Note	Archived equipment will remain active in subscribed account

Number	9.2.2.12
Requirement	User will be able to create an event in equipment calendar in department account
Note	Event will be visible in sharing and subscribed accounts

Number	9.2.2.13
Requirement	User will be able to add/edit maintenance type only for equipment created in department account
Note	• Changes in equipment maintenance will appear in subscribed account • Admin and equipment owner only

Number	9.2.2.14
Requirement	User will be able to edit maintenance template only for equipment created in department account
Note	• Changes in maintenance template will appear in subscribed account • Admin and equipment owner only

Number	9.2.2.15
Requirement	User will be able to perform maintenance event only for equipment created in department account
Note	Maintenance event will appear in subscribed account

Number	9.2.2.16
Requirement	User will be able to edit a maintenance event only for equipment that was created in department account

Note	Changes in maintenance event will be appear in subscribed account
------	---

Number	9.2.2.17
Requirement	User will be able to see maintenance log
Note	

9.2.3 On Subscribing Lab account

Number	9.2.3.1
Requirement	User will not be able to edit equipment subscribed from department
Note	

Number	9.2.3.2
Requirement	User will be able to view the equipment location's name of equipment subscribed from department
Note	Equipment location's name is not linked to the equipment location

Number	9.2.3.3
Requirement	User can upload an attachment to equipment subscribed from department
Note	Attachment will not appear in sharing, department or another subscribed account

Number	9.2.3.4
Requirement	User can tag equipment subscribed from department
Note	Tag will not appear in sharing, department or another subscribed account

Number	9.2.3.5
Requirement	User will be able to delete equipment subscribed from department
Note	Deleted equipment will remain available in sharing, department or another subscribed account

Number	9.2.3.6
Requirement	User will be able to archive equipment subscribed from department
Note	Archived equipment will remain active in sharing, department or another subscribed account

Number	9.2.3.7
Requirement	User will be able to add event on the calendar of equipment subscribed from department
Note	Event will not appear in sharing, department or another subscribed account

Number	9.2.3.8
Requirement	User will not be able to add/edit maintenance type for equipment subscribed from department account
Note	

Number	9.2.3.9
Requirement	User will not be able to edit maintenance template for equipment subscribed from department account
Note	

Number	9.2.3.10
Requirement	User will not be able to perform a maintenance event for equipment subscribed from department account
Note	

Number	9.2.3.11
Requirement	User will not be able to edit a maintenance event for equipment subscribed from department account
Note	

Number	9.2.3.12
Requirement	User will be to add subscribed equipment to schedule
Note	The equipment event created from the schedule will be visible in department and subscribed accounts

9.3 Shared protocols

Number	9.3.1
Requirement	User will be able to share a protocol to the department protocols directory
Note	- By clicking on the 'Share to repository' button and then choosing the department directory from the 'Share protocol' popup - The protocol will be added to the department protocols directory - Available only for administrators - If the 'Repositories sharing' option is marked on the account settings

Number	9.3.2
Requirement	User will be able to access the department protocols directory
Note	By clicking on the 'Add from protocols directory' button and then choosing the department repository

Number	9.3.3
Requirement	User will be able to unshare a protocol from the department directory
Note	- By clicking on the 'Share/remove from the repository' button and then un-marking the department directory from the 'Share protocol' popup - The protocol will be removed from the department protocols directory

Number	9.3.4
Requirement	User will be able to view shared protocols
Note	By clicking on them from the department protocols directory

Number	9.3.5
Requirement	User will be able to print a shared protocol
Note	

Number	9.3.6
Requirement	User will be able to copy a shared protocol to his account
Note	- By clicking on the 'Copy protocol' button - The protocol will be added to the lab account - The 'Copy protocol' button will be replaced with a 'Copied' button

10. Labguru import/export (IMP)

10.1 Section Removed

10.2 Import/Update

10.2.1 General

10.2.1.1 Import page

Number	10.2.1.1.1
Requirement	User will be able to import items using a template
Note	By downloading import template file

Number	10.2.1.1.2
Requirement	User will be able to upload template file to Labguru
Note	

Number	10.2.1.1.3
Requirement	User will be notified when imported used invalid template
Note	Invalid templates can be templates of other collections or modules and templates that were modified by user after downloading

10.2.1.2 Excel Template

Number	10.2.1.2.1
Requirement	User will be required to enter Name field
Note	Item's name must be unique when unique name field is enabled

Number	10.2.1.2.2
Requirement	User will not be able to update 'Created at' field
Note	When filling the field, the system will ignore it

Number	10.2.1.2.3
Requirement	User will be able to select the owner from pre-defined list
Note	If left empty, the owner of the item will be the user who performed the import

Number	10.2.1.2.4
Requirement	User will be required to enter custom fields defined as mandatory
Note	

Number	10.2.1.2.5
Requirement	User will be able to select values from list of fields for custom fields from Pre-defined list or LOV type
Note	

Number	10.2.1.2.6
Requirement	User will be able to see the excel template columns separated by different colors for: • Pre-defined list • Custom fields • Parent collections • Sequence info • Stock info
Note	

Number	10.2.1.2.7
Requirement	User will be able to update an existing item/stock
Note	• By entering an existing ID in 'System ID' field

- | | |
|--|---|
| | <ul style="list-style-type: none"> For a new item, 'System ID' field should be empty |
|--|---|

10.2.1.3 Import Summary page

Number	10.2.1.3.1
Requirement	User will be able to view Import info summary: - Import status (in progress or completed) - Import id - File name - Started at - Import at - Items not imported (if file contain error) - Items imported successfully - Import errors (in a spreadsheet element)
Note	

Number	10.2.1.3.2
Requirement	User will be able to edit the import for failed items
Note	When import has errors, a spreadsheet will appear in import summary page showing only failed items

Number	10.2.1.3.3
Requirement	User will be able to see for each failed item the reason it failed
Note	Reason for failure will be added into column A in spreadsheet

Number	10.2.1.3.4
Requirement	User will be able to re-import failed items from spreadsheet
Note	By clicking the 'I'm done, Re-import now' button

Number	10.2.1.3.5
Requirement	User will be able to comment on the import summary page
Note	

Number	10.2.1.3.6
Requirement	User will be able to access the items that were imported successfully
Note	

10.2.2 Import in collections

10.2.2.1 Excel Template

Number	10.2.2.1.1
Requirement	Requirement removed

Note	
------	--

Number	10.2.2.1.2
Requirement	User will be required to enter 'Manufacturer' and 'Catalog no.' or 'Produced by' fields
Note	- Only when importing consumables or items from collection with purchasable attributes - Catalog number must be unique

Number	10.2.2.1.3
Requirement	User will be able to create tags for new items
Note	

Number	10.2.2.1.4
Requirement	User will be able to import compounds with an existing structure
Note	- By adding a value in 'Structure (SMILES)' or 'Structure (InChI)' fields in an import template - Duplicates will be linked to their original compound

Number	10.2.2.1.5
Requirement	User will not be able to update auto calculated fields in collections
Note	

10.2.2.2 Import Collection Summary page

Number	10.2.2.2.1
Requirement	User will be able to edit imported item
Note	

Number	10.2.2.2.2
Requirement	User will be able to search imported items by name and sysID
Note	

Number	10.2.2.2.3
Requirement	User will be able to filter imported items by any combination of columns in display
Note	

Number	10.2.2.2.4
Requirement	User will be able to tag imported items

Note	By selecting items and clicking on ‘Tag’
------	--

Number	10.2.2.2.5
Requirement	User will be able to filter imported items by tags
Note	

Number	10.2.2.2.6
Requirement	User will be able to select which columns to view on the index page
Note	

Number	10.2.2.2.7
Requirement	User will be able to sort items by any column in display in ascending/descending order
Note	

Number	10.2.2.2.8
Requirement	User will be able to delete selected imported items
Note	

Number	10.2.2.2.9
Requirement	User will be able to archive imported item
Note	

10.2.3 Import Stocks

Number	10.2.3.1
Requirement	Requirement removed
Note	

Number	10.2.3.2
Requirement	User will be able to select Privacy (public or private)
Note	If no option selected, stock privacy will be public

Number	10.2.3.3
Requirement	User will be able to select Stock type from drop-down list
Note	If no option selected, stock type will be tube

Number	10.2.3.4
Requirement	User will be able to select Stock color from drop-down list
Note	

Confidential

Number	10.2.3.5
Requirement	User will be able to select Concentration/Volume/Weight units from drop-down list
Note	

Number	10.2.3.6
Requirement	User will be able to create multiple identical stocks per item
Note	By using the Stock Count field

Number	10.2.3.7
Requirement	User will be able to select Stock owner
Note	If left empty, the owner will be the user who performed the import

Number	10.2.3.8
Requirement	User will be able to select Stored/frozen by
Note	If left empty, the Stored/Frozen by will be the user who performed the import

Number	10.2.3.9
Requirement	User will be required to enter either name of box in 'Box name' field (new/existing) or select storage location from drop-down list to place the stock
Note	- If box name does not exist in account, a new box will be created with that name - If creating a new box, the value under 'Storage location' will be used as the location of the new box - If creating a new box, and leaving Storage location empty, the new box will be created as an unstored box - If entering both name of existing box and storage location, stocks will be created in the box

Number	10.2.3.10
Requirement	User will be able to enter Box dimensions (#rows or #columns)
Note	- Relevant only if box is created via import - If left empty, box dimensions will be 10 ×10

Number	10.2.3.11
Requirement	User will be able to select Stock position in a box
Note	- Enter numeric or alphanumeric position (for example 1 or A1) - Position must be valid and available in the selected box - If left empty, stock will be placed in the first available slot in the box

Number	10.2.3.12
Requirement	User will be required to enter either Inventory item name or SysID
Note	• If entering both, item name will be used • If there are multiple items with the entered name, stock will be associated to the item that was created earlier • If entering item name that does not exist in collection, a new item will be created

Number	10.2.3.13
Requirement	User will be able to add 'Box location in Rack-cells'
Note	• By entering cell number • If left empty but selecting rack type as storage location, stock will be placed in the first available cell in the rack • If entering both name of existing box and cell number in 'Box location in Rack-cells', stocks will be created in the box

10.2.4 Import Equipment

Number	10.2.4.1
Requirement	Requirement removed
Note	

Number	10.2.4.2
Requirement	User will be able to add Manufacturer
Note	• Creating a new manufacture will add it in the companies page • If manufacturer already exist it will add it to the company

Number	10.2.4.3
Requirement	User will be able to create tags for new equipment
Note	

Number	10.2.4.4
Requirement	User will be able to select to locate the equipment into a storage location from drop-down list
Note	

10.3 Export

Number	10.3.1
Requirement	User will be able to export an excel file
Note	Export file will be received by email

Number	10.3.2
Requirement	Requirement removed
Note	

Number	10.3.3
Requirement	Requirement removed
Note	

10.4 Bulk-Update

Number	10.4.1
Requirement	User will be able to bulk-update collection/stock/equipment using export-import
Note	• By clicking Export and edit the excel or Edit the data • Export file will be received by email • Do not delete existing data from file • Large exports are created as CSV files and should be copied to an empty import template

Number	10.4.2
Requirement	User will be able to bulk-update collection/stock/equipment using Edit
Note	By clicking Edit on the toolbar

Number	10.4.3
Requirement	User will be able to bulk update collection using edit (full collection/selected items)
Note	By clicking on 'Edit' on the toolbar

11. Section removed

12. Labguru Workflow Editor module^(WFE)

12.1 General

Number	12.1.1
Requirement	User with permission to access the Workflow Editor and Dashboards will be able to do it from the account drop down menu
Note	Admins and members with a role that has "can access Workflow Editor" and Dashboards checked

Number	12.1.2
Requirement	User will be able to go back to flows index page from any page
Note	By clicking on the Labguru logo

Number	12.1.3
Requirement	User will be able to move to any existing flow page in the account from any flow page
Note	

Number	12.1.4
Requirement	User will be able to logout from the Editor and be directed back to the Labguru account
Note	

Number	12.1.5
Requirement	User will be able to view JSON samples of different entities in the system
Note	

Number	12.1.6
Requirement	User will be able to see the number of exceptions occurred in account since last time visiting the 'Expectations' page
Note	

12.2 Flows index page

Number	12.2.1
Requirement	Upon entering the Workflow Editor, user will see a list of all created workflows
Note	

Number	12.2.2
Requirement	User will see for each flow the following columns: • Id

	• Title • Status • Run counter • Last run • Description • Trigger • Created at
Note	status can be either "Active" or "Inactive"

Number	12.2.3
Requirement	User will be able to sort the list by any of the columns
Note	By default, all columns are visible

Number	12.2.4
Requirement	User will be able to choose which columns can be viewed in the list
Note	By default, all columns are visible

Number	12.2.5
Requirement	User will be able to enter each of the flow
Note	User will be directed to the 'Editor' tab in the flow page

Number	12.2.6
Requirement	User will be able to add a new workflow
Note	

Number	12.2.7
Requirement	User will be able re-register all webhooks
Note	

Number	12.2.8
Requirement	User will be able to de-activate all flows
Note	

12.3 Flow page - Editor tab

Number	12.3.1
Requirement	User will be able to add/change trigger
Note	

Confidential

Number	12.3.2
Requirement	User will be able to add a new following step after every step in the flow
Note	

Number	12.3.3
Requirement	User will be able to add a new branch to steps that already have a following step
Note	

Number	12.3.4
Requirement	User will be able to save every step
Note	

Number	12.3.5
Requirement	User will be able to see steps API docs
Note	

Number	12.3.6
Requirement	User will be able to delete each step
Note	User will be asked to confirm once after trying to delete

Number	12.3.7
Requirement	User will see logs of last 3 un-archived runs
Note	Will appear in a panel to the right of the steps display

Number	12.3.8
Requirement	User will be able to re-run one of the last 3 unarchived runs
Note	

Number	12.3.9
Requirement	User will be able to archive one of the last 3 unarchived runs
Note	

Number	12.3.10
Requirement	User will be able to view all variables created as part of the flow
Note	Divided in to variables created from pre-built steps (Variables) and from scripter steps (Generic steps)

Number	12.3.11
Requirement	User will be able to set the flow to be active or not active

Note	Inactivated flows will not start when trigger is happening
------	--

Number	12.3.12
Requirement	User will be able to set the flow to be synced or async
Note	Async flows will run immediately without queuing order

Number	12.3.13
Requirement	User will be able to re-register the webhook
Note	

Number	12.3.14
Requirement	User will be able to export the flow as a yaml file
Note	

Number	12.3.15
Requirement	User will be able to import a flow from a yaml file
Note	User will need to select trigger after uploading the file to complete the import

Number	12.3.16
Requirement	User will be able to open the WF API documentation page
Note	

12.4 Flow page - Runs tab

Number	12.4.1
Requirement	User will be able to switch the list view between only active/ only archived/ both
Note	

Number	12.4.2
Requirement	User will be able to re-run every run
Note	

Number	12.4.3
Requirement	User will be able to archive/un-archive every run
Note	

Number	12.4.4
Requirement	User will be able to refresh the page
Note	

Number	12.4.5
Requirement	User will be able to select how many runs to present on page
Note	Selection between 7,30 or 300 runs by moving between pages

Number	12.4.6
Requirement	User will be able to filter list of run to show only successful runs
Note	

Number	12.4.7
Requirement	User will be able to filter list of runs by run date
Note	

12.5 Triggers

Number	12.5.1
Requirement	User will be able to set webhook based trigger
Note	

Number	12.5.2
Requirement	User will be able to set periodic trigger
Note	By using cron schedule expressions

Number	12.5.3
Requirement	User will be able to set external trigger
Note	When setting external trigger, user will get url to be used when calling for the activation of the flow

Number	12.5.4
Requirement	User will be able to set manual trigger
Note	

12.6 Steps

12.6.1 Conditions

Number	12.6.1.1
Requirement	User will be able to add condition steps, selecting between numerical, text, exists? and between operator types
Note	

Number	12.6.1.2
Requirement	User will be able to set the condition on ‘Stop’ or ‘Continue’ if condition is met
Note	

Number	12.6.1.3
Requirement	User will be able to add parameters and operator for condition
Note	Based on type of condition selected

Number	12.6.1.4
Requirement	User will be able to set flow to archive when condition is not met
Note	

12.6.2 Scripter

Number	12.6.2.1
Requirement	User will be able to add scripter step to flow
Note	

Number	12.6.2.2
Requirement	User will be able to name scripter step
Note	

Number	12.6.2.3
Requirement	User will be able to select scripter language
Note	Selection between: Python, Ruby

Number	12.6.2.4
Requirement	User will be able to connect scripter to ‘Production’ or ‘Tester’ code templates
Note	

Number	12.6.2.5
Requirement	User will be able to type code
Note	

Number	12.6.2.6
Requirement	Requirement removed
Note	

Number	12.6.2.7
Requirement	User will be able to name generic variables

Note	
------	--

Number	12.6.2.8
Requirement	User will be able to select code language
Note	Supported languages: Ruby, Python

Number	12.6.2.9
Requirement	User will be able to make API commands to Labguru
Note	- By using get_api, post_api, put_api commands - Commands should work with all supported languages

Number	12.6.2.10
Requirement	User will be able to add to the run's log
Note	By using the log() command

12.7 Code templates

Number	12.7.1
Requirement	User will be able to write code template in Ruby or Python
Note	

Number	12.7.2
Requirement	User will be able to write templates as 'Production' or 'Tester'
Note	Each option is available in its own window, under 'Code templates' menu

Number	12.7.3
Requirement	User will be able to save changes made in code templates
Note	

12.8 Exceptions

Number	12.8.1
Requirement	User will see 200 recently occurred exceptions in a list
Note	For each exception: Flow, Error and date occurred

Number	12.8.2
Requirement	User will be able to go directly to flow page by clicking on flow's name

Note	
------	--

Number	12.8.3
Requirement	User will be able to mark exception as read
Note	

Number	12.8.4
Requirement	User will be able to mark all exceptions on page as read
Note	

Number	12.8.5
Requirement	User will be able to adjust email notification setting on exceptions
Note	Selection between immediate notifications, daily report or weekly report

12.9 Flow page - Variables tab

Number	12.9.1
Requirement	User will be able to add variable
Note	- By adding key and value - Variables can be used in the flow, in scripter steps using <i>variable(key)</i> command

Number	12.9.2
Requirement	User will be able to edit variable
Note	By changing key, value or both

Number	12.9.3
Requirement	User will be able to remove variable
Note	

13. Labguru Audit trail module^(ADT)

13.1 General Audit Trail page

Number	13.1.1
Requirement	User will be able to view every action in Labguru as: - Creation of a new entity - Updating an entity - Deletion of an entity
Note	

Confidential

Number	13.1.2
Requirement	User will be able to filter actions by: • Date • Active member • Teams • Former member • Action
Note	

Number	13.1.3
Requirement	User will be able to view more audit entries
Note	• By clicking on the double down arrow button • By default, the list will show 25 entries • Every click will add 25 more entries to list

Number	13.1.4
Requirement	User will be able to show/hide system calls
Note	By clicking on ‘Include system calls’

Number	13.1.5
Requirement	User will be able to refresh page to show new actions
Note	By clicking on ‘Refresh

Number	13.1.6
Requirement	User will be able to print current audit trail view into pdf
Note	The pdf will include only audit trail entries that are on screen

Number	13.1.7
Requirement	User will be able to view the required entries information: • Timestamp • User ID • Action type • Event description • Object • Context
Note	Entry made by auto-save action will be displayed in ‘Action’ column

Number	13.1.8
Requirement	User will not be able to edit an audit trail entry
Note	

Number	13.1.9
Requirement	User will be able to clear filters to reset the audit trail view
Note	

13.2 Diff column

Number	13.2.1
Requirement	When editing spreadsheet, user will see sheet, position, old value and new value
Note	

Number	13.2.2
Requirement	When editing input field in form, user will see field name, old value and new value
Note	

Number	13.2.3
Requirement	User will be able to see Reason for Change value when editing form element in experiment
Note	

Number	13.2.4
Requirement	When duplicating row in form element, user will see the new duplicated fields data
Note	Each input field in a new duplicate row will have a unique name (field name plus a serial identifier)

Number	13.2.5
Requirement	User will be able to see the changed values when editing sample element
Note	Name and ID of stock/item will be displayed User will be able to see if item/stock is added or removed

Number	13.2.6
Requirement	User will be able to see the changed values when editing equipment element
Note	Name and ID of stock/item will be displayed User will be able to see if equipment is added or removed

14. Components

14.1 DNA viewer

14.1.1 General

Number	14.1.1.1
Requirement	User will be able to select if to view the sequence as circular or linear
Note	Default is 'Circular'

Number	14.1.1.2
Requirement	User will be able to view sequences properties
Note	

Number	14.1.1.3
Requirement	User will be able to show/hide features
Note	By default, features will be shown

Number	14.1.1.4
Requirement	User will be able to show/hide primers
Note	By default, primers will be shown

Number	14.1.1.5
Requirement	User will be able to show/hide cut sites
Note	By default, one-cutters will be shown

Number	14.1.1.6
Requirement	User will be able to show/hide ORF translations
Note	By default, ORF translations will be hidden

Number	14.1.1.7
Requirement	User will be able to switch between 'Locked' and 'Edit' modes
Note	By default, sequence will be locked

Number	14.1.1.8
Requirement	User will be able to search for a sequence or an object
Note	

Number	14.1.1.9
Requirement	User will be able to undo/redo last action
Note	

Confidential

Number	14.1.1.10
Requirement	User will be able to view the nucleotide positions (Base numbers)
Note	

Number	14.1.1.11
Requirement	User will see features and primers names on map
Note	When hovering over the object

Number	14.1.1.12
Requirement	User will be able to select part of the sequence
Note	• By Clicking left click and dragging • Selection can be done from the map side or the sequence side • Selection will be highlighted

Number	14.1.1.13
Requirement	User will be able to select an object
Note	• By clicking on object in map side or sequence side • Selection will be highlighted

Number	14.1.1.14
Requirement	User will be able cancel selection
Note	By left-clicking anywhere in map

Number	14.1.1.15
Requirement	user will be able to view the Tm for a selected sequence
Note	

Number	14.1.1.16
Requirement	User will be able to view the length of the whole sequence
Note	

Number	14.1.1.17
Requirement	User will be able to visualize digestion with restriction enzymes
Note	

Number	14.1.1.18
Requirement	The user will be able to export the plasmid map and sequence
Note	Export options are as JSON, FASTA or gb files

14.1.2 Actions menu

Number	14.1.2.1
Requirement	User will be able to copy the sequence to clipboard using the following options: • Copy sequence • Copy reverse complement • Copy AA translation • Copy reverse complement AA translation • Copy complement
Note	

Number	14.1.2.2
Requirement	User will be able to add a new feature
Note	Only in ‘Edit’ mode

Number	14.1.2.3
Requirement	User will be able to add a new primer
Note	Only in ‘Edit’ mode

Number	14.1.2.4
Requirement	User will be able to add a new translation
Note	

Number	14.1.2.5
Requirement	User will be able to add a new reverse translation
Note	

Number	14.1.2.6
Requirement	User will be able to change sequence to lower/upper case
Note	

Number	14.1.2.7
Requirement	User will be able to edit feature/primer
Note	Only in ‘Edit’ mode

Number	14.1.2.8
Requirement	User will be able to delete feature/primer
Note	Only in ‘Edit’ mode

Number	14.1.2.9
Requirement	User will be able to inverse selection
Note	

14.1.3 Editing sequence

Number	14.1.3.1
Requirement	User will be able to cut part of sequence
Note	Only in 'Edit' mode

Number	14.1.3.2
Requirement	User will be able to add nucleotides to the sequence
Note	Only in 'Edit' mode

Number	14.1.3.3
Requirement	User will be able to replace nucleotides in the sequence
Note	Only in 'Edit' mode

Number	14.1.3.4
Requirement	User will be able to add a copied sequence to the sequence
Note	Only in 'Edit' mode

14.1.4 Quick Search in viewer

Number	14.1.4.1
Requirement	User will see number of hits/number of matching objects
Note	

Number	14.1.4.2
Requirement	User will be able to scroll between results
Note	Current selected results will be highlighted

Number	14.1.4.3
Requirement	User will be able to select DNA or AA search
Note	

Number	14.1.4.4
Requirement	User will be able to select literal or ambiguous search
Note	

Number	14.1.4.5
Requirement	User will be able to highlight all search results
Note	

Number	14.1.4.6
Requirement	User will be able to switch to expanded search menu mode
Note	

14.2 EasyEdit - Desktop application

14.2.1 EasyEdit setup and access

Number	14.2.1.1
Requirement	User will be able to install EasyEdit
Note	By downloading the 'Easy Edit' app for mac or for windows

Number	14.2.1.2
Requirement	User will be able to access EasyEdit menu
Note	By clicking on EasyEdit icon in system tray (Windows) and menu bar items (MacOS)

Number	14.2.1.3
Requirement	User will be able to quit EasyEdit
Note	

Number	14.2.1.4
Requirement	User will be able to auto-update EasyEdit to a new version
Note	By receiving a notification to download a new version of EasyEdit and accepting the update

14.2.2 EasyEdit connection

Number	14.2.2.1
Requirement	User will be required to authenticate EasyEdit with a Labguru account
Note	- Log in with Labguru credentials in Labguru web app. Account must match EasyEdit entered host. - Each first sign-in to an account creates a Labguru account folder with the following sub-folders: 'Backup', 'Uploader' and 'Local Edit'.

Number	14.2.2.2
Requirement	User will be able to specify the host of their Labguru account
Note	

Number	14.2.2.3
Requirement	User will be able to reset the host to default
Note	

Number	14.2.2.4
Requirement	User will be able to open Labguru web app
Note	By clicking 'Open Labguru'

Number	14.2.2.5
Requirement	User will be able to force connection between EasyEdit and Labguru account
Note	• By clicking 'Connect' • Available if a connection has not yet been established between EasyEdit and a signed-in Labguru account

Number	14.2.2.6
Requirement	User will be able to restart connection between EasyEdit to Labguru account
Note	• By clicking 'Reconnect' • Available if EasyEdit is already connected to a Labguru account

14.2.3 Editing and uploading attachments

Number	14.2.3.1
Requirement	User will be able to access EasyEdit folder
Note	• By clicking 'Open Folder' the current logged in Labguru account folder opens. • Pop-up message with instructions about the folders is presented

Number	14.2.3.2
Requirement	User will be able to add equipment
Note	• By clicking 'Add Instrument' and choosing an equipment from the associated account equipment list • A folder with the equipment name will be created in the EasyEdit folder

Number	14.2.3.3
Requirement	User will be able attach files to an equipment
Note	• By adding files into the equipment's folder • Attachments can be viewed in the instrument's output files and in the Dashboard 'unattached files' or 'unattached images' • Files uploaded to subfolders in the folder are not uploaded

Number	14.2.3.4
Requirement	User will be able to upload files to Labguru
Note	• By adding files into 'uploader' folder. • Newly uploaded files will be in the ELN under the dashboard 'unattached images' or 'unattached files'.

	Files uploaded to subfolders in the folder are not uploaded
--	---

Number	14.2.3.5
Requirement	User will be able to edit the attachment
Note	- By clicking on 'Easy Edit' in the attachment element file and locally editing the file. - EasyEdit opens attachments locally, creating a 'local-edit' copy and a session to track file changes - Download pop-up for EasyEdit appears when not installed or opened - Unable to edit signed witness experiment/protocol/document - File is protected for 60 minutes from being edited by another user - Saving the section disconnects the file from EasyEdit - When editing multiple files at the same time, changes will only be saved in one file

Number	14.2.3.6
Requirement	User will be able to save changes into Labguru attachments
Note	By saving the changes in the locally edited attachment. Upon saving, EasyEdit re-uploads the file into Labguru unless: - User disconnected from account or EasyEdit - Connection time limit expired (one hour) and another user connected to the file

Number	14.2.3.7
Requirement	User will be able to see changes he made in excel elements that were added to the page, in the file itself
Note	

14.2.4 Logs and Information

Number	14.2.4.1
Requirement	User will be able to access log files folder
Note	Log files are found in 'EasyEdit config' folder and can accessed by clicking on 'Logs folder'

Number	14.2.4.2
Requirement	User will be able to access connected file information
Note	By clicking 'Connected Files,' user can view: Current and past connections of files with EasyEdit

	Current file time session
--	---

Number	14.2.4.3
Requirement	User will be able to view the recent activity log
Note	- By clicking 'Logs' - Log contains details about system events and activities

Number	14.2.4.4
Requirement	User will be able to view EasyEdit product information
Note	

Number	14.2.4.5
Requirement	User will be able to view EasyEdit license information
Note	

Number	14.2.4.6
Requirement	User will be able to control notifications
Note	- By clicking 'Disable notifications' or 'Enable Notifications' - Windows / MacBook notifications indicate file connection status and upload success

14.2.5 EasyEdit Geneious Prime plug-in

Number	14.2.5.1
Requirement	User will be required to download EasyEdit plugin for Geneious Prime
Note	Download available from Geneious Prime software

Number	14.2.5.2
Requirement	User will be required to install EasyEdit plugin in Geneious Prime
Note	Labguru Sync Folder is created in Geneious Prime

Number	14.2.5.3
Requirement	User will be able to uninstall EasyEdit plugin in Geneious Prime
Note	

Number	14.2.5.4
Requirement	User will be able to add plasmid items to Labguru
Note	By adding Geneious Prime files to Labguru Sync folder

	• EasyEdit must be open and connected to a Labguru account • Only Geneious Prime file (.g) can be uploaded to Labguru • Geneious Prime files are added to 'Plasmids' collection
--	---

14.3 Small Molecule Drawer^(SMD)

14.3.1 General

Number	14.3.1.1
Requirement	User will be able to draw a compound in 2D
Note	

Number	14.3.1.2
Requirement	User will be able to clear the drawing area
Note	By clicking 'Clear canvas' in the top toolbar

Number	14.3.1.3
Requirement	User will be able to import compounds/reactions from file
Note	By uploading from a file, dragging, or entering manually input Supported file types: • DL Molfile (v2000 and v3000) • RXN file format • InChI String • ChemAxon Extended SMILES • ChemAxon Extended CML file formats

Number	14.3.1.4
Requirement	User will be able to save compounds/reactions to molecular file
Note	Supported file types: • DL Molfile (v2000 and v3000) • RXN file format • InChI String • ChemAxon Extended SMILES • ChemAxon Extended CML file formats • PNG • SVG

Number	14.3.1.5
Requirement	User will be able to copy, paste and cut structures from and on the canvas
Note	

Confidential

Number	14.3.1.6
Requirement	User will be able to undo/redo last action
Note	

Number	14.3.1.7
Requirement	User will be able to Aromatize/De-aromatize a structure
Note	

Number	14.3.1.8
Requirement	User will be able to optimize the layout of a structure on the canvas
Note	By clicking 'Layout' tool

Number	14.3.1.9
Requirement	User will be able to improve structure appearance
Note	Clicking the 'Clean Up' tool will adjust bond lengths and angles in the structure

Number	14.3.1.10
Requirement	User will be able to calculate R/S and E/Z configurations
Note	By clicking 'Calculate CIP'

Number	14.3.1.11
Requirement	User will be able to check a structure for errors
Note	By clicking 'Check structure'

Number	14.3.1.12
Requirement	User will be able to view the calculated values of compounds displayed on the canvas: • Chemical Formula • Molecular Weight • Exact mass • Elemental analysis
Note	Molecular weight and exact mass decimal places are editable

Number	14.3.1.13
Requirement	User will be able to add/remove explicit hydrogens
Note	

Number	14.3.1.14
Requirement	User will be able to view structure in three-dimensional
Note	

Confidential

Number	14.3.1.15
Requirement	User will be able to view structure in fullscreen mode
Note	

Number	14.3.1.16
Requirement	User will be able to adjust the zoom display
Note	

Number	14.3.1.17
Requirement	User will be able to navigate on the Canvas using a hand tool
Note	

Number	14.3.1.18
Requirement	User will be able to partially/fully delete any element of the structure
Note	• By clicking 'Erase' • Deleting one atom/bond by selecting 'Rectangle'/'Lasso' options + 'Erase' • Deleting full structure by selecting 'Fragment' or 'Lasso'/'Rectangle'+ ctrl click + 'Erase' • Deleting multiple atoms/bonds/structures or other elements 'Lasso'/'Rectangle'+ shift click+ 'Erase'

Number	14.3.1.19
Requirement	User will be able to partially/fully select elements to edit the structure
Note	By selecting one of the options: • Rectangle selection • Lasso selection • Fragment selection

Number	14.3.1.20
Requirement	User will be able to change structure display vertical/horizontal flip
Note	By clicking vertically or horizontally flip using the following options: • 'Fragment' option • 'Lasso'/'Rectangle'+ ctrl click

Number	14.3.1.21
Requirement	User will be able to rotate selected objects on the canvas
Note	By selecting the object and clicking on the rotating tool

14.3.2 Drawing Bonds

Number	14.3.2.1
Requirement	User will be able to add a bond type to a structure
Note	From the drop-down list/ by right-clicking: • Single • Double • Triple • Dative • Hydrogen • Single up • Single down • Single Up/Down • Double Cis/Trans

Number	14.3.2.2
Requirement	User will be able to edit an existing bond type
Note	• From drop down list/ by right-clicking • User will be able to change the bond type by clicking on it (single/ double/ triple)

Number	14.3.2.3
Requirement	User will be able to draw/add a carbon chain to the structure
Note	

14.3.3 Drawing Atoms

Number	14.3.3.1
Requirement	User will be able to add a specific atom to the structure
Note	From the atoms toolbar/ Periodic table or by double-click (or right-click) an atom to open its properties dialogue

Number	14.3.3.2
Requirement	User will be able to change an atom or its properties
Note	By left-clicking and edit

Number	14.3.3.3
Requirement	User will be able to increase/decrease the charge of an atom in a structure
Note	By selecting the 'Charge Plus' or 'Charge Minus' tool

14.3.4 Marking S-Groups

Number	14.3.4.1
Requirement	The user will be able to define a repeating unit in a structure as a group
Note	- By clicking 'S-Group' properties tools and selecting multiple group type - Repeat count range from 1 to 200

14.3.5 Template toolbar

Number	14.3.5.1
Requirement	The user will be able to add rings to the structure
Note	User will be able to select the connection as single bond, bond-to-bond or a fused structure

Number	14.3.5.2
Requirement	The user will be able to add a structure from the custom template library
Note	By clicking 'Structure library' and selecting structure from template library

Number	14.3.5.3
Requirement	User will be able to add functional groups to the structure
Note	- By clicking 'Structure library' and then 'Functional groups' tab - Selection from predefined functional groups

Number	14.3.5.4
Requirement	User will be able to expand functional groups
Note	By right-clicking on the functional group and expand abbreviation

Number	14.3.5.5
Requirement	User will be able to remove functional groups
Note	By right-clicking on the functional group and remove abbreviation

15. **Requests**

15.1 Labguru Request Templates module

15.1.1 Request Templates index page

Number	15.1.1.1
Requirement	User will be able to add a new request template
Note	

Confidential

Number	15.1.1.2
Requirement	User will be able to filter request templates by tags
Note	

Number	15.1.1.3
Requirement	User will be able to search for a request template by name
Note	

Number	15.1.1.4
Requirement	User will be able to change the number of request templates in a single page
Note	User will see 20 per page by default and can change to 5, 50,100, and 200 documents per page

Number	15.1.1.5
Requirement	User will be able to sort template requests by any column in display, by ascending/descending order
Note	

Number	15.1.1.6
Requirement	User will be able to filter by any combination of columns in display
Note	

Number	15.1.1.7
Requirement	User will be able to choose which columns to view in index page
Note	- By choosing headers from list of columns option - User's preference will be saved

Number	15.1.1.8
Requirement	User will be able to change the width of the column
Note	By dragging and dropping the column borders

Number	15.1.1.9
Requirement	User will be able to see list of active template requests that user has access to view
Note	

Number	15.1.1.10
Requirement	User will be able to view archived request templates

Note	
------	--

15.1.2 Request Templates show page

Number	15.1.2.1
Requirement	User will be able to change request template name
Note	

Number	15.1.2.2
Requirement	User will be able to tag a template request
Note	

Number	15.1.2.3
Requirement	User will be able to select request type
Note	

Number	15.1.2.4
Requirement	User will be able to open a new request from request template
Note	

Number	15.1.2.5
Requirement	User will be able to sign request template
Note	When another user is editing the request template page, the sign function will be disabled, and a message will appear on hover indicating the page is being edited

Number	15.1.2.6
Requirement	User will be able to duplicate template request
Note	

Number	15.1.2.7
Requirement	Request template owner will be able to set edit permission for other member/teams
Note	By default, only admins, owner will have edit permissions

Number	15.1.2.8
Requirement	Request template owner will be able to set view permission for other member/teams
Note	By default, request template is visible to all members of account

Confidential

Number	15.1.2.9
Requirement	User will be able to link request template to other pages
Note	• By clicking the linked resources button on the 'More' menu • By clicking the 'Add link' button on the 'Linked resources' section on the bottom of the page

Number	15.1.2.10
Requirement	User will be able to view version history of a request template
Note	

Number	15.1.2.11
Requirement	User will be able to add a new section to request template
Note	By clicking on the 'add section' link and choosing protocol to load from

Number	15.1.2.12
Requirement	User will be able to add a comment to request template
Note	

Number	15.1.2.13
Requirement	User will be able to download all attachments added to request template
Note	

Number	15.1.2.14
Requirement	User will be able to archive request template
Note	

Number	15.1.2.15
Requirement	User will be able to restore archived request template
Note	

Number	15.1.2.16
Requirement	User will be able to print request template to pdf file
Note	

15.1.3 Actions in Request templates section

Number	15.1.3.1
Requirement	User will be able to add to request template section the following elements:

	• Text • Spreadsheet • Steps • Samples • Plate • Equipment • Form • Sketch • Compound • Reaction • Code • Page break • Attachments
Note	

15.2 Labguru Request module

15.2.1 Request Index Page

Number	15.2.1.1
Requirement	User will be able to start a new request by selecting request template
Note	By clicking 'New request' and selecting a request template opening the drop-down menu

Number	15.2.1.2
Requirement	User will be able to create a new request template
Note	By clicking 'New request' and creating a new template request

Number	15.2.1.3
Requirement	User will be able to filter request by tags
Note	

Number	15.2.1.4
Requirement	User will be able to search for a request by name
Note	

Number	15.2.1.5
Requirement	User will be able to change the number of request in a single page
Note	User will see 20 per page by default and can change to 5, 50, 100, and 200 documents per page

Confidential

Number	15.2.1.6
Requirement	User will be able to sort requests by any column in display, by ascending/descending order
Note	

Number	15.2.1.7
Requirement	User will be able to filter by any combination of columns in display
Note	

Number	15.2.1.8
Requirement	User will be able to choose which columns to view in index page
Note	• By choosing headers from list of columns option • User's preference will be saved

Number	15.2.1.9
Requirement	User will be able to change the width of the column
Note	By dragging and dropping the column borders

Number	15.2.1.10
Requirement	User will be able to see list of active requests that user has access to view
Note	

15.2.2 Request Show Page

Number	15.2.2.1
Requirement	User will be able to change request name
Note	

Number	15.2.2.2
Requirement	User will be able to tag a request
Note	

Number	15.2.2.3
Requirement	Requirement removed
Note	

Number	15.2.2.4
Requirement	User will be able to change request priority
Note	• Urgent, high, medium or low

Confidential

	As default, the priority is low
--	---

Number	15.2.2.5
Requirement	User will be able to print request to pdf file
Note	

Number	15.2.2.6
Requirement	User will be able to archive request
Note	

Number	15.2.2.7
Requirement	User will be able to restore archived request
Note	

Number	15.2.2.8
Requirement	User will be able to link request to other pages
Note	

Number	15.2.2.9
Requirement	User will be able to add a new section to request
Note	By clicking on the ‘add section’ link and choosing protocol to load from

Number	15.2.2.10
Requirement	User will be able to add a comment to request
Note	

Number	15.2.2.11
Requirement	User will be able to download all attachments added to request
Note	

Number	15.2.2.12
Requirement	User will be able to delete request
Note	- Only requestor/system admin - Completed request cannot be deleted

Number	15.2.2.13
Requirement	User will be able to duplicate request
Note	

15.2.3 Navigating Request Flow

Number	15.2.3.1
Requirement	Requester will be able to submit draft request to Pending for review
Note	• New status will be Pending • Request will be uneditable to requester • Request will be editable to requester team leader/admin • Request will be assigned to team leader • Requester team leader will receive email notification • Requester will be able to select a team leader as approver • When the requestor is part of one team, the team leader will be assigned by default as approver • When the requester is part of more than one team, team leader can be selected from drop down list

Number	15.2.3.2
Requirement	Requester will be required to select executing team and due date
Note	

Number	15.2.3.3
Requirement	Requester team leader will be able to reject pending request
Note	• New status will be Draft • Request will be editable to requester/admin • Request will be assigned to requester • Requester will receive email notification

Number	15.2.3.4
Requirement	Requester team leader will be able to approve pending request
Note	• New status will be Approved • Request will be uneditable to requester team leader • Request will be editable to executing team leader/admin • Request will be assigned to executing team leader • Executing team leader will receive email notification

Number	15.2.3.5
Requirement	Executing team leader will be able to reject the request
Note	• New status will be Draft • Request will be uneditable to executing team leader • Request will be editable to requester/admin • Request will be assigned to requester

Confidential

	Requester will receive email notification
--	---

Number	15.2.3.6
Requirement	Executing team leader will be able to accept the request
Note	- New status will be Awaiting assignment - New field 'Performed by' is added

Number	15.2.3.7
Requirement	Executing team leader and admin member will be able to select a performer from executing team members
Note	- New status will be Assigned - Request will be editable to Performer/admin - Executing team leader will remain the assignee - Performer will receive email notification

Number	15.2.3.8
Requirement	Performer will be able to reject the request
Note	- New status will be awaiting assignment - Request will be uneditable to performer - Request will be editable to executing team leader/admin - Executing team leader will receive email notification

Number	15.2.3.9
Requirement	Performer will be able to start the request
Note	New status will be in progress

Number	15.2.3.10
Requirement	Performer user will be able to complete the request
Note	- New status will be completed - Request will be uneditable to performer

15.2.4 Actions in Request templates section

Number	15.2.4.1
Requirement	User will be able to add to requests section the following elements: - Text - Spreadsheet - Steps - Samples - Plate - Equipment - Sketch

	• Compound • Reaction • Code • Page break
Note	

16. Batch Management

16.1 Labguru Processes module

16.1.1 Processes index page

Number	16.1.1.1
Requirement	User will be able to add a new process
Note	

Number	16.1.1.2
Requirement	User will be able to view list of processes
Note	• User will be able to enter process page by clicking on its name • User will be able to enter project page by clicking on default project link • User will be able to enter member page by clicking on member link

Number	16.1.1.3
Requirement	User will be able to add new batch record for each process
Note	By clicking 'Add batch record'

16.1.2 Processes show page

16.1.2.1 General

Number	16.1.2.1.1
Requirement	User will be able to add a new batch record
Note	

Number	16.1.2.1.2
Requirement	User will be able to view process details: • ID • Item type • Created by • Created at • Owner • Updated at

Confidential

	• Naming convention • Indexed at • Threshold • Automated reports • Default project • Allow project select • Process description • Catalog number • Custom fields
Note	• User will be able to enter member page by clicking on created by/ created at • User will be able to enter project page by clicking on default project link

Number	16.1.2.1.3
Requirement	User will be able to edit process
Note	

Number	16.1.2.1.4
Requirement	User will be able to tag the process
Note	By clicking on 'add tag' in 'Tags' side box

Number	16.1.2.1.5
Requirement	User will be able to link process to other pages
Note	By clicking on 'add link' in 'Linked Resources' side box

Number	16.1.2.1.6
Requirement	User will be able to add attachment to process page
Note	By clicking on 'Upload a file' in 'Attachments' side box

Number	16.1.2.1.7
Requirement	User will see list of attached images/files
Note	

Number	16.1.2.1.8
Requirement	User will be able to download a file attached to process
Note	

Number	16.1.2.1.9
Requirement	User will be able to view an image attached to process in full size
Note	

Number	16.1.2.1.10
Requirement	User will be able to remove a file attached to process
Note	

Number	16.1.2.1.11
Requirement	User will be able to add a comment
Note	

Number	16.1.2.1.12
Requirement	User will be able to add a task to process
Note	

16.1.2.2 SOP tabs

Number	16.1.2.2.1
Requirement	User will be able to view the SOP used in process
Note	User will be able to enter SOP page by clicking on its name

16.1.2.3 Protocols tabs

Number	16.1.2.3.1
Requirement	User will be able to view protocols used in process
Note	

Number	16.1.2.3.2
Requirement	User will be able to add protocols from drop down list
Note	By clicking on ‘Manage protocols’

Number	16.1.2.3.3
Requirement	User will be able to remove process protocols
Note	- By clicking the trash icon/‘Manage protocols’ - Users will not be able to remove SOP protocols that are part of the process

16.1.2.4 Batch Records tabs

Number	16.1.2.4.1
Requirement	User will be able to view batch records created in process

Note	User will be able to enter batch record page by clicking on its name
------	--

Number	16.1.2.4.2
Requirement	User will be able to schedule an activity for a batch record
Note	

16.1.2.5 Schedule activity page

Number	16.1.2.5.1
Requirement	User will be required to name the activity
Note	

Number	16.1.2.5.2
Requirement	User will be required to select activity type
Note	From dropdown list of protocols associated with process

Number	16.1.2.5.3
Requirement	User will be required to select start date and time for the activity
Note	

Number	16.1.2.5.4
Requirement	User will be able to assign the activity to another members
Note	- By default the activity will be signed to the user himself - By typing the first 3 letters of member's name

Number	16.1.2.5.5
Requirement	User will be able to add remarks to the activity
Note	

16.1.2.6 Bill of Materials (BOM) & Equipment tab

16.1.2.6.1 Recipe material and equipment tables

Number	16.1.2.6.1.1
Requirement	User will be able to view list of material items of the process in required materials table
Note	

Number	16.1.2.6.1.2
--------	--------------

Confidential

Requirement	User will be able to view cost of process
Note	Cost value was calculated in bill of material step (refer section 16.1.3.4.7)

Number	16.1.2.6.1.3
Requirement	User will be able to view available quantity for production
Note	Available quantity is the total material of the item (units are identical as size units)

Number	16.1.2.6.1.4
Requirement	User will be able to view capacity for production
Note	Based on available quantity and required quantity values

Number	16.1.2.6.1.5
Requirement	User will be able to sort the required materials in ascending or descending order by any column header
Note	- By clicking on a header name - By clicking on the kendo menu and select one of the sorting options

Number	16.1.2.6.1.6
Requirement	User will be able to filter required materials by any combination of multiple columns
Note	

Number	16.1.2.6.1.7
Requirement	User will be able to add and remove columns from the display
Note	- By clicking on the kendo menu -> 'Columns' and checking the wanted columns - User preferences will be saved

Number	16.1.2.6.1.8
Requirement	User will be able to adjust the width of the columns
Note	- By dragging and dropping the columns' borders - Auto-adjust by double clicking the border

Number	16.1.2.6.1.9
Requirement	User will be able to view list of required equipment of the process
Note	

Number	16.1.2.6.1.10
Requirement	User will be able to view equipment's usage cost
Note	

Number	16.1.2.6.1.11
Requirement	User will be able to sort the required equipments in ascending or descending order by any column header
Note	• By clicking on a header name • By clicking on the kendo menu and select one of the sorting options

Number	16.1.2.6.1.12
Requirement	User will be able to filter required equipments by any combination of multiple columns
Note	

Number	16.1.2.6.1.13
Requirement	User will be able to add and remove columns from the display in required equipments table
Note	• By clicking on the kendo menu -> 'Columns' and checking the wanted columns • User preferences will be saved

Number	16.1.2.6.1.14
Requirement	User will be able to adjust the width of the columns in required equipments table
Note	• By dragging and dropping the columns' borders • Auto-adjust by double clicking the border

16.1.2.6.2 Notification section

Number	16.1.2.6.2.1
Requirement	User will be able to view number of batches that can be produced based on BoM recipe
Note	In 'Current Production Capacity'

Number	16.1.2.6.2.2
Requirement	User will be able to view estimated batch cost for a single batch production
Note	• For materials, equipments and both material& equipment

Confidential

	When a recipe includes multiple types of currency, the cost will be separated by each currency type
--	---

Number	16.1.2.6.2.3
Requirement	User will be able to view list of materials limiting production
Note	- When material quantity is below the threshold value specified for batch production - The production capacity at risk warning will display the number of materials with low stock levels

Number	16.1.2.6.2.4
Requirement	The user will be able to view cost breakdown for each currency type based on the following categories: - Material - Collection
Note	User will be able to select any combination of multiple currency

Number	16.1.2.6.2.5
Requirement	User will be able to view list of unmaintained equipments
Note	The production capacity at risk warning will display the number of equipment needing maintenance

Number	16.1.2.6.2.6
Requirement	The user will be able to view cost breakdown by equipment for each currency type
Note	User will be able to select any combination of multiple currency

16.1.2.7 Reports tab

Number	16.1.2.7.1
Requirement	User will be able to view reports created in process
Note	- User will be able to enter report page by clicking on its name - User will be able to enter batch record page by clicking on its link

Number	16.1.2.7.2
Requirement	User will be able to edit process
Note	By clicking 'Manage report template'

16.1.3 Add/Edit Process
16.1.3.1 Processes info step

Number	16.1.3.1.1
Requirement	User will be required to name the process
Note	Name must be unique

Number	16.1.3.1.2
Requirement	User will be required to select item type from dropdown list
Note	

Number	16.1.3.1.3
Requirement	User will be required to select process owner from dropdown list
Note	By default the user himself will be the owner

Number	16.1.3.1.4
Requirement	User will be able to add process description
Note	

Number	16.1.3.1.5
Requirement	User will be able to add catalog number associated with produced item
Note	

Number	16.1.3.1.6
Requirement	User will be able to save/ move to the next step
Note	

16.1.3.2 Batch record info step

Number	16.1.3.2.1
Requirement	User will be able to set naming conversation (Test-%ID-%d-%m-%Y-%H:%M) for the batch records
Note	Name format can include the following: • Test • Year (%Y) • Month (%m) • Day (%d)

Confidential

	• Hour (%H) • Minute (%M) • Second (%S) • Batch Record ID (%ID)
--	--

Number	16.1.3.2.2
Requirement	User will be able to add custom fields to process
Note	By clicking ‘Add new row’

Number	16.1.3.2.3
Requirement	User will be able to set a custom field as: • Regular field (text or numbers) • Only numbers • URL • Date • Date and Time • Pre-defined list • Text area • Pre-defined list with multiple select • Simple long text area • Account LOV
Note	

Number	16.1.3.2.4
Requirement	User will be able to mark a custom field as required
Note	

Number	16.1.3.2.5
Requirement	User will be able to delete custom field
Note	

16.1.3.3 Work plan & Activities

Number	16.1.3.3.1
Requirement	User will be able to select a project for the process
Note	• By clicking ‘Default project’ • By typing the first 3 letters of project’s name

Number	16.1.3.3.2
Requirement	User will be able to enable/disable project selection
Note	When enabled, the user can choose a project other than the default

Number	16.1.3.3.3
Requirement	User will be able to select SOP for the process
Note	By typing the first 3 letters of SOP's name

Number	16.1.3.3.4
Requirement	User will be able to mark SOP as required
Note	

16.1.3.4 Bill of Materials step

Number	16.1.3.4.1
Requirement	User will be able to view list of materials required for process
Note	Process required materials are inventory items that included in the protocols and are part of the process's SOP

Number	16.1.3.4.2
Requirement	User will be able to sort the materials in ascending or descending order by any column header
Note	• By clicking on a header name • By clicking on the kendo menu and selecting one of the sorting options

Number	16.1.3.4.3
Requirement	User will be able to filter materials by any combination of multiple columns
Note	

Number	16.1.3.4.4
Requirement	User will be able to add and remove columns from the display
Note	• By clicking on the kendo menu -> 'Columns' and checking the wanted columns • User preferences will be saved

Number	16.1.3.4.5
Requirement	User will be able to adjust the width of the columns
Note	• By dragging and dropping the columns' borders • Auto-adjust by double clicking the border

Number	16.1.3.4.6
Requirement	User will be able to add the required quantity and usage units
Note	• Required quantity is numeric field • User will be able to select volume/weight units from drop down list

Number	16.1.3.4.7
Requirement	User will be able to view calculated cost
Note	• Required quantity, usage units, price and size are required for the calculation • The size and price of are purchasable attributes of inventory item • Usage unit type must align with the size units in inventory • When the cost cannot be calculated, "Uncalculable will be displayed

Number	16.1.3.4.8
Requirement	User will be able to edit the following material values: • Required quantity • Usage units • Price and currency • Size and units
Note	• User will be able to select usage units, currency and size units from drop down list

Number	16.1.3.4.9
Requirement	The user will be able to specify low production capacity threshold
Note	Integer numbers only

Number	16.1.3.4.10
Requirement	The user will be able to see estimated batch cost
Note	When a recipe includes multiple types of currency, the cost will be separated by each currency type

16.1.3.5 Process reports

Number	16.1.3.5.1
Requirement	User will be able to design a template for a report
Note	By enabling 'Automated reports'

Number	16.1.3.5.2
Requirement	User will be required to name the report type
Note	

Number	16.1.3.5.3
Requirement	User will be required to add report layout

Note	By selecting from list of existing reports or creating a new report
------	---

Number	16.1.3.5.4
Requirement	User will be able to set naming conversation (Test-%ID-%d-%m-%Y-%H:%M) for the batch records
Note	Name must be unique and defined to include: • Test • Year (%Y) • Month (%m) • Day (%d) • Hour (%H) • Minute (%M) • Second (%S) • Batch Record ID (%ID)

Number	16.1.3.5.5
Requirement	User will be able to customize the report by adding the following types of sections: • Protocol
Note	• In protocol type, user will be able to select an existing protocol and data/section from created activities, to be copied to the generated report when triggered by a selected action

Number	16.1.3.5.6
Requirement	User will be able to add/remove new sections
Note	By clicking ‘Add new row or trash icon

Number	16.1.3.5.7
Requirement	User will be able to disable automated reports
Note	

16.1.3.6 Equipment Requirements

Number	16.1.3.6.1
Requirement	User will be able to view list of equipment required for process
Note	

Number	16.1.3.6.2
Requirement	User will be able to sort the required equipment in ascending or descending order by any column header
Note	• By clicking on a header name

Confidential

	By clicking on the kendo menu and selecting one of the sorting options
--	--

Number	16.1.3.6.3
Requirement	User will be able to filter equipment by any combination of multiple columns
Note	

Number	16.1.3.6.4
Requirement	User will be able to add and remove columns from the display
Note	- By clicking on the kendo menu -> 'Columns' and checking the wanted columns - User preferences will be saved

Number	16.1.3.6.5
Requirement	User will be able to add uses per step
Note	- When usage rate input is per use - Total operating time field will be uneditable

Number	16.1.3.6.6
Requirement	User will be able to add total operating time
Note	- When usage rate input is per hour - Total operation time dropdown list values: seconds/ minutes/hours/days/weeks - Uses per step field will be uneditable

Number	16.1.3.6.7
Requirement	User will be able to edit usage rate
Note	

Number	16.1.3.6.8
Requirement	User will be able to view calculated cost
Note	- Uses per step/Total operation time is required - When the cost cannot be calculated, "Uncalculable will be displayed

Number	16.1.3.6.9
Requirement	The user will be able to view maintenance status: maintained/not maintained/ not set
Note	

Number	16.1.3.6.10
Requirement	The user will be able to view estimated equipment usage cost
Note	When a recipe includes multiple types of currency, the cost will be separated by each currency type

16.2 Labguru Batch Records module

16.2.1 Batch Records index page

Number	16.2.1.1
Requirement	User will be able to add a new batch records
Note	User will be able to select an existing process or create a new process

Number	16.2.1.2
Requirement	User will be able to view list of batch records
Note	• User will be able to enter batch records page by clicking on its name • User will be able to enter project page by clicking on its link • User will be able to enter batch folder by clicking on its link • User will be able to enter process by clicking on its link

16.2.2 Batch Records show page

16.2.2.1 General

Number	16.2.2.1.1
Requirement	User will be able to view batch record details: • ID • Item type • Process • Project • Created by • Created at • Owner • Updated at • Status • Indexed at • Uuid • Storage • Description • Customize fields
Note	• User will be able to enter member page by clicking on created by/ owner

Confidential

	- User will be able to enter project page by clicking on project link - User will be able to enter process page by clicking on process link
--	--

Number	16.2.2.1.2
Requirement	User will be able to edit batch record
Note	

Number	16.2.2.1.3
Requirement	User will be able to tag the batch record
Note	By clicking on ‘add tag’ in ‘Tags’ side box

Number	16.2.2.1.4
Requirement	User will be able to link batch record to other pages
Note	By clicking on ‘add link’ in ‘Linked Resources’ side box

Number	16.2.2.1.5
Requirement	User will be able to add attachment to batch record page
Note	By clicking on ‘Upload a file’ in ‘Attachments’ side box

Number	16.2.2.1.6
Requirement	User will see list of attached images/files
Note	

Number	16.2.2.1.7
Requirement	User will be able to download a file attached to batch record
Note	

Number	16.2.2.1.8
Requirement	User will be able to view an image attached to batch record in full size
Note	

Number	16.2.2.1.9
Requirement	User will be able to remove a file attached to batch record
Note	

Number	16.2.2.1.10
Requirement	User will be able to add a comment
Note	

Number	16.2.2.1.11
Requirement	User will be able to add a task to batch record
Note	

Number	16.2.2.1.12
Requirement	User will be able to delete batch record
Note	

16.2.2.2 SOPs tab

Number	16.2.2.2.1
Requirement	User will be able to view the SOPs used in batch record
Note	User will be able to enter SOP page by clicking on its name

Number	16.2.2.2.2
Requirement	User will be able to create activities from SOP protocols
Note	• By clicking ‘Start SOP’ • When batch record status is ‘Pending’ • User will be required to enter Folder name • User will not be able to edit Project name • New status will be ‘In progress’

16.2.2.3 Used Materials tab

Number	16.2.2.3.1
Requirement	User will be able to view list of stock used to produce the batch
Note	

Number	16.2.2.3.2
Requirement	User will be able to sort stocks in ascending or descending order by any column header
Note	• By clicking on a header name • By clicking on the kendo menu and select one of the sorting options

Number	16.2.2.3.3
Requirement	User will be able to filter stocks by any combination of multiple columns
Note	

Number	16.2.2.3.4
Requirement	User will be able to add and remove columns from the display

Confidential

Note	- By clicking on the kendo menu -> 'Columns' and checking the wanted columns - User preferences will be saved
------	---

Number	16.2.2.3.5
Requirement	User will be able to adjust the width of the columns
Note	- By dragging and dropping the columns' borders - Auto-adjust by double clicking the border

Number	16.2.2.3.6
Requirement	User will be able to print used material PDF report
Note	- By clicking 'Generate report' - User will be able to select which columns to print from the index options (up to 8 columns)

Number	16.2.2.3.7
Requirement	User will be able to export used material report
Note	- By clicking 'Export' - User will be notified when the export is done

16.2.2.4 Stocks tab

Number	16.2.2.4.1
Requirement	User will be able to view list of stocks' batch record
Note	

Number	16.2.2.4.2
Requirement	User will be able to add new stocks to batch record
Note	- By clicking 'Add stock' - Limited to 5000 stocks

Number	16.2.2.4.3
Requirement	User will be able to edit selected stocks
Note	- By marking up to 100 stocks and clicking on 'Edit selected' - By clicking on the pencil icon to edit one stock

Number	16.2.2.4.4
Requirement	User will be able to tag stock/s
Note	By selecting stocks and clicking on 'Tag'

Number	16.2.2.4.5
Requirement	User will be able to mark stock/s as consumed

Note	• By selecting stocks and clicking on ‘Mark as consumed’ • By clicking on the stamp icon to consume one stock
------	--

Number	16.2.2.4.6
Requirement	User will be able to dispatch stock/s
Note	By selecting stocks and clicking on ‘Dispatch’

16.2.2.5 Used Equipment tab

Number	16.2.2.5.1
Requirement	User will be able to view list of equipment used in the batch record
Note	Based on adding a usage record to equipment in the batch activity

Number	16.2.2.5.2
Requirement	User will be able to sort used equipments in ascending or descending order by any column header
Note	• By clicking on a header name • By clicking on the kendo menu and select one of the sorting options

Number	16.2.2.5.3
Requirement	User will be able to filter used equipments by any combination of multiple columns
Note	

Number	16.2.2.5.4
Requirement	User will be able to add and remove columns from the display
Note	• By clicking on the kendo menu -> 'Columns' and checking the wanted columns • User preferences will be saved

Number	16.2.2.5.5
Requirement	User will be able to adjust the width of the columns
Note	• By dragging and dropping the columns' borders • Auto-adjust by double clicking the border

Number	16.2.2.5.6
Requirement	User will be able to print used equipment PDF report
Note	• By clicking ‘Generate report’ • User will be able to select which columns to print from the index options (up to 8 columns)

Number	16.2.2.5.7
Requirement	User will be able to export used equipment report
Note	• By clicking ‘Export’ • User will be notified when the export is done

16.2.2.6 Activities

Number	16.2.2.6.1
Requirement	Users will be able to view a list of protocols available for scheduling in the batch record
Note	By clicking ‘Protocols’ tab

Number	16.2.2.6.2
Requirement	Users will be able to schedule an activity
Note	• Refer section 16.1.2.5 • ‘Select activity’ field is not editable

Number	16.2.2.6.3
Requirement	User will be able to perform an activity
Note	• From list of activities scheduled for today in ‘Today’s Activities’ tab • New status will be ‘View’ • User will be able to enter an activity page by clicking on ‘View’

Number	16.2.2.6.4
Requirement	User will be able to re-schedule an activity
Note	• From list of planned activities in ‘Future Activities’ tab • Schedule and activities pop up will be displayed

Number	16.2.2.6.5
Requirement	User will be able to view activities history list
Note	User will be able to enter an activity history page by clicking on ‘View’

16.2.2.7 Datasets tab

Number	16.2.2.7.1
Requirement	User will be able to see vectors of datasets associated with the batch
Note	Each vector is linked to the original dataset of the vector

Number	16.2.2.7.2
Requirement	User will be able to export the list of vectors to a CSV file
Note	• By clicking the "Export to CSV" button • Export will be sent in email

Number	16.2.2.7.3
Requirement	User will be able to export the list of vectors to an Excel file
Note	• By clicking the "Export to Excel" button • Export will be sent in email • Export will include separate sheets for each dataset type, with each sheet displaying all the vectors associated with that specific dataset

16.2.2.8 Certificates tab

Number	16.2.2.8.1
Requirement	User will be able to print basic batch record report
Note	By clicking 'Generate report' in 'Basic Report' tab

Number	16.2.2.8.2
Requirement	User will be able to select which sections to print
Note	• By selecting 'Sections,' the user can choose either all sections or specific elements. • By default, all sections are selected

Number	16.2.2.8.3
Requirement	User will be able to fit content to page
Note	

Number	16.2.2.8.4
Requirement	User will be able to add header & footer to the report
Note	• By marking 'Set header&footer' • User can choose elements from a predefined list and position them on the left, right, or center of the page

Number	16.2.2.8.5
Requirement	User will be able to view custom batch record report
Note	• By clicking report's name link • Based on process settings

17. Customer Portal

17.1 General

17.1.1 Customer Portal First Login

Applicable to the following: Super Admin, Admin and users

Number	17.1.1.1
Requirement	User will be required to accept invitation when invited to customer portal
Note	

Number	17.1.1.2
Requirement	User will be required to set a strong password
Note	Password must be at least 11 characters long and may include a mix of lowercase letters, uppercase letters, numbers, and special characters

Number	17.1.1.3
Requirement	User will be required to confirm password
Note	

17.1.2 Customer Portal Login/logout

Number	17.1.2.1
Requirement	User will be required to enter email address and password to sign in
Note	

Number	17.1.2.2
Requirement	User will be able to reset his password in case it was forgotten
Note	- By clicking the 'Forgot password' link - Email with a link to reset password will be sent to user

Number	17.1.2.3
Requirement	User will be able to logout
Note	- By clicking sign out - Automatic time-out will occur after 30 minutes

17.2 User roles and access controls

17.2.1 Super Administrator index page

Number	17.2.1.1
Requirement	Super admin will be able to view list of admins
Note	- Details will include: - Admin ID - Name - Email - Super admin/ Admin - Invitation sent/ accepted at - Invited by - Clicking on the admin's line will redirect the super admin user to the admin show page

Number	17.2.1.2
Requirement	Super admin will be able to add admin users
Note	- By clicking 'Invite new admin' in the actions menu - Required fields: email and name - Optional field: bio/notes

Number	17.2.1.3
Requirement	Super admin will be able to edit admin users
Note	

Number	17.2.1.4
Requirement	Super admin will be able to delete admin users
Note	Confirmation notification will be displayed

Number	17.2.1.5
Requirement	Super admin will be able to sort the admin list
Note	By ID

Number	17.2.1.6
Requirement	Super admin will be able to change the number of the admin list displayed per page
Note	By selecting the number of items displayed per page (12, 24, 48 and 72)

17.2.2 Super Administrator show page

Number	17.2.2.1
Requirement	Super admin will be able to view admin details: • Admin ID • Name • Email • Bio notes • Super admin/non super admin • Invitation sent/ accepted at • Invited by
Note	Clicking on 'Invited by' link redirects to admin show page

Number	17.2.2.2
Requirement	Super admin will be able to edit admin users
Note	

Number	17.2.2.3
Requirement	Super admin will be able to delete admin users
Note	Confirmation notification will be displayed

17.2.3 Customers index page

Number	17.2.3.1
Requirement	Admin will be able to view all customers
Note	• Default view is alphabetically (ascending) • Details will include: ○ Customer ID ○ Name ○ Notes/Bio • Clicking on the customer's line will redirect to the customer show page

Number	17.2.3.2
Requirement	Admin will be able to edit customer
Note	

Number	17.2.3.3
Requirement	Admin will be able to delete customer
Note	• Confirmation notification will be displayed • Associated users will be removed • Associated samples will be removed

Number	17.2.3.4
Requirement	Admin will be able to sort the customer list
Note	By ID and groups users

Number	17.2.3.5
Requirement	Admin will be able to change the number of customers list displayed per page
Note	By selecting the number of items displayed per page (12, 24, 48 and 72)

Number	17.2.3.6
Requirement	Admin will be able to add new customer
Note	• By clicking ‘Create new customer’ • Name is required fields • Optional field: bio/notes

17.2.4 Customers show page

17.2.4.1 General

Number	17.2.4.1.1
Requirement	Admin will be able to view customer details: • Customer ID • Name • Notes • Created/ updated at • Groups • Users
Note	

Number	17.2.4.1.2
Requirement	Admin will be able to edit customer
Note	

Number	17.2.4.1.3
Requirement	Super admin will be able to delete customer
Note	• Confirmation notification will be displayed • Associated users will be removed • Associated samples will be removed

Number	17.2.4.1.4
Requirement	Admin will be able to add new group
Note	• By clicking ‘Add new group’ in the actions menu • Group name is required

Number	17.2.4.1.5
Requirement	Admin will be able to invite new user
Note	• By clicking 'Invite new user' in the actions menu • Required fields: email and name • Optional field: bio/notes • Refer 17.1.1 for login password requirements

Number	17.2.4.1.6
Requirement	Admin will be able to view list of groups
Note	• In 'Groups' tab • Details will include: ID and name

Number	17.2.4.1.7
Requirement	Admin will be able to view list of users
Note	• In 'Users' tab • Details will include: ○ ID ○ Name ○ Email ○ Invitation sent/ accepted by ○ Invited by

17.2.4.2 Groups page

Number	17.2.4.2.1
Requirement	Admin will be able to view groups details: • ID • Name • Customer • Users associated
Note	• By clicking group name

Number	17.2.4.2.2
Requirement	Admin will be able to delete group
Note	• Confirmation notification will be displayed • Associated samples will be deleted • Users assigned to group will not be able to access group data

Number	17.2.4.2.3
Requirement	Admin will be able to edit group name
Note	

Number	17.2.4.2.4
Requirement	Admin will be able to assign user to group/s
Note	• By selecting from list of uses • User can assign a single user/ multiple users by clicking 'Assign'/'Assign and Continue' • User that is not part of a group cannot view data related

Number	17.2.4.2.5
Requirement	Admin will be able to remove assigned user from group
Note	• By clicking unlink icon • Confirmation notification will be displayed

17.2.4.3 Users page

Number	17.2.4.3.1
Requirement	Admin will be able to view user details: • ID • Name • Email • Bio notes • Invitation sent/ accepted by • Invited by • Customer • Groups
Note	By clicking user name

Number	17.2.4.3.2
Requirement	Admin will be able to delete user
Note	Confirmation notification will be displayed

Number	17.2.4.3.3
Requirement	Admin will be able to edit user
Note	• Editable fields: name and bio • Email is non editable

Number	17.2.4.3.4
Requirement	Admin will be able to assign group to user
Note	• By selecting from list of groups • User can assign a single group/ multiple groups by clicking 'Assign'/'Assign and Continue' • User that is not part of a group cannot view data related

Number	17.2.4.3.5
Requirement	Admin will be able to remove assigned group from user
Note	• By clicking unlink icon • Confirmation notification will be displayed

17.3 Sample management

17.3.1 Samples index page

Number	17.3.1.1
Requirement	User will be able to view list of samples
Note	• Only samples assigned to customer and group they are part of • Clicking on the group/customer link will redirect the user to its show page • Notification of new/update samples will be highlighted

Number	17.3.1.2
Requirement	User will be able export samples to CSV file
Note	

Number	17.3.1.3
Requirement	User will be able to search sample
Note	By sample name and Labguru custom field text type

Number	17.3.1.4
Requirement	User will be able to sort samples
Note	By sample ID

Number	17.3.1.5
Requirement	User will be able to change the number of samples list displayed per page
Note	By selecting the number of items displayed per page (12, 24, 48 and 72)

Number	17.3.1.6
Requirement	User will not be able to delete/edit samples
Note	Also admin and super admin cannot delete samples

17.3.2 Samples show page

Number	17.3.2.1
Requirement	User will be able to view samples details
Note	Clicking on the group link will redirect to its show page

Number	17.3.2.2
Requirement	User will be able to download a file attached to sample
Note	

17.3.3 Shared Samples

Number	17.3.3.1
Requirement	User be able to share samples to the customer portal
Note	• A custom collection is created for sharing through the customer portal • Group filed is required for sharing the samples (cp_group_id) • User will receive a daily email notification (at 16:00 UTC) for new sample • Only users that are part of the group can view the samples

Number	17.3.3.2
Requirement	Admin user be able to select collection fields to share in the customer portal
Note	By selecting from prefixed custom fields list (with cp_)

Change Log

Version	Author	Date	Details	Approved by
5.0.1	Rotem Miara	16-DEC-2019	Initial Document	Gal Haase
5.0.8	Gal Haase	13-FEB-2020	Update to 5.0.8	Rotem Miara
5.0.9	Rotem Miara	25-MAY-2020	Update to 5.0.9	Gal Haase
5.0.10	Rotem Miara	02-SEP-2020	Update to 5.0.10	Gal Haase
5.0.11	Gal Haase	09-SEP-2020	Update to 5.0.11	Rotem Miara
5.0.12	Rotem Miara	21-OCT-2020	Update to 5.0.12	Gal Haase
5.0.13	Rotem Miara	17-DEC-2020	Update to 5.0.13	Gal Haase
5.14	Rotem Miara	20-JAN-2021	Update to 5.14	Gal Haase
5.15	Rotem Miara	22-FEB-2021	Update to 5.15	Gal Haase
5.16	Rotem Miara	08-MAR-2021	Update to 5.16	Gal Haase
5.17	Rotem Miara	19-MAY-2021	Update to 5.17	Gal Haase
5.18	Gal Haase	24-MAY-2021	Update to 5.18	Rotem Miara
5.19	Gal Haase	11-JUN-2021	Update to 5.19	Rotem Miara
5.20	Gal Haase	09-NOV-2021	Update to 5.20	Rotem Miara
5.21	Gal Haase	30-MAR-2022	Update to 5.21	Rotem Miara
5.22	Gal Haase	31-JUL-2022	Update to 5.22	Talia Hadad
5.23	Gal Haase	11-AUG-2022	Update to 5.23	Talia Hadad
5.24	Gal Haase	26-FEB-2023	Update to 5.24, Labguru version 5.96	Talia Hadad
5.25	Carmel Oren	31-JUL-2023	Update to 5.25, Labguru version 6.5	Gal Haase
5.26	Carmel Oren	29-Nov-2023	Update to 5.26, Labguru version 6.9	Gal Haase
5.27	Carmel Oren	28-Jan-2024	Update to 5.27, Labguru version 6.12	Gal Haase
5.28	Carmel Oren	15-May-2024	Update to 5.28, Labguru version 6.16	Gal Haase
5.29	Carmel Oren	14-Aug-2024	Update to 5.29, Labguru version 6.21	Gal Haase
5.30	Carmel Oren	13-Nov-2024	Update to 5.30, Labguru version 6.26	Gal Haase
5.31	Carmel Oren	03-Feb-2025	Update to 5.31, Labguru version 6.28	Gal Haase
5.32	Carmel Oren	21-May-2025	Update to 5.32, Labguru version 6.33	Gal Haase
5.33	Carmel Oren	12-AUG-2025	Update to 5.33, Labguru version 6.37	Gal Haase

Appendixes

Name	Version
Appendix A - User permission matrix	3.0